

Jul 10, 2019

## 08 Exploration of expression condition

DOI

[dx.doi.org/10.17504/protocols.io.5cfg2tn](https://dx.doi.org/10.17504/protocols.io.5cfg2tn)



Tjusls China<sup>1</sup>

<sup>1</sup>Tianjin University



**Tjusls China**

Tianjin University

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.5cfg2tn](https://dx.doi.org/10.17504/protocols.io.5cfg2tn)

**Protocol Citation:** Tjusls China 2019. 08 Exploration of expression condition. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5cfg2tn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 10, 2019

**Last Modified:** July 10, 2019

**Protocol Integer ID:** 25703



## Materials

### MATERIALS

⊗ IPTG **Bio Basic Inc. Catalog #IB0168.SIZE.100g**

⊗ SDS-PAGE Tricine Loading Buffer **Boster Bio Catalog #AR1143**

⊗ LB medium **Catalog #/**

### Before start

Set the gradient of condition to explore how to express it best. For example, we often use 0.5mM IPTG, 16°C/0.5mM IPTG, 37°C/1mM IPTG, 16°C/1mM IPTG, 37°C as different conditions.



- 1 Transform the plasmid into bacteria used to express target protein(e.g. E.coli BL21(DE3)).
- 2 Take monoclonal in the culture plate into LB tube and cultivate in shaking incubator overnight(10-12h) to activate bacteria.
- 3 Test the OD600 number of bacteria, then pipet 5-10ul into each new 5 mL LB tube. Don't forget to add antibiotic into tubes and mark them.
- 4 Cultivate in shaking incubator for 3-4 hours until the OD600 of bacteria range from 0.6 to 0.8.
- 5 Pipet 200ul bacterial liquid as uninduced sample, and take another 600ul to mix with 400ul 50%glycerol to store. Then add inducer IPTG into each tube in different concentration, and incubate at 16°C for 16 hours or at 37°C for 4 hours shaking at 200-300rpm.

600 µL
- 6 After cultivating, pipet 200ul for each as induced sample.  
(The method to make samples
  - 1) Centrifuge the taken bacterial liquid at 12,000rpm for 3 minutes
  - 2) Drop the supernatant and resuspend the precipitate using 100ul ddw.
  - 3) Pipet 50ul resuspending liquid to mix with 10ul 6XSDS Loading buffer
  - 4) Boil it in dry bath at 100°C for 10 minutes)

Use SDS-PAGE to check whether the target protein express or not and what the most suitable condition for its expression is.

200 µL