



Sep 21, 2025

1000x streptomycin stock (50 mg/ml)

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5k325v1b/v1>

Jordan Ward¹

¹UC, Santa Cruz



Jordan Ward

UC, Santa Cruz

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5k325v1b/v1>

Protocol Citation: Jordan Ward 2025. 1000x streptomycin stock (50 mg/ml). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8epv5k325v1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 19, 2025



Last Modified: September 21, 2025

Protocol Integer ID: 227729

Keywords: 1000x streptomycin stock

Funders Acknowledgements:

NIH

Grant ID: R35GM158317

Abstract

Our standard protocol for making streptomycin stock solutions



- 1 Weigh 5 grams of streptomycin (VWR/Avantor 97061-528) and add to a small beaker.
- 2 Add a magnetic stir bar and 50 ml of milliQ water and put on a magnetic stir plate. Mix until streptomycin is dissolved.
- 3 Put a sterile 0.22 μm syringe filter onto a 50 ml sterile syringe and use it to sterilize streptomycin stock, injecting it into a sterile 50 ml conical tube.
- 4 Make 1 ml aliquots and store at -20°C . Aliquots are good for at least one year.