

Jan 28, 2026

10X Genomics 5' scRNA-seq with sgRNA direct capture

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge1xe7v47/v1>

Boxun Li¹, Charles Gersbach¹

¹Duke University

Gersbach Lab



Andrea R Daniel

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge1xe7v47/v1>

Protocol Citation: Boxun Li, Charles Gersbach 2026. 10X Genomics 5' scRNA-seq with sgRNA direct capture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ge1xe7v47/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 26, 2026

Last Modified: January 28, 2026

Protocol Integer ID: 241578

Keywords: 10X Genomics, single-cell RNA-seq, neuron, astrocyte, sgRNA direct capture, primary astrocytes with sgrna direct capture, sgrna direct capture, seq of neuron mono, 10x genomic, mouse primary astrocyte, primary astrocyte, cell rna, neuron mono, seq

Funders Acknowledgements:

NIH

Grant ID: HG012053

Disclaimer

This protocol is adapted from the work of Boxun Li in the Gersbach lab at Duke University.

Abstract

This protocol describes 10X Genomics 5' HT v2 single-cell RNA-seq of neuron mono-culture or co-culture with mouse primary astrocytes with sgRNA direct capture.

- 1 Grow and dissociate neuron and astrocyte monocultures and neuron-astrocyte co-cultures (according to previous protocols Ngn2 neuron differentiation [Ref. 1], and Neuron dissociation [Ref. 2]).

For Perturb-seq, two different protocols are used.

1. For high-MOI experiments (i.e., the putative regulatory element (pRE) screens): iPSCs are seeded at 300k-500k/mL for 2mL/well on 6-well plates. Lentivirus expressing the sgRNA library is added at high MOI while cells are in suspension. iPSCs are transitioned to pre-differentiation as detailed in protocol Ngn2 neuron differentiation (Ref. 1) and pre-differentiated and matured as normal. Neurons are dissociated on Day 14 of maturation for scRNA-seq directly.
 2. For low-MOI experiments (i.e., the promoter screens): iPSCs are pre-differentiated as normal as detailed in protocol Ngn2 neuron differentiation (Ref. 1). On the day of reseedling (Maturation Day 0), pre-differentiated cells are seeded as detailed in the protocol. Lentivirus expressing the sgRNA library is added at low MOI while cells are in suspension. Neurons are dissociated for Fluorescence Activated Cell Sorting (FACS) first to isolate EGFP+ cells (expressed by lentivirus) before moving on to scRNA-seq.
- 2 For promoter Perturb-seq only: FACS isolate EGFP+ cells on a Sony Cell Sorter SH800Z.
- 3 Perform single-cell RNA-seq:
 1. Single cell suspensions are prepared according to the steps above and protocol Neuron dissociation (Ref. 2).
 2. Concentrations of single-cell suspensions are adjusted based on post-filtration counts and loaded onto a 10x Genomics Chromium Next GEM 5' HT v2 Chip N with 5' CRISPR kit to allow sgRNA capture, following the manufacturer's protocol (10X Genomics, PN-1000375, PN-1000451). Loading volumes are calculated by assuming effective concentrations ~1.2-fold lower than measured, to target ~20,000 cell recovery per lane with ~20% overloading (e.g., treating a cell suspension concentrated at 1000 cells/ μ L as 800 cells/ μ L).
 3. Single-cell gene expression and sgRNA expression libraries are prepared using the 10x Genomics Next GEM 5' HT v2 and 5' CRISPR kits following the manufacturer's protocols (10X Genomics, User Guide CG000512 Rev B, PN-1000374, PN-1000375, PN-1000451, PN-1000215).
 4. Final sequencing libraries are quality checked based on concentration (Qubit 4 Fluorometer; Qubit dsDNA High Sensitivity Assay Kit, Invitrogen Q32854) and fragment-size distribution (Agilent TapeStation 2200 or 4200; DNA D1000 and D5000 assays, Agilent 5067-5582, 5067-5583, 5067-5588 5067-5589) and are sequenced using Illumina Nextseq 2000, NovaSeq 6000 or NovaSeq X Plus sequencing to target



depths of $\sim 5 \times 10^4$ and $\sim 5 \times 10^3$ reads/cell for the gene expression and sgRNA expression libraries, respectively.

Protocol references

1. Li B, Gersbach CA. Excitatory neuron differentiation from inducible Ngn2 iPSCs (WTC11).
<https://www.protocols.io/view/excitatory-neuron-differentiation-from-inducible-n-ewov1kqxpgr2/v1>
2. Li B, Gersbach CA. Neuron dissociation protocol for single-cell sequencing applications.
<https://www.protocols.io/view/neuron-dissociation-protocol-for-single-cell-seque-j8nlk1or5g5r/v1>