

Apr 18, 2024

Version 1

# 10x Protocols: Chromium Next GEM Single Cell 5' -- University of Minnesota TMCsTMCs (CG000331 Rev E) V.1

DOI

[dx.doi.org/10.17504/protocols.io.5qpvokn89l4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvokn89l4o/v1)

Laura J Niedernhofer<sup>1</sup>, David A Bernlohr<sup>1</sup>

<sup>1</sup>University of Minnesota, Minneapolis, MN

Cellular Senescence Net...



Allie Abdellatif

MIBAM, University of Minnesota, Minneapolis, MN

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.5qpvokn89l4o/v1>

**Protocol Citation:** Laura J Niedernhofer, David A Bernlohr 2024. 10x Protocols: Chromium Next GEM Single Cell 5' -- University of Minnesota TMCsTMCs (CG000331 Rev E). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvokn89l4o/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 02, 2024

**Last Modified:** April 18, 2024

**Protocol Integer ID:** 97664

**Keywords:** minnesota genomics center, 10x genomic, chromium next gem single cell, barcoded cdna, 10x protocol, cdna, university of minnesota tmcstmc, university of minnesota tmc, nuclei isolation, minnesota tmcstmc, cell, minnesota tmc, following protocol, protocol, nuclei, microfluidic

**Funders Acknowledgements:**

NIH

Grant ID: 5U54AG076041-03

NIH

Grant ID: 5U54AG079754-02

## Abstract

Following single-cell dissociation or nuclei isolation, the Next GEM 5' assay uses microfluidics to partition and assign cell- or nuclei-specific barcodes to transcript cDNA at the 5' end. Barcoded cDNA is prepared with adaptors for sequencing by synthesis (SBS). The following protocol has been used at the University of Minnesota TMCs in collaboration with the University of Minnesota Genomics Center. The following protocol has been adapted from protocols developed by 10x Genomics and Illumina to be used at the University of Minnesota TMCs in collaboration with the University of Minnesota Genomics Center. These protocols are owned by their respective companies and are subject to periodic revision.


## Troubleshooting



## Tissue Preparation

- 1 Complete single cell or nuclei isolation prior to starting this protocol

## Library Preparation

- 2  CG000331\_Rev-E.pdf 4.3MB

### Note

Sequence with the read format 29,10,10,89

### Note

Sequencers used at UMN Genomics Center:

- Illumina NextSeq 2000
- Illumina NovaSeq 6000
- Illumina NovaSeq X Plus

## FASTQ Generation

- 3 BCL data from Illumina sequencer is demultiplexed and converted into FASTQ format using bcl2fastq version 2.20.0