



Jun 26, 2023

10x Protocols: Tissue Fixation & Dissociation for Chromium Fixed RNA Profiling -- University of Minnesota TMCs (CG000553 Rev A)

DOI

dx.doi.org/10.17504/protocols.io.q26g7y2pkgwz/v1

IOx Genomics¹, Laura Niedernhofer²

¹IOx Genomics, info@ioxgenomics.com; ²University of Minnesota Medical School, Minneapolis, MN USA

Cellular Senescence Net...



Mickayla DuFresne-To

University of Minnesota - Twin Cities

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7y2pkgwz/v1>

Protocol Citation: IOx Genomics, Laura Niedernhofer 2023. 10x Protocols: Tissue Fixation & Dissociation for Chromium Fixed RNA Profiling -- University of Minnesota TMCs (CG000553 Rev A). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q26g7y2pkgwz/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 03, 2023

Last Modified: June 26, 2023

Protocol Integer ID: 81379

Keywords: 10x, scRNAseq, Fixed, Dissociation, Fix-Dissociation, UMN, UMinnesota, University of Minnesota, preparation for 10x chromium flex scrnaseq, chromium fixed rna profiling, dissociation for chromium fixed rna profiling, 10x chromium flex scrnaseq, 10x genomics protocol, preparation for chromium single cell expression, 10x genomic, chromium single cell expression, chromium single cell, 10x protocol cg000553, tissue fixation, 10x protocol, rna, genomic, dissociation of whole tissue, whole tissue, single cell, 10x

Abstract

10x Genomics protocol where you fix the whole tissue and then dissociate it in preparation for Chromium Single Cell Expression flex

Protocol ID# (CG) and Revision letter provided here:

10x protocol CG000553, Revision A-- Fixation and dissociation of whole tissue; preparation for 10x Chromium Flex scRNAseq

Note: These protocols may not be the current version offered by the company but were used to produce the specific datasets connected to them. Please review the company support websites for the most recent versions of the protocols prior to starting your experiment.

Troubleshooting



1



10x Tissue Fixation+Dissociation, R...