

Aug 11, 2023

10x Protocols: Visium v1 FFPE Library Construction -- University of Minnesota TMCs (CG000407 Rev D)

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm2rpng3p/v1>

10x Genomics¹, Laura J Niedernhofer², David A Bernlohr²

¹10x Genomics, info@10xgenomics.com; ²University of Minnesota, Minneapolis, MN USA

Cellular Senescence Net...



Mickayla DuFresne-To

University of Minnesota - Twin Cities

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm2rpng3p/v1>

Protocol Citation: 10x Genomics, Laura J Niedernhofer, David A Bernlohr 2023. 10x Protocols: Visium v1 FFPE Library Construction -- University of Minnesota TMCs (CG000407 Rev D). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvm2rpng3p/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 05, 2023

Last Modified: August 11, 2023

Protocol Integer ID: 81460

Keywords: Spatial, Visium, UMN, UMinnesota, University of Minnesota, 10x Genomics, 10x, Niedernhofer, Bernlohr, FFPE, protocols from 10x genomic, 10x genomic, 10x genomics support website, visium spatial gene expression, 10x protocol for tissue preparation guideline, 10x protocol for imaging guideline, 10x protocol cg000407, visium v1 ffpe library construction, 10x protocol, genomic, ffpe sample, tissue preparation guideline, gene expression, cytassist component, 10x, chemistry on ffpe sample, recent versions of the protocol

Funders Acknowledgements:

NIH

Grant ID: 1U54AG07041-01

NIH

Grant ID: 1U54AG079754-01

Abstract

Protocols from 10x Genomics for Visium Spatial Gene Expression v1 chemistry on FPPE samples (without the CytAssist component).

1. Protocol ID# (CG) and Revision used:

10x protocol CG000407, **Revision D** -- Library Construction

2. Additional protocols added for guidance:

10x protocol for imaging guidelines

10x protocol for tissue preparation guidelines

Note: These protocols may not be the current version offered by 10x Genomics but were used to produce the specific datasets connected to them. Please review the 10x Genomics support website for the most recent versions of the protocols prior to starting your experiment.



1 **10x protocol CG000407, Revision D (Library construction)**

 CG000407_VisiumSpatialGeneExpr...

2 **Additional Protocols/Resources**

 10x Manual Visium FFPE Imaging G...

 10x Manual Visium Spatial FFPE Tis...

<https://www.10xgenomics.com/support/spatial-gene-expression-ffpe>