


Jul 10, 2019

## 14 SDS-PAGE

DOI

[dx.doi.org/10.17504/protocols.io.49qgz5w](https://dx.doi.org/10.17504/protocols.io.49qgz5w)



Tjusls China<sup>1</sup>

<sup>1</sup>Tianjin University



Tjusls China

Tianjin University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.49qgz5w>

**Protocol Citation:** Tjusls China 2019. 14 SDS-PAGE. protocols.io <https://dx.doi.org/10.17504/protocols.io.49qgz5w>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** July 08, 2019

**Last Modified:** July 10, 2019

**Protocol Integer ID:** 25616

**Keywords:** sd

## Guidelines

	SDS-PAGE Separation of gel samples	Scope of application
	6 % sample gel	50-150kD
	8 % sample gel	30-90kD
	10 % sample gel	20-80kD
	12 % sample gel	12-60kD
	15 % sample gel	10-40kD

## Materials

### MATERIALS

⊗ TEMED **Bio-Rad Laboratories Catalog #1610801**

⊗ 40g Acrylamide/ Bisacrylamide (29:1); Premixed Powder **G-Biosciences Catalog #786-506**

⊗ Coomassie brilliant blue G-250 **Bio Basic Inc. Catalog #CB0038.SIZE.100g**

⊗ SDS, 10%(w/v) solution **Bio Basic Inc. Catalog #SD8118.SIZE.100ml**

⊗ Tris pH 6.8

⊗ 10 % Ammonium persulfate **Catalog #/**

⊗ 1M Tris(PH=8.8) **Catalog #/**

## Troubleshooting

## Safety warnings

! Be careful of 30 % Acr-Bis M and TEMED owing to their potential harm.



## Before start

In general, 15 % sample gel is used, so here is only the composition of 15 % sample gel.


### **Composition of the 15 % separated of the gel sample(4mL):**

	Composition	volume (mL)
	Distilled water	0.4
	30 % Acr-Bis (29 : 1)	2
	1M Tris pH=8.8	1.5
	10 % SDS	0.04
	10 % Ammonium persulfate	0.04
	TEMED	0.002

### **Composition of concentrated sample gel(1mL):**

	Composition	volume (mL)
	Distilled water	0.7
	30 % Acr-Bis (29 : 1)	0.165
	1M Tris pH=6.8	0.125
	10 % SDS	0.01
	10 % Ammonium persulfate	0.01
	TEMED	0.001

- 1 Mount the gels in the vertical electrophoresis apparatus.
- 2 Load samples under the cathode buffer. Apply 10µL sample volumes to  $0.7 \times 5$  mm sample wells.

 10 µL
- 3 Set running conditions appropriate to your type of gel. We usually set 200V to run at the beginning, after 10 minutes change the voltage into 300V for another 30 minutes.
- 4 Protein can be visualized directly in the gel by Coomassie staining or silver staining. We usually choose Coomassie staining, then transfer the gel to the decolorizing solution and remove the dye solution.
- 5 Recording with gel imager.