

Aug 13, 2019

## 15 Determination of Enzyme Activity

DOI

[dx.doi.org/10.17504/protocols.io.6gthbwn](https://dx.doi.org/10.17504/protocols.io.6gthbwn)

Tjusls China<sup>1</sup>

<sup>1</sup>Tianjin University



Tjusls China

Tianjin University

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.6gthbwn](https://dx.doi.org/10.17504/protocols.io.6gthbwn)

**Protocol Citation:** Tjusls China 2019. 15 Determination of Enzyme Activity. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6gthbwn>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 13, 2019

**Last Modified:** August 13, 2019

**Protocol Integer ID:** 26867



## Materials

### MATERIALS

✕ Ultrasonic Processor **Cole-Parmer Catalog #UX-04714-52**

✕ DMSO **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8418**

✕ General 96-well plates (Black) **Catalog #/**






✕ Infinite M1000 Pro Automatic Microplate Reader **Catalog #/**

✕ Multi-channel adjustable pipette **Catalog #/**

✕ Fluorescent Probe(CDC-1) **Catalog #/**

✕ Target Enzyme(beta-lactamase) **Catalog #/**

## System setup – protein concentration

- 1 Soak the 96-well plates in 75% ethanol and put the container in ultrasonic cleaner for 30min to 1 hour, then use ddH<sub>2</sub>O to wash these plates several times. Put these clean plates in drying oven at 55°C.  
 00:30:00 ~  01:00:00  
 55 °C
- 2 Dilute the enzyme using its buffer. There we pipet 1 µL protein stock solution in 1mL buffer and mix gently. Then pipet 100 µL protein solution and mix with 400µL buffer each time, in order to dilute it as a 5-time gradient.
- 3 Pipet 94µL protein solution into 9 wells in plate, usually choosing B2-D4 area, to set 3 parallel controls. Pipet 94 µL buffer without protein into 3 wells as negative controls. Then add 6µL fluorescent substrate into wells.
- 4 Set up the program in Infinite M1000 Pro Automatic Microplate Reader. Shake for 10 sec at 654 rpm.  
 00:00:10  654 rpm  
Kinetic Cycle (to read fluorescent intensity each cycle)  
Fluorescent measure, 75 cycle, 10sec for each cycle.
- 5 Put the plate in Microplate reader, and click Start button.
- 6 When the facility ends testing, save data and import it into GraphPad Prism Software. Use “nonlinear fit” – “straight line” and compare R<sup>2</sup> of lines under different concentrations to pick up the best linear fit one, whose R<sup>2</sup> is most close to 1 .
- 7 Take this concentration as standard value, then set up parallel gradient of its 2x, 0.5x, 0.25x, etc. Repeat step 3-5.
- 8 Calculate the ratio of emission( $rE = Q_0/Q_m$ ,  $Q_0$  means the maximum fluorescent intensity of each reaction under different protein concentrations,  $Q_m$  means the maximum fluorescent intensity of all reactions under different protein concentrations). Use GraphPad Prism Software to calculate EC80 value. Set log(concentration of protein) as X, the rate of emission as Y. Use “nonlinear fit” – “log(agonist) vs. response—Find ECanything”, input 80 as the value of F parameter.
- 9 Usually we use the EC80 value as suitable protein concentration, and it can be adjusted according to the actual situation.

## System setup – buffer

- 10 Design experimental groups with the “N+(N-1) principle”.  
Since we use PBS as our protein buffer, and class B beta-lactamases are depend on  $\text{Zn}^{2+}$ , so we choose the concentration of NaCl, the concentration of  $\text{ZnCl}_2$ , and pH, as variables.
- 11 Repeat step 3-5 to measure.
- 12 When the facility ends testing, save data and import it into GraphPad Prism Software. Use “nonlinear fit” – “straight line” to calculate the initial velocity of each reaction a.k.a. its slope value. Choose the condition with higher initial velocity.

## Kinetic Constant Measurement

- 13 Dilute protein again with the ensured most suitable solution into proper concentration.
- 14 Dilute the fluorescent substrate as 2-time gradient for 8 groups.
- 15 Repeat step 3-5 to measure.
- 16 When the facility ends testing, save data and import it into GraphPad Prism Software. Use “nonlinear fit” – “straight line” to calculate the initial velocity of each reaction a.k.a. its slope value.
- 17 Use “nonlinear fit” – “Michaelis-Menten” to fit Michaelis plot of this beta-lactamase. At the same time the software will calculate kinetic constants  $K_m$ ,  $V_{max}$  automatically.
- 18 Dilute protein as 2-time gradient for several groups. Repeat step 3-5 to measure. Take FI as Y, [S] as X, then use “nonlinear fit” – “straight line” to calculate fluorescent calibration value.
- 19 Calculate  $k_{cat}$  value.  $K_{cat} = V_{max}/[E]$ .