

Apr 30, 2025

16S and 18S rRNA mock community construction

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5r19dg1b/v1>

Rob Lampe^{1,2}, Ariel Rabines^{1,2}, Bryce Ellman^{1,2}, Hong Zheng^{1,2}, Andrew E Allen^{1,2}

¹Scipps Institution of Oceanography; ²J. Craig Venter Institute

A.E. Allen Lab



Rob Lampe

Woods Hole Oceanographic Institution

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5r19dg1b/v1>

Protocol Citation: Rob Lampe, Ariel Rabines, Bryce Ellman, Hong Zheng, Andrew E Allen 2025. 16S and 18S rRNA mock community construction. protocols.io <https://dx.doi.org/10.17504/protocols.io.8epv5r19dg1b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 10, 2024

Last Modified: April 30, 2025

Protocol Integer ID: 98029

Keywords: mock communities, marine microorganisms, metabarcoding, sequencing, standards, external standards, amplicon sequencing, DNA, rRNA gene, rRNA gene sequences from diverse marine microorganism, rRNA mock community construction, rRNA gene sequence, rRNA mock community construction mock community, diverse marine microorganism, rRNA, mock community, amplicon, DNA for methods development, DNA

Funders Acknowledgements:

National Oceanic and Atmospheric Administration

Grant ID: NA19NOS4780181

National Science Foundation

Grant ID: OCE-2224726

Simons Foundation

Grant ID: 970820

Abstract

Mock communities are defined mixtures of cells or DNA for methods development, testing, or validation. Here we present mock communities based on previously used 16S and 18S rRNA gene sequences from diverse marine microorganisms. Their use as controls for amplicon sequencing enables detection of biased or aberrant runs.

Attachments



[16S Sequences.csv](#)

43KB



[18S Sequences.csv](#)

30KB



[18S Mock Community](#)

[R...](#)
15KB



[16S Mock Community](#)

[R...](#)
17KB

Materials

- Stock plasmids containing individual 16S or 18S sequences (see spreadsheet)
- PCR reagents: NEB OneTaq 2X Master Mix with Standard Buffer (Cat no. M0482)
- PCR and gel clean up kits: Zymo DNA Clean & Concentrator (D4013), NucleoSpin Gel and PCR Clean-up XS kit (Cat. no. 740611), Zymoclean Gel DNA Recovery Kit (Cat. no. D4001)
- Qubit dsDNA, High Sensitivity (Q32851 or Q32854)

The kits listed here simply describe the ones used for the creation of our mock communities and may be substituted for others that serve the same role.



Obtain plasmids

- 1 The stock plasmids were synthesized by Twist with the Kan High Copy vector and rehydrated to 40 ng/μL in low-EDTA TE Buffer. Use 1:10 dilutions of the stock plasmids (~4 ng/μL) for the first PCR.

First PCR of inserts

- 2 Set up  25 μL PCR reactions with M13 primers and NEB OneTaq 2X Master Mix with Standard Buffer:

Reagent	Volume (μL)
10 uM M13 Forward Primer	1
10 uM M13 Reverse Primer	1
Master Mix	12.5
Template DNA	1
Water	10.5

- 3 Run PCR with the following conditions:

Step	Temp (°C)	Time
Initial Denaturation	94	30 seconds
25 cycles	94	30 seconds
	55	1 minute
	68	2 minutes
Final Extension	68	5 minutes

- 4 Run PCR products on a gel to confirm successful amplification. Most samples are cleaned up with the Zymo DNA Clean-up and Concentrator Kit following the manufacturer's instructions. Perform the final elution in  50 μL of elution buffer and prepare 1:10 dilutions in  50 μL with molecular-grade H₂O for the next PCR.



18S Sample 14 (3323) and the archael sequences (6445, 6446, 6447) are gel extracted with either the Macherey-Nagel NucleoSpin Gel and PCR Clean-up XS kit (Cat. no. 740611) or Zymoclean Gel DNA Recovery Kit (Cat. no. D4001)

Second PCR of inserts

- 5 Set up PCR reactions as previously described with the following primers and annealing temperatures depending on the sequence type:

Type	Forward Primer	Reverse Primer	Annealing Temp (°C)
16S Bacteria	27F	1492R	55
16S Archaea	20F	1392R	54
18S	EukA	EukB	54

Primer	Sequence (5' - 3')
27F	AGAGTTTGATCMTGGCTCAG
1492R	TACGGYTACCTTGTTACGACTT
20F	AGTTTGATCMTGGCTCAG
1392R	ACGGGCGGTGTGTRC
EukA	AACCTGGTTGATCCTGCCAGT
EukB	GATCCTTCTGCAGGTTACCTAC

27F/1492R: Heuer, Holger, et al. "Analysis of actinomycete communities by specific amplification of genes encoding 16S rRNA and gel-electrophoretic separation in denaturing gradients." *Applied and environmental microbiology* 63.8 (1997): 3233-3241.

20F/1392R: Lane, D. J. "16S/23S rRNA sequencing." *Nucleic acid techniques in bacterial systematics* (1991).

EukA/EukB: Medlin, Linda, et al. "The characterization of enzymatically amplified eukaryotic 16S-like rRNA-coding regions." *Gene* 71.2 (1988): 491-499.

- 6 Perform clean-ups again with the Zymo DNA Clean-up and Concentrator Kit. Perform the final elution in  50 μ L . Here the archael sequences are gel extracted again.



Quantification, dilution, and mixing

- 7 Quantify the cleaned-up PCR product in duplicate with the Qubit dsDNA High Sensitivity (HS) assay kit.
- 7.1 (optional) Perform sanger sequencing with the second PCR primers for validation. This was performed with the construction of our mock communities.
- 8 Using the recipe spreadsheets, enter the concentrations, prepare dilutions, and mix accordingly. The expected ratios are also included in these spreadsheets. The final mock communities are also then serially diluted to approximately 0.01 ng/μL.
- 9 (optional) Confirm successful amplification following the PCR conditions and primers used for amplicon library preparation with both the diluted PCR product and mixtures. This was performed with the construction of our mock communities.

Protocol references

Parada AE, Needham DM, Fuhrman JA. 2016. Every base matters: assessing small subunit rRNA primers for marine microbiomes with mock communities, time series and global field samples. *Environ Microbiol* 18:1403-1414.

Yeh Y-C, McNichol J, Needham DM, Fichot EB, Berdjeb L, Fuhrman JA. 2021. Comprehensive single-PCR 16S and 18S rRNA community analysis validated with mock communities, and estimation of sequencing bias against 18S. *Environ Microbiol* 23:3240-3250.

Acknowledgements

We are grateful to Dr. Jed Fuhrman and members of his laboratory who previously shared the mock communities that the ones used here are based on.