

Jul 04, 2025

## 18% SPRI reagent

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l21mr4g1y/v1>

Nadin Rohland<sup>1</sup>

<sup>1</sup>Harvard Medical School



Nadin Rohland

Harvard Medical School

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.4r3l21mr4g1y/v1>

**Protocol Citation:** Nadin Rohland 2025. 18% SPRI reagent. protocols.io  
<https://dx.doi.org/10.17504/protocols.io.4r3l21mr4g1y/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** July 03, 2025



**Last Modified:** July 04, 2025

**Protocol Integer ID:** 221619

**Keywords:** ampure xp bead, reagent, brew reagent

## Abstract

This protocol explains how to prepare 18% SPRI beads that can be used instead of AMPure XP beads.

Note that as with any home-brew reagent, batches can show different behavior and should be tested (with a DNA ladder) prior to use.

## Guidelines

This reagent, if stored in the refrigerator and covered from light, can last for up to 2 years (and potentially longer) in our experience.

Keep the stock at 4C in the dark and make smaller working aliquots after the beads are well suspended.



## Materials

Magnetic racks for different tubes/plates:

- 50ml magnetic separation rack (e.g. NEB: S1507S or V&P Scientific: VP 772FB-2)
- 1.5/2.0ml magnetic separation rack (e.g. DynaMag-2 Magnet; ThermoFisher Scientific: 12321D)
- 96-well plate separation rack (ring magnet works well for automation, e.g. 96S Super Magne;, Alpaqua Engineering: A001322)
- 96-well plate separation rack (bar magnet or other magnets for 96 well plates work well for manual multichannel pipetting, e.g. Invitrogen: 12027 or Maxi Scientific: 4E-MRK100296)

1 Liter bottle, e.g. Nalgene Square PET Media Bottles with Closure, ThermoFisher Scientific: 342040-1000

|  | item                                                             | vendor            | part number            |
|--|------------------------------------------------------------------|-------------------|------------------------|
|  | Sera-Mag Speedbead carboxylate -modified [E3] magnetic particles | Cytiva            | 6515210505 0250 (15ml) |
|  | PEG-8000                                                         | Millipore-Sigma   | 89510                  |
|  | 5M NaCl                                                          | Millipore-Sigma   | S5150                  |
|  | 1M Tris-HCl                                                      | Invitrogen        | 15568025               |
|  | 0.5M EDTA                                                        | Millipore-Sigma   | 03690                  |
|  | water                                                            | Millipore-Sigma   | W4502                  |
|  | optional: GeneRuler Ultra Low Range DNA ladder                   | Thermo Scientific | FERSM1211              |

Prepare TE for bead washing and resuspension

|  | reagent     | volume | final concentration |
|--|-------------|--------|---------------------|
|  | 1M Tris-HCl | 10     | 10 mM               |
|  | 0.5M EDTA   | 2      | 1 mM                |
|  | water       | 988    |                     |
|  | TOTAL       | 1000   |                     |



## Safety warnings

! Mix Sera-Mag Speedbeads well before washing and SPRI reagent before using or aliquoting.

## Before start

This protocol is for 750ml 18% SPRI reagent (using a single bottle of 15ml of the Sera-Mag Speedbeads).  
If you make less, proportionally downscale.



## Wash beads

- 1 For a total of 750ml 18% SPRI reagent, use the entire 15ml bottle of Sera-Mag Speed beads carboxylate-modified and vortex vigorously to remove thick aggregates of beads; transfer into a 50ml Falcon tube.
- 2 Add 20 ml TE to rinse the storage bead bottle, mix and transfer to the Falcon tube; mix well by vortexing, then briefly centrifuge tube to collect liquid from inside the lid (e.g. by swinging your arm or short low centrifugation)
- 3 Place tube in magnet to collect beads on tube wall, this takes about 10 minutes or until liquid is clear
- 4 Remove supernatant with serological pipet and add between 20 and 35 ml TE, mix well to wash beads more by vortexing.
- 5 Repeat steps 3 and 4 for a total of 2 washes, and resuspend washed beads in 15ml TE.

## preparation of 18% SPRI reagent

- 6 Weigh out 135 g of PEG-8000 into a 1 Liter bottle.
- 7 To the PEG-8000 add below reagents

|  | reagent                  | volume in ml | final concentration  |
|--|--------------------------|--------------|----------------------|
|  | 5M NaCl                  | 150          | 1 M                  |
|  | 1M Tris-HCl              | 7.5          | 10 mM                |
|  | 0.5M EDTA                | 1.5          | 1                    |
|  | washed beads from Step 5 | 15           | 0.1% solid particles |

and fill up with water to 750 ml.

Mix well and let stand until all PEG particles are dissolved; intermittent shaking is helpful

- 8 Label with date and store in a dark box (or cover with aluminum foil) in the refrigerator.



- 9 This reagent is comparable to AmPure XP beads and can be used as recommended for the original product.  
(The recommended volume is 1.8x for standard PCR cleanup to retain dsDNA > 100bp and clean away shorter primers and adapters etc)

### Optional: preparing of 38% SPRI reagent

- 10 If a SPRI reagent is needed that allows to retain shorter fragments and to work in smaller volumes (e.g. to fit in a specific plate for automation), prepare a 38% SPRI reagent by following the above, but use 285 g of PEG-8000.  
The 38% SPRI reagent can be used in equal volume to the cleanup PCR reactions from e.g. ancient DNA library preparations.

### Recommended: testing of new reagent batch

- 11 Use 5  $\mu$ l (2.5  $\mu$ g total) of GeneRuler Ultra Low Range DNA ladder  
Add 95  $\mu$ l water  
add 180  $\mu$ l 18%SPRI (1.8x as recommended for AMPure XP)  
(or 100  $\mu$ l 38% SPRI)  
mix by vortexing or pipetting  
incubate 5 min at room temperature  
spin and place on appropriate magnet for 4 minutes or until liquid is clear (it takes longer for the more viscous 38%SPRI)  
wash beads while tubes remain on magnet with fresh 80% ethanol twice  
spin, put back on magnet and remove remaining ethanol with pipet  
air dry for 1min at room temperature  
elute in 20  $\mu$ l TE by resuspending the beads by pipetting  
spin, place on magnet and let separate for 1min  
collect elution into new tube
- 12 visualize on BioAnalyzer or TapeStation together with diluted original ladder of comparable starting amount  
18% SPRI: the 100 bp band should be less strong and the shorter fragments should be gone  
38% SPRI: the 50 bp should be as strong as before the cleanup and 35 bp should be more faint, 25 bp band should be gone.



## Protocol references

DeAngelis MM, Wang DG, Hawkins TL. Solid-phase reversible immobilization for the isolation of PCR products. *Nucleic Acids Res.* 1995 Nov 25;23(22):4742-3. doi: 10.1093/nar/23.22.4742. PMID: 8524672; PMCID: PMC307455.

Rohland N, Reich D. Cost-effective, high-throughput DNA sequencing libraries for multiplexed target capture. *Genome Res.* 2012 May;22(5):939-46. doi: 10.1101/gr.128124.111. Epub 2012 Jan 20. PMID: 22267522; PMCID: PMC3337438.

## Acknowledgements

I'd like to thank Matthias Meyer and his team at MPI-EVA for developing the 38% SPRI version.