

Dec 09, 2025

Version 1

# 1D Batch Processing and SAND Using NMRPipe V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3xen1g25/v1>

Christopher Esselman<sup>1</sup>, Arthur Edison<sup>1</sup>

<sup>1</sup>University of Georgia



Christopher Esselman

UGA

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.kqdg3xen1g25/v1>

**Protocol Citation:** Christopher Esselman, Arthur Edison 2025. 1D Batch Processing and SAND Using NMRPipe. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg3xen1g25/v1> Version created by **Christopher Esselman**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** November 30, 2025

**Last Modified:** December 09, 2025

**Protocol Integer ID:** 233829

**Keywords:** NMRPipe, SAND, NMRbox, nmr data, 1d batch processing, proton, sand, process, 1d

**Funders Acknowledgements:**

Network for Advanced NMR/NSF

Grant ID: 1946970

NIH/NIGMS R35

Grant ID: GM145348

## Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Abstract

This protocol describes how to batch process and SAND one-dimensional (1D) proton ( $^1\text{H}$ ) NMR data.


## Before start

This protocol assumes:

- You have a NMRbox account
- Directions for creating an account can be found at: <https://nmrbox.nmrhub.org/pages/getting-started>
- This protocol uses NMRPipe. Further Demo Data and processing tips can be found at: <https://www.ibbr.umd.edu/nmrpipe/demo.html>

## Acquire Necessary Scripts

- 1 Sign in to your NMRbox account.
- 2 Open the terminal, type `ssh`, and press Enter.




- 3 Paste the command provided below into the terminal window and press Enter.

### Command

Retrieve sand\_demo1 files from NMRPipe Website

```
wget https://www.ibbr.umd.edu/nmrpipe/sand_demo1.tar
```


- 4 Open the File Manager and extract sand\_demo1.tar to the home directory by right-clicking on sand\_demo1.tar and selecting "Extract Here."



- 4.1 The contents can also be extracted using the terminal by executing the command below:

```
tar -xvf sand_demo1.tar
```

- 5 Navigate to the extracted sand\_demo1 directory, then copy and paste each ".com" file into the directory that contains your unprocessed data.



- 5.1 The ".com" files can also be copied using the terminal:


```
cp ./*.com ~/destination_folder
```

## Creating Metadata table

- 6 Create a copy of meta.tab and place it into the directory containing your unprocessed data.
- 6.1 Update meta.tab to match the data within your unprocessed directory. Specifically, change the "DIR\_NAME" column (the second column of numbers) to match the directories of the spectra you wish to batch process. You should also change the "TITLE" column (the third column) to describe the spectra you are processing, ensuring that there are no spaces in the title name.

```
Open ▾ + meta.tab
~/sand_demo1
1
2 #
3 # Metadata for spike-in urine series.
4 # Required: INDEX DIR_NAME TITLE
5 # Optional: P0 P1
6
7 VARS INDEX DIR_NAME TITLE P0 P1
8 FORMAT %3d %10s %s %s %s
9
10 1 1 spike0-urine540ul-buf60ul-spike0ul Auto 0.0
11 2 2 spike1-urine540ul-buf60ul-spike20ul Auto 0.0
12 3 3 spike2-urine540ul-buf60ul-spike40ul Auto 0.0
13 4 4 spike3-urine540ul-buf60ul-spike60ul Auto 0.0
14 5 5 spike4-urine540ul-buf60ul-spike80ul Auto 0.0
15 6 6 ref-d200ul-buf60ul-spike20ul Auto 0.0
16 7 7 buf-d20540ul-buf60ul-spike0ul Auto 0.0
```

6.2 Modifying meta.tab by hand is cumbersome, but you can use the script below to automate the process.



1. Place the metaTabCreate.sh script in the same directory as the meta.tab file you just moved.
2. Ensure the file is executable by right-clicking on metaTabCreate.sh, choosing "Properties", selecting "Permissions," and checking "Allow this file to run as a program".
3. Using a terminal, navigate (cd) to the directory containing the files and enter the command

```
./metaTabCreate.sh
```

4. Provide a title for the dataset, ensuring you do not use spaces in the name.

6.3 Ensure that the rows in the meta.tab file contain metadata only for the files you intend to process. You should delete any rows corresponding to folders that do not contain data or should not be processed.

## Perform Processing

- 7 Using a new terminal window, navigate (cd) into the directory containing the unprocessed data, ".com" files, and meta.tab.
- 8 Enter the command

```
./all.com
```

The data processing is complete once a graphical image of the spectra appears on the screen. The final processed data will be located in the newly created folder named interp.

- 9 To view the spectra at any time after processing, open a terminal, navigate to the directory containing the unprocessed and processed data, and enter the command

```
specView.tcl -in interp/*.ft1
```

## Optimizing Processing

- 10 Sometimes, the automatic zero-order phasing does not phase the spectra perfectly. If this is the case, please follow the steps below.
- 11 Visualize the spectra to identify which spectra are not phased correctly. Then, reprocess the chosen spectra with 0.0 zero-order phasing and 0.0 first-order phasing by changing the "P0" values for those spectra in meta.tab from Auto to 0.0 and then reexecuting

```
./all.com
```


- 12 After processing is complete, open a new terminal, navigate to the interp folder, and enter the command

```
nmrDraw
```


Note that nmrDraw functions best when used with a three-button mouse. Further directions on spectral processing and use of NMRPipe can be found at:

<https://www.ibbr.umd.edu/nmrpipe/index.html>


- 12.1 Select the spectrum you want to phase by first right-clicking on "File" and then left-clicking on "Select File".




12.2 Use the top sliding bar to scroll through the spectra, and double left-click on the specific spectrum you wish to phase.



### 12.3 Use the sliders to phase the spectra by left-clicking and dragging them.



### 12.4 After the phasing looks correct, the values shown in the red box below will be the values used to replace the zero-order and first-order phase values for that spectrum in the meta.tab file.





- 13 After you have changed the phasing values in the meta.tab file, reexecute

```
./all.com
```

## SAND

- 14 Download the run\_sand and copy\_csv\_sand.sh files from the FLAnalyzer GitHub repository.

<https://github.com/edisonomics/FLAnalyzer/tree/main/sand>

- 15 Place the run\_sand file into your scratch directory on NMRbox. Your scratch directory is located at the path

```
/scratch/xxxx/username
```

, where xxxx is the number where your NMRbox username is stored (e.g., 0030). You can find this number by opening the terminal on NMRbox and typing

```
pwd
```

; the number appears immediately before your username.

An example output of executing the pwd command is: /home/nmrbox/0030/cesselman  
The corresponding scratch folder for this user is: /scratch/0030/cesselman

Some users do not have a number before their name.

An example output of executing the pwd command would be: /home/nmrbox/cesselman  
The corresponding scratch folder for this user would be: /scratch/cesselman

- 15.1 Ensure that the file is executable by using the same method described in Step 6.2.

- 16 Copy the interp folder created during batch processing into the scratch directory.

- 17 Edit the run\_sand script as follows:

```
1#!/bin/bash
2
3# Set the source directory on your scratch disk
4scratch_directory=""
5interp_directory=""
6
7# Set the destination directory on your local machine
8local_directory=""
9
10# Define the keyword to search for
11keyword="Done"
12
13# Iterate through the files in the source directory
14for file in "$interp_directory"/*; do
15    if [ -f "$file" ]; then
16        # Extract the file name without extension
17        filename=$(basename -- "$file")
18        filename_noext="${filename%.*}"
19
20        sand1D.com -in "$file" -x1 10ppm -xn -0.5ppm -simAll
21
22
23        # Copy the resulting folder to the local directory
24        mkdir "$local_directory/$filename_noext/"
25        cp -r "$scratch_directory/sand0001/"* "$local_directory/$filename_noext/"
26        rm -r "$scratch_directory/sand0001"
27    fi
28done
```

Line 4: Add the path to the scratch directory (e.g.,  
scratch\_directory="/scratch/0030/cesselman").

Line 5: Add the path to the interp folder in the scratch directory (e.g.,  
interp\_directory="/scratch/0030/cesselman/interp").

Line 8: Add the path to the local folder where SAND results will be saved. Ensure that this  
local folder already exists.  
(e.g., local\_directory="/home/nmrbox/0030/cesselman/Documents/").

Line 20: Change the ppm values to cover the full chemical shift range in which peaks  
appear  
(e.g., change 10ppm and -0.5ppm to values that suit your analysis)

Completion: After editing the file, save and close it.

- 18 On NMRbox, a single user can access up to three Virtual Machines (VMs) simultaneously. To reduce the time needed for SANDING, you should use three VMs. Repeat steps 14 through 17 on the different VMs. In the interp folder for each VM, delete any files that should not be SANDED by that specific VM. For example, if you want to SAND 140 spectra, keep the spectra in the scratch folder as follows: 1-48 on the first VM, 49-96 on the second VM, and 97-140 on the third VM.
- 19 On each VM, open the terminal, navigate to the scratch directory, and execute

```
./run_sand
```

## View SAND Results

- 20 In the local directory where the SAND results are saved, navigate to the folder of an individual Sanded spectrum. You can use the command below to assess the SAND results.

```
specView.tcl -in ft/*.ft1 sim/0001.ft1 simAll/1/*.ft1
```

## Gather SAND CSV Files

- 21 It is helpful to gather all the SAND tables into a single directory, and this is required for use with FLAnalyzer. Download the `copy_csv_sand.sh` file from the FLAnalyzer GitHub (Step 14) and ensure the permissions are set to allow the script to run as an executable.
- 22 Edit the `copy_csv_sand.sh` script as follows:

```
1#!/bin/bash
2
3# Define the output directory where files will be copied
4output_dir=""
5
6# Loop through each initial directory
7for initial_dir in /Change_to_path_destination_directory/*; do
8    if [ -d "$initial_dir" ]; then
9        # Get the name of the initial directory
10        dir_name=$(basename "$initial_dir")
11
12        # Find the 0001.ft1 file inside each initial directory
13        file_to_copy="$initial_dir/csv/0001.csv"
14
15        # Check if the file exists
16        if [ -f "$file_to_copy" ]; then
17            # Copy the file to the output directory with the desired name
18            cp "$file_to_copy" "$output_dir/$dir_name.csv"
19            echo "Copied $file_to_copy to $output_dir/$dir_name.csv"
20        else
21            echo "File $file_to_copy not found."
22        fi
23    fi
24done
```

Line 4: Add the absolute path to the directory where the CSV files will be copied.



Line 7: Change this to the absolute path of the destination directory for the SAND output folders, as specified in Line 8 of the run\_sand script (Step 17).

23 Execute the script by entering

```
./copy_csv_sand.sh
```

## Protocol references

1. Delaglio, F. et al. NMRPipe: a multidimensional spectral processing system based on UNIX pipes. *J. Biomol. NMR.* **6**, 277-293 (1995).
2. Wu, Y. et al. SAND: Automated Time-Domain Modeling of NMR Spectra Applied to Metabolite Quantification. *Anal Chem* **96**, 1843-1851 (2024).

## Acknowledgements

Thanks to Laura Morris for reviewing the protocol.