



Mar 23, 2023

Version 2

2-step PCR mixture and conditions (Barcoded-head primers for seqs pooling) V.2

DOI

dx.doi.org/10.17504/protocols.io.j8nlkky3xl5r/v2

Yin-Tse Huang¹, Tsu-Chun Hung¹

¹Kaohsiung Medical University



Yin-Tse Huang

Kaohsiung Medical University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkky3xl5r/v2>

Protocol Citation: Yin-Tse Huang, Tsu-Chun Hung 2023. 2-step PCR mixture and conditions (Barcoded-head primers for seqs pooling). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkky3xl5r/v2> Version created by **Yin-Tse Huang**

Manuscript citation:

Herbold CW, Pelikan C, Kuzyk O, Hausmann B, Angel R, Berry D, Loy A. 2015. A flexible and economical barcoding approach for highly multiplexed amplicon sequencing of diverse target genes. *Front. Microbiol.* [Internet] 6:731. Available from: <http://dx.doi.org/10.3389/fmicb.2015.00731>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 23, 2023

Last Modified: September 13, 2023

Protocol Integer ID: 79299

Keywords: pcr mixture, step pcr mixture, 2x supergreen pcr master mix, head primers for seq, pcr, head primer, primer, seq

Abstract

PCR mixture and condition (2X SUPERGREEN PCR MASTER MIX)

Troubleshooting

- 1 Wear glove, clean up the working bench w. 1% bleach

For 1' PCR head-primers




- 2 Prepare 1' PCR master mixutre for **head-primers (prepare 1.2X of solutions for pipetting error if needed)**

PCR mixture for head-primers for each reaction

	A	B	C	D
	Component	Volume	Volume (1.2X)	Final conc.
	Forward Primer (10 µM)	1.6 µl	1.9 µl	1 µM
	Reverse Primer (10 µM)	1.6 µl	1.9 µl	1 µM
	2X Supergreen PCR Master Mix	7.8 µl	9.4 µl	-
	ddH2O	4.1 µl	4.9 µl	-
	Total volume	15 µl	18 µl	-

Note


Negative control ALWAYS NEEDED! For example, if you have 5 PCR reactions to run, prepare master mixture for 6 reactions (5 DNA template + 1 negative control).

- 3 Mix the 1' PCR master mixture gently by pippeting. Quick spin the tube.
- 4 Transfer  15 µL 1' PCR master mixutre in 8-strip PCR tubes.
- 5 Add  0.6 µL DNA template in 8-strip PCR tubes, resulting in a  15.6 µL reaction mixture for 1' PCR.





Note

Negative control contains only  15 μL master mixture but not DNA template

- 6 Mix the reaction mixture gently by tapping the tubes. Quick spin the tubes.
- 7 Carry out PCR using the following condition:

1' PCR condition for **head-primers**

	A	B	C	D
	Step	Temp	Sec	Cycle
	<i>Initial denaturation</i>	95 °C	180	
	<i>Denaturation</i>	98 °C	30	20-25 cycles
	<i>Annealing</i>	60-66 °C varied (b)	30	
	<i>Extension</i>	72 °C	180	
	<i>Final extension</i>	72 °C	210	
	<i>Preservation</i>	Preservation	4 °C	∞

b. Annealing varied, **60-66°C** is working; Refer to 1' PCR primers for annealing temperature

c. 1kb ~ 1min extension; enough time allow full extension of sequence

7.1 1' hear-primers used in Huang lab


	A	B	C	D
	Name	Sequence	Tm°C	CG%
	NS1B1ngs_H1	GCTATGCGCGAGCTGCcctngttgatyctg ccagt	71.7	60
	LR5_H1	GCTATGCGCGAGCTGCTcctgagggaaac ttcg	70.2	60.6
	EF1-526F_H1	GCTATGCGCGAGCTGCgtcgtgytytygg hcaygt	71	59.3
	EF1-2218R_H1	GCTATGCGCGAGCTGCatgacaccracrg cracrgtytg	72.2	60.3

- 8 Carry out **electrophoresis** for inspection of DNA products



Gel before markdown

- 9 Markdown wells and upload the pictures to the Lab Google drive



Marked gel picture go to the Lab Google drive



For 2' PCR barcoded-head primers







- 10 Prepare 2' PCR master mixutre for **barcoded-primers** (prepare 1.2X of solutions for pipetting error if needed)

PCR mixture for barcoded-primers for each reaction **(NO PRIMERS at this point!!)**

	A	B	C	D
	Component	Volume	Volume (1.2X)	Final conc.
	2X Supergreen PCR Master Mix	10.75 µL	12.9 µL	-
	ddH2O	10.75 µL	12.9 µL	-
	Total volume	21.5 µL	25.8 µL	-

Note

Negative control ALWAYS NEEDED! For example, if you have 5 PCR reactions to run, prepare master mixture for 6 reactions (5 DNA template + 1 negative control).

- 11 Mix the 2' PCR master mixture gently by pippeting. Quick spin the tube.
- 12 Transfer  21.5 µL of the 2' PCR master mixture to 8-strip PCR tubes.
- 13 Add  2.5 µL **pre-mixed barcoded-head primers** (Forward + Reverse) to each PCR tubes.
- 14 Add  1 µL of **1' PCR product as template**, resulting in  25 µL reaction mixture for 2' PCR. 
- Negative control** contains only  24 µL master mixture and premixed barcoded-head primers but not DNA template
- 15 Mix gently by tapping the tubes. Quick spin the tubes.

16 Carry out 2' PCR using the following condition:


2' PCR condition for **barcoded-head primers**

	A	B	C	D
	Step	Temp	Sec	Cycle
	<i>Initial denaturation</i>	98 °C	30	
	<i>Denaturation</i>	98 °C	15	10-15 cycles
	<i>Annealing</i>	64-68 °C varied (a)	15	
	<i>Extension</i>	72 °C	30 (b)	
	<i>Final extension</i>	72 °C	210	
	<i>Preservation</i>	Preservation	4 °C	∞

a. Annealing varied, **65 °C** is working based on test on 220531; Refer 2' PCR primers for annealing temperature


b. 1kb ~ 1min extension; enough time allow full extension of sequence

17 Carry out **electrophoresis** for inspection of DNA products



Gel before markdown

18 Markdown wells and upload the pictures to the Lab Google drive



Marked gel picture go to the Lab Google drive