



Apr 17, 2025

2D and 3D Neuron/Microglia co-culture systems

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4zkzjvo5/v1>

Lucas Blasco-Agell^{1,2}, Veronica Testa^{1,2}, Valentina Baruffi^{1,2}, Antonella Consiglio^{1,2}, Miquel Vila³

¹Department of Pathology and Experimental Therapeutics, Bellvitge University Hospital- IDIBELL, 08908 Hospitalet de Llobregat, Spain;

²Institute of Biomedicine of the University of Barcelona (IBUB), Carrer Baldori Reixac 15-21, Barcelona 08028, Spain.;

³VHIR-CIBERNED-ASAP

Vilalab Public



Miquel Vila

VHIR-CIBERNED-ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4zkzjvo5/v1>



Protocol Citation: Lucas Blasco-Agell, Veronica Testa, Valentina Baruffi, Antonella Consiglio, Miquel Vila 2025. 2D and 3D Neuron/Microglia co-culture systems. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4zkzjvo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 08, 2025

Last Modified: April 17, 2025

Protocol Integer ID: 126365

Keywords: microglia, microglia interaction, confocal microscopy, 3d neuron, ubiquitous phosphoglycerate kinase, encoding green fluorescent protein, control of the ubiquitous phosphoglycerate kinase, green fluorescent protein, gfp fluorescence, cell, ipsc, neuron

Funders Acknowledgements:

ASAP

Grant ID: ASAP-020505

Abstract

To visualize iPSC-derived hMG within 2D and 3D co-culture systems, iPSCs are transduced with a lentiviral vector (LV) encoding Green Fluorescent Protein (GFP) under the control of the ubiquitous Phosphoglycerate Kinase (PGK) promoter. After transduction, cells are detached and GFP fluorescence is assessed 72 hours later. Neuron–microglia interactions are examined in a 2D co-culture model using iPSC-derived vmDAn, and confocal microscopy is used to capture images.

- 1 In order to visualize iPSCs-derived hMG within 2D/3D co-cultures, iPSCs are transduced with a lentiviral vector (LV) expressing Green Fluorescence Protein (GFP) under the Phosphoglycerate Kinase (PGK) ubiquitous promoter.
- 2 iPSCs are detached and 10^6 iPSCs are infected in 100 μ L of mTeSR™ within a 15 mL Falcon tube with a mean of infection (MOI) of 10.
- 3 Cells are incubated under normal conditions for 1 hour and plated in a 6 wp with 1 mL of mTeSR ON. After 16-18 hours, 1 mL of mTeSR is added.
- 4 Fluorescence is checked after 72 hours under an inverted microscope. To study neuron/microglia interactions, we generate 2D co-culture system with iPSC-derived vmDAn. First, mature GFP-iPSC-derived hMG are conditioned with medium from vmDAn without B27 vitA.
- 5 One week later, hMG are detached with accutase, counted and plated on top of day 42 vmDAn on a human astrocyte feeder layer, at a ratio of 1 hMG per every 2 neurons.
- 6 Astrocytes feeder layer is employed for every analysis where we need to assess neuronal health. Co-culture medium consists of RPMI with Neurobasal 50/50% supplemented with 1% Glut, 1% P/S, M-CSF, IL-34, BDNF, GDNF, TGF- β 3, cAMP and AA.
- 7 Half of the medium is changed every 2-3 days. vmDAn and hMG are co-cultured for 1 week in hypoxia conditions, before being used for experiments. For culture treatments, 0.5 μ g/mL of NM is added to co-cultures for 8 hours, and after 2 hours IVM 3 μ M is added, for a total of 6 hours. In case of PFF, 10 μ g/mL are added for 8 hours.
- 8 Images are taken under a 63x objective employing a confocal microscopy. vmDAn morphology and neurite complexity are assessed by ICC staining of TH and analyzed by Sholl analysis using the Simple Neurite Tracer plugin from ImageJ. hMG morphology is analyzed for circularity employing ImageJ.
- 9 Reconstruction of vmDAn and hMG are performed with IMARIS. Surface of vmDAn is generated with a Smooth of 0.103 μ m and a Background subtraction of 5 μ m.
- 10 Surface of hMG was generated with a Smooth of 0.103 μ m and a Background subtraction of 3 μ m. For the 3D co-culture system, one vmO is manually placed in a well of a U 96-wp with 50,000 GFP-hMG and incubated under normal conditions in an orbital shaker for 24 hours.
- 11 Then, every organoid is collected and plated on low-adhesion 6-wp. Coculture medium is replaced every 2-3 days. vmO and hMG are co-cultured for 1 week in normal incubation conditions in agitation, before being used for experiments.



- 12 100 ng/mL of LPS (Sigma-Aldrich) is added for 72 hours to co-cultures to study hMG morphological changes. Reconstruction of hMG is done with IMARIS.