

Apr 12, 2024

3,3-Diaminobenzidine Tetrahydrochloride (DAB) Immunohistochemistry on Mouse Brain Tissue Sections

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl497zjgo5/v1>

madalynn.erb Erb¹

¹Van Andel Research Institute

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl497zjgo5/v1>

Protocol Citation: madalynn.erb Erb 2024. 3,3-Diaminobenzidine Tetrahydrochloride (DAB) Immunohistochemistry on Mouse Brain Tissue Sections. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl497zjgo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: April 08, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 98141

Keywords: ASAPCRN, immunohistochemistry on mouse brain tissue section, floating mouse brain section, mouse brain tissue section, mouse brain section, mouse brain, diaminobenzidine tetrahydrochloride, dab, immunohistochemistry, mouse

Abstract

This protocol details the staining for free floating mouse brain sections.

Materials

 Pyrex® Staining Dish **Ted Pella Inc. Catalog #36754-60**

 8-Section Staining Nets **Ted Pella Inc. Catalog #36154-64**

 VECTASTAIN® Elite® ABC-HRP Kit, Peroxidase (Standard) **Vector Laboratories Catalog #PK-6100**

 Pierce™ DAB Substrate Kit **Thermo Scientific Catalog #34002**

 Fisherbrand™ Superfrost™ Plus Microscope Slides **Fisher Scientific Catalog #12-550-15**

 Entellan **Electron Microscopy Sciences Catalog #14802**

Recipes

10X TBS

	A	B
	1M Tris	500 ml
	5M NaCl	300 ml
	H ₂ O	200 mL

1X TBS + 0.1% triton → 1 liter

	A	B
	10X TBS	100 ml
	Triton X	1 ml
	Water	899 ml

TBS → 1 liter

	A	B
	1M Tris	50 ml
	5M NaCl	30ml
	Water	919 ml

TBST (0.1% Tween-20) → 1 liter



	A	B
	10X TBS	100 ml
	Tween-20	1ml
	Water	899 ml

0.2M Phosphate Buffer

	A	B
	Na ₂ HPO ₄	22.72 g
	NaH ₂ PO ₄	5.52 g
	H ₂ O	1000 mL

0.1M Phosphate Buffer

	A	B
	0.2M PB	500 mL
	H ₂ O	500 mL

Cresyl Violet (500 mL)

	A	B
	H ₂ O	500 mL
	Cresyl Violet Acetate	0.5g
	100% Glacial Acetic Acid	1.25mL

Warm to 60°C and stir until completely dissolved filter with whatman paper (wrap in foil to protect from light)

Quenching solution→100ml



	A	B
	30% hydrogen peroxide (Sigma H1009)	1 ml
	Methanol	99 ml

ABC mix→60ml

	A	B
	TBS	60 ml
	Reagent A	24 drops
	Reagent B	24 drops

DAB Solution→30ml

	A	B
	PBS	30 ml
	(DAB reagent 1)	12 drops buffer pH = 7.5
	(DAB reagent 2)	24 drops DAB
	(DAB reagent 3)	12 drops hydrogen peroxide

vortex between reagents

Blocking solution→100ml

	A	B
	NGS	10 ml
	TBS + triton	90 ml + 0.1%

Primary Antibodies:

	A	B	C	D
	Target	Species	Conc	Manufacturer



	A	B	C	D
	Iba1	Rabbit	1:1000	WAKO 019-19741
	GFAP	Mouse	1:1000	Sigma G3893
	pSer129-Synuclein	Rabbit	1:1000	Abcam ab51253
	TH	Rabbit	1:2000	Novus Biological N300109
	pSer202, pThr205 Tau (AT8)	Mouse	1:1000	Thermofisher MN1020
	GFP	Rabbit	1:1000	Thermofisher A-11122

Secondary Antibodies:

	A	B	C	D	E
	Target	Species	Conc	Manufacturer	
	Mouse	Goat	1:500	BA-9200-1.5	Biotinylated
	Rabbit	Goat	1:500	BA-1000-1.5	Biotinylated



Day 1

2d 1h 35m

- 1 Staining protocol for 35µm free floating mouse brain sections.
- 2 Staining is performed in glass staining dishes (Pyrex 36754-60) using 8-section staining nets (Ted Pella 36154-64) .
 - 2.1 The sections can be transferred between wells using a paint brush.
 - 2.2
 - Volume of solution:
 1. 20 mL - 30 mL for antibodies or ABC mix .
 2. 50 mL for washing.
- 3 Gently rock the plates during the washing and incubation steps.
- 4 Wash sections: 3 times (5 min per wash) in Tris Buffer Saline (TBS) at Room temperature .
 - 4.1 Wash sections: Wash for 00:05:00 in Tris Buffer Saline (TBS) at Room temperature (1/3).
 - 4.2 Wash sections: Wash for 00:05:00 in Tris Buffer Saline (TBS) at Room temperature (2/3).
 - 4.3 Wash sections: Wash for 00:05:00 in Tris Buffer Saline (TBS) at Room temperature (3/3).
- 5 Quenching: Incubate sections in Quenching Solution (0.3% H₂O₂ in Methanol) for 00:10:00 at 4 °C .
- 6 Wash sections: 3 short rinses in TBS then 2 times 5min in TBS at Room temperature .



5m



5m



5m



10m





6.1 Then wash for  00:05:00 in TBS at  Room temperature (1/2).

5m



6.2 Then wash for  00:05:00 in TBS at  Room temperature (2/2).

5m



7 Block sections in blocking solution (10% Normal Goat Serum (NGS) + 0.1% Triton-X in TBS) for  01:00:00 at  Room temperature .

1h



8 Primary Antibodies: Dilute antibodies in 5% NGS + 0.1% Triton-X in TBS.

8.1 Incubate in primary antibody solution for  48:00:00 at  4 °C .

2d



8.2 Cover dishes with parafilm for this step.

Day 4

1d 3h 47m

9 Wash 3 times (10 min each wash) in TBST (0.1% Tween 20 in TBS) at  Room temperature .



9.1 Wash for  00:10:00 in TBST (0.1% Tween 20 in TBS) at  Room temperature (1/3).

10m



9.2 Wash for  00:10:00 in TBST (0.1% Tween 20 in TBS) at  Room temperature (2/3).

10m

9.3 Wash for  00:10:00 in TBST (0.1% Tween 20 in TBS) at  Room temperature (3/3).

10m



10 Secondary Antibodies:



10.1 Dilute secondary antibodies in TBST.



10.2 Biotinylated goat anti-rabbit IgG (1:500 dilution) or Biotinylated goat anti-mouse IgG (1:500 dilution).

10.3 Incubate for 02:00:00 at Room temperature or 24:00:00 at 4 °C .

1d 2h



11 Wash sections 3 times (10min each wash) in TBST at Room temperature .

11.1 Wash for 00:10:00 in TBST at Room temperature (1/3).

10m

11.2 Wash for 00:10:00 in TBST at Room temperature (2/3).

10m

11.3 Wash for 00:10:00 in TBST at Room temperature (3/3).

10m

12 Prepare ABC solution using VECTASTAIN® Elite® ABC-HRP Kit, Peroxidase (Standard) (PK-6100).

12.1 ABC mix: 5 mL TBS + 2 drops of reagent A (vortex) + 2 drops reagent B (vortex) let stand at Room temperature for 00:30:00 before using.

30m



12.2 Vortex between reagents (after adding reagent A and after adding reagent B).



12.3 Make 30 mL ABC solution per dish.

13 Incubate sections in ABC solution for 01:00:00 at Room temperature .

1h



14 Wash sections 3 times (10min each wash) in TBST at Room temperature .



14.1 Wash for 00:10:00 in TBST at Room temperature (1/3).

10m



14.2 Wash for 00:10:00 in TBST at Room temperature (2/3).

10m



14.3 Wash for 00:10:00 in TBST at Room temperature (3/3).

10m



15 Wash sections 2 times (5 min each wash) in TBS at Room temperature .

15.1 Wash sections for 00:05:00 in TBS at Room temperature (1/2).

5m



15.2 Wash sections for 00:05:00 in TBS at Room temperature (2/2).

5m



16 Prepare DAB solution using DAB Substrate Kit, Peroxidase (HRP), with Nickel, (3,3'-diaminobenzidine) (SK-4100).



1. 5 mL TBS + 2 drops of Reagent 1 (Buffer pH=7.5) + 4 drops Reagent 2 (DAB) + 2 drops Reagent 3 (H₂O₂).
2. Vortex between reagents (after reagent 1, after reagent 2 and after reagent 3).
3. Make 30 mL DAB solution per dish.

17 Revelation with DAB chromogen (until cells are visible):

2m

Note

typically about 00:02:00 for TH.

- All dishes must be developed for the same amount of time – do them together, ask a friend for help if you have more than 2 dishes of tissue.

18 Perform 3 short rinses in TBS at Room temperature .



18.1 Perform rinse in TBS at Room temperature (1/3).





18.2 Perform rinse in TBS at Room temperature (2/3).



18.3 Perform rinse in TBS at Room temperature (3/3).



19 Wash sections 5 time (5 min each wash) in TBS at Room temperature .



19.1 Wash sections for 00:05:00 in TBS at Room temperature (1/5).

5m



19.2 Wash sections for 00:05:00 in TBS at Room temperature (2/5).

5m



19.3 Wash sections for 00:05:00 in TBS at Room temperature (3/5).

5m



19.4 Wash sections for 00:05:00 in TBS at Room temperature (4/5).

5m



19.5 Wash sections for 00:05:00 in TBS at Room temperature (5/5).

5m



20 Mount sections:

20.1 Use a petri dish filled with PBS to mount the sections on superfrost plus slides (Fisher Scientific 12-550-15).

- This is easiest to do using a medium sized paint brush.

20.2 After mounting let slides dry at Room temperature (00:02:00 - 00:10:00) then dip in miliQ H₂O.

12m



20.3 Dry slides Overnight at Room temperature before moving to the next step.

5m



Day 5 or 6

9h 18m

21 Rehydrate slides in H₂O for ⌚ 00:02:00 - ⌚ 00:03:00 at 🌡 Room temperature .

5m



22 Dehydrate slides in sequential EtOH washes (70%, 95%, 100% EtOH) for ⌚ 00:05:00 each wash at 🌡 Room temperature .

5m



23 Clear tissue for ⌚ 00:05:00 in fresh Xylene I then incubate in fresh Xylene for ≥ ⌚ 00:05:00 .

10m



- a. Keep slides in Xylene until coverslipping.
- b. Perform Xylene steps and coverslipping in a chemical fume hood.

24 Use Entellan mounting medium (Electron Microscopy Sciences 14802) to coverslip slides.



24.1 Apply one line of Entellan across the middle of the slide.

24.2 Cover with glass coverslip and gently press down with forceps.

24.3 Let slides dry in the fume hood for ⌚ 01:00:00 – then dry on the bench ⌚ Overnight at 🌡 Room temperature .

1h 5m



24.4 Gently remove excess Entellen with a razor blade.

Optional Cresyl Violet Protocol (NISSL stain)

10m

25 After rehydration and before ethanol steps (dehydration).

26 Incubate slides in fresh cresyl violet for ⌚ 00:10:00 at 🌡 Room temperature .

10m



27 Proceed with dehydration steps normally.

28 Must stain with cresyl violet for TH stereology in substantia nigra.

Secondary Antibodies

29

	A	B	C	D	E
	Target	Species	Conc	Manufacturer	
	Mouse	Goat	1:500	BA-9200-1.5	Biotinylated
	Rabbit	Goat	1:500	BA-1000-1.5	Biotinylated