

Oct 08, 2025

.3% Triton X-100 in PBS

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbw9dngpk/v1>

isabel.dilandro¹

¹Claremont McKenna College

Isabel DiLandro



isabel.dilandro

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbw9dngpk/v1>

Protocol Citation: isabel.dilandro 2025. .3% Triton X-100 in PBS. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbw9dngpk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: October 07, 2025

Last Modified: October 08, 2025

Protocol Integer ID: 229257

Keywords: pb

Abstract

.3% Triton X-100 in PBS

Materials

1L beaker

PBS

Triton X-100

1000 μ L micropipette

Stir plate

Stir bar

1L bottle



- 1 Add 100 mL of PBS to your 1 L beaker and place beaker on stir plate. Add a stir bar.
- 2 Pipette 3000 μ L of Triton X-100 into the beaker and stir until completely dissolved. It may help to break up the Triton X-100 by mixing larger pieces using the micropipette before activating the stir bar.
- 3 Once completely incorporated, add the remaining 900 mL of PBS to the 1 L beaker and stir until clear and no visible Triton X-100 remains.
- 4 Pour the solution into the 1 L bottle and store in the refrigerator.