

Feb 25, 2020

'Uniclear' water-based brain clearing for light sheet imaging of electrode tracks

DOI

<https://dx.doi.org/10.17504/protocols.io.zndf5a6>

Susu Chen¹, Karel Svoboda¹


¹Janelia Research Campus

Optical Clearing of Tissue



Susu Chen

Janelia Research Campus



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.zndf5a6>

Protocol Citation: Susu Chen, Karel Svoboda 2020. 'Uniclear' water-based brain clearing for light sheet imaging of electrode tracks. protocols.io <https://dx.doi.org/10.17504/protocols.io.zndf5a6>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 01, 2019

Last Modified: February 25, 2020

Protocol Integer ID: 21925

Keywords: brain clearing for light sheet imaging, whole mouse brain clearing, based brain clearing, lightsheet microscope, light sheet imaging, cleared brain, electrode track, uniclear procedure, refractive index matching, robust specimen, zeiss z1 lightsheet, brain

Abstract

UniClear procedure for whole mouse brain clearing and refractive index matching. The advantages of this method are that the method is water-based and produces a mechanically robust specimen. Cleared brains can be imaged with lightsheet microscopes (e.g. Zeiss Z1 lightsheet).



Materials

MATERIALS

⊗ DMSO Merck MilliporeSigma (Sigma-Aldrich) Catalog #472301

⊗ OptiPrep™ Density Gradient Medium Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1556)

⊗ 2-Propanol (99.5 %) Merck MilliporeSigma (Sigma-Aldrich) Catalog #278475

2-methyl-2-butanol

<https://www.sigmaaldrich.com/catalog/product/sial/152463?lang=en@ion=US>

2-propanol

<https://www.sigmaaldrich.com/catalog/product/sial/278475?lang=en@ion=US>

Sodium dodecyl sulfate solution (SDS), 10% in H₂O

<https://www.sigmaaldrich.com/catalog/product/sigma/71736?lang=en@ion=US>

OptiPrep™ Density Gradient Medium

<https://www.sigmaaldrich.com/catalog/product/sigma/d1556?lang=en@ion=US>

Dimethyl sulfoxide (DMSO)

<https://www.sigmaaldrich.com/catalog/product/sial/276855?lang=en@ion=US>

Nycodenz

<https://www.progen.com/nycodenz.html>

Preparation of SBiP buffer: 500ml

*** Mix SBiP at 4C till fully dissolved, kept at 4C.

At room temperature, SBiP will get activated and emulsified for delipidation. Use each batch within a month for best effect.

- Ice cold H₂O 350 ml
- 50mM Na₂HPO₄ 2 ml
- 4% SDS (in H₂O, pH7.4) 10 ml
- 2-methyl-2-butanol 80 ml
- 2-propanol 40 ml

Preparation of Opti-prep refractive index matching solution

- 60% (w%) Opti-prep solution



- 20% (w%) DMSO
- 20% (w%) PBS
- 20 mM Tris
- add overall 1.12g/ml Nycodenz (dissolve gradually by adding ~1/10th of Nycodenz powder at a time until fully mixed, heating up and stirring can help to dissolve faster)

References:

Economo, Clark, et al., A platform for brain-wide imaging and reconstruction of individual neurons.


Elife.2016 Jan 20;5:e10566. doi: 10.7554/eLife.10566


Winnubst et al., Reconstruction of 1,000 Projection Neurons Reveals New Cell Types and Organization of Long-Range Connectivity in the Mouse Brain

Cell.2019 Sep 19;179(1):268-281.e13. doi: 10.1016/j.cell.2019.07.042. Epub 2019 Sep 5.



- 1 Perfuse mouse with cold 4% PFA. Post fix brain for at least  24:00:00 at room temperature (RT). Wash brain in PBS to remove fixative. Fixed brain samples can be kept in PBS at  4 °C before further processing.
- 2 Delipidate brain with SBiP buffer:  20 mL per adult mouse brain, shaking (70 RPM) at  37 °C . On the first day: change buffer at  03:00:00 ,  06:00:00 ,  12:00:00 and leave over night. Starting on the second day: change SBiP buffer every  24:00:00 for consecutive 14 days , shaking (70 RPM) at  37 °C .
- 3 After delipidation, wash brain in PBS for at least  12:00:00 before further processing. Sample can be stored in PBS at RT before refractive index matching.
- 4 Refractive index matching: first day in  20 mL 50% PBS + 50% opti-prep mixed solution, rotating at RT for  24:00:00 . On second day transfer the brain to  40 mL 100% opti-prep solution, rotating at RT for  24:00:00 . Index-matched brain should appear transparent.
- 5 For light-sheet imaging (e.g. Zeiss Z1) mount brain vertically.





Example of a mounted brain (fluorescent yellow color is DNA dye YO-PRO, which is optional). For electrode localization we use autofluorescence.

To image tissue autofluorescence use 488 nm laser. To image CM-Dil tracks use 561 nm laser. Brain can be imaged with a voxel size of $1.22\ \mu\text{m} \times 1.22\ \mu\text{m} \times 8\ \mu\text{m}$ on Z1 with a 5x objective.