



Dec 13, 2024

4-Gene Standard Uniplex PCR Protocol for Antibiotic Resistance Genes and Determinants *sul1*, *erm(B)*, *ctx-m₃₂*, and *int1*



DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvme6jbg3p/v1>

Justine C. Condon¹, Lisa M. Durso¹

¹USDA-ARS



Justine C. Condon

USDA-ARS

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvme6jbg3p/v1>

Protocol Citation: Justine C. Condon, Lisa M. Durso 2024. 4-Gene Standard Uniplex PCR Protocol for Antibiotic Resistance Genes and Determinants *sul1*, *erm(B)*, *ctx-m₃₂*, and *int1*. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvme6jbg3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 24, 2024

Last Modified: December 13, 2024

Protocol Integer ID: 108736

Keywords: Polymerase Chain Reaction (PCR), Gel electrophoresis, 4-Gene Standard Uniplex PCR, antibiotic resistance gene, ARG, antibiotic resistance determinant, *su1*, *erm(B)*, *ctx-M-*, *ctx-M-32*, *int1*, sulfonamide resistance, erythromycin resistance, cefotaxime resistance, extended spectrum beta-lactamase resistance, integrase, antibiotic resistance gene, antibiotic resistance determinant, antibiotic resistance, integrase gene *int1*, pcr assay, gene standard uniplex pcr protocol, gel electrophoresis, assayed drug class, general screening tool, assay, determinants *su1*

Funders Acknowledgements:

USDA-ARS

Grant ID: NP 212 Soil and Air

Disclaimer

The use of trade, firm, or corporation names in this publication (or page) is for the information and convenience of the reader. Such use does not constitute an official endorsement or approval by the United States Department of Agriculture or the Agricultural Research Service of any product or service to the exclusion of others that may be suitable.

Abstract

A suite of four antibiotic resistance determinants were selected as a general screening tool for agricultural and environmental samples. It included antibiotic resistance genes from three commonly assayed drug classes (*su1*, *erm(b)* and *ctx-m*) and one determinant often associated with human-impacted samples (the integrase gene *int1*). This protocol describes the correct procedures on how to set-up and run the PCR assay and gel electrophoresis.

Image Attribution

Stock Photo

Guidelines

Scope

- Polymerase Chain Reaction (PCR) is a technique that enables researchers to rapidly and efficiently amplify (copy) a specific region of DNA (Deoxyribonucleic Acid) or RNA (Ribonucleic Acid) into millions of copies.
- A simple PCR reaction consists of target DNA, a set of synthetic oligonucleotide primers that flank the target DNA sequence, a thermostable DNA polymerase (usually Taq polymerase), and nucleotides.
- Thermal cycling exposes these contents to repeated cycles of heating and cooling that promotes DNA replication. There are three steps involved in the process:

1. denaturation
2. annealing
3. elongation

This assay is for the detection of *sul1*, *erm(B)*, *ctx-m-32*, and *int1*. In June 2020, the forward and reverse primers for *int1* were flip-flopped from the Hardwick orientation and are now based according to Dr. Lana Castleberry's recommendations and her use of the primers in the 4G multiplex assay. See 'References' section for references.

Responsibility

This SOP applies to all staff members and students. These individuals must be knowledgeable about the requirements set forth within this document. The lab manager or designee shall ensure that all staff and students know the proper techniques.

Disclaimer

The use of trade, firm, or corporation names in this publication (or page) is for the information and convenience of the reader. Such use does not constitute an official endorsement or approval by the United States Department of Agriculture or the Agricultural Research Service of any product or service to the exclusion of others that may be suitable.



Materials

Reagents and Materials

-  JumpStart™ REDTaq® ReadyMix™ Reaction Mix **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P0982**
- Forward and Reverse Primers – See **Table 2** (Step 4.3)
- PCR Water (HyClone SH30538.01 Sigma, St. Louis, MO)
- Positive Controls – gBlocks– See **Appendix** (Steps 31-34)
-  Agarose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9539**
-  SYBR SAFE DNA stain **Invitrogen - Thermo Fisher Catalog #S33102**
-  DirectLoad™ PCR 100 bp Low Ladder **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D3687**
- 0.2 mL x 8 PCR reaction strips (T320-2N, Simport, Quebec, Canada)
- 2 mL sterile centrifuge tube (02.681.299 Fisher Scientific, Pittsburg, PA)
- Aerosol barrier no retention pipette tips (USA Scientific, Orlando, FL)
- 1X SB buffer (20X Sodium Borate (**SB**) Buffer - SB20-4, FasterBetter Media, Hunt Valley, MD)

Equipment

- PCR Prep Cabinet (3620804 Labconco, Kansas City, MO)
- Pipettors for various volumes from 1-1000µl (Gilson, Middleton, WI)
- Thermal cycler (T100, BioRad, Hercules, CA)
- Mini Centrifuge (C1301P, Labnet International, Big Flats, NY)
- 1 Liter or 500ml sterile Pyrex glass bottle
- Gel casting tray, gel combs, and gel boat
- Electrophoresis machine (PowerEase 300W, Life Technologies, Carlsbad, CA)
- Gel imaging apparatus (Gel Doc XR+, BioRad, Hercules, CA)
- Lighter

Safety warnings

-  Wear proper personal protective equipment (PPE) when working through PCR: nitrile gloves and lab coat.

When working with the hot agarose, wear thermal hot-pad gloves.

Procedure - PCR

15m

- 1 If frozen, thaw all reagents and samples and keep on an ice block or ice for stability. If kept in the refrigerator, keep there until ready for use, and keep on an ice block or ice for stability when in use. 
- 2 Label tubes for the samples and controls.
- 3 Prepare the PCR prep area for setting up a PCR reaction to eliminate contamination – UV 0.2mL 8-strip PCR tubes, PCR water, and microcentrifuge tube for  00:15:00 . 15m 
- 4 Combine the Master Mix components for the number of reactions (allow for 10% extra Master Mix) using the recipe from **Table 1** ( [go to step #4.2](#)) with the correct primer sets (**Table 2** ( [go to step #4.3](#))) for 20µL reactions.
- 4.1 Before adding to the Master Mix - mix the primers by aspirating and expelling with the pipette tip.  

Note

***Beware, overmixing can degrade the DNA.**

4.2

A	B	C	D
Component	Description	Volume/RX N (µl)	Final Concentration
Primer 1	Forward Primer (100µM)	0.4	2µM
Primer 2	Reverse Primer (100µM)	0.4	2µM
TAQ	JumpStart RedTaq ReadyMix Reaction Mix (Cat.#P0982)	10	1X

	A	B	C	D
	PCR Water	PCR Water (HyClone SH30538.01)	4.2	N/A
	Sample	Template (& positive control)	5.0	

Table 1. PCR Master Mix recipe. 20µL total reaction volume.

4.3



	A	B	C	D	E	F	G
	Gene	Type	Sequence 5' to 3'	T _M (°C)	TC Conditions	Amplicon Size (bp)	Sequence Reference
	<i>sul1</i>	Forward	GACGAGATT GTGCGGTTC TT	64	1 cycle at 94°C for 2 min; 35 cycles at 94°C for 30s, 64°C for 30 s, 72°C for 2 min; 1 cycle at 72°C for 5 min.	185	Szczepanowski et al., 2009
		Reverse	GAGACCAAT AGCGGAAGC C				
	<i>erm(B)</i>	Forward	GATACCGTT TACGAAATT GG	58	1 cycle at 94°C for 2 min, 35 cycles at 94°C for 30 s, 58°C for 30 s, 72°C for 2 min; 1 cycle at 72°C for 5 min.	364	Chen et al, 2007
		Reverse	GAATCGAGA CTTGAGTGT GC				
	<i>ctx-m-32</i>	Forward	CGTCACGCT GTTGTTAGG AA	63	1 cycle at 94°C for 2 min; 35 cycles at 94°C for 30 s, 63°C for 30 s, 72°C for 2 min; 1 cycle at 72°C for 5 min.	156	Szczepanowski et al., 2009
		Reverse	CGCTCATCA GCACGATAA				

	A	B	C	D	E	F	G
		e	AG				
	<i>intl1</i>	Forward	ACATGCGTG TAAATCATC GTCG	60	1 cycle at 94°C for 2 min; 35 cycles at 94°C for 30 s, 60°C for 30 s, 72°C for 2 min; 1 cycle at 72°C for 5 min	473	Hardwick et al., 2008
		Reverse	CTGGATTTC GATCACGGC ACG				Castleberry 2018

Table 2. Primer sequences, melting temperature, and thermal cycling conditions.

5 Vortex the completed Master Mix for at least  00:00:20 .

20s



6

Note

Check the pipette's volume as you go along, making sure it is exactly at the correct volume

7 Carefully pipette  5 µL of PCR Water to the non-template control (NTC) tube.

8 Add  15 µL of Master Mix to each labeled, 0.2mL PCR tubes.



9 Next, gently vortex to mix, then pipette  5 µL of sample ( 5 µL of working stock; or  1 µL straight DNA +  4 µL PCR water) into its same labeled tube.



Note

Use a new pipette tip for each sample.



9.1 To make a working stock dilution, see **SOP: IDT Primer and gBlock Hydration and Aliquots**.

10 Add  5 μL of positive control (mix first by aspirating and expelling with the pipette tip) to its same labeled tube.



11 Centrifuge the strip tubes so that the entire volume is at the bottom.



Note

Make sure the centrifuge is balanced with both sides containing tubes.

12 Load the tubes into the Thermocycler, keeping them in order.

13 Find the programmed gene conditions listed under '4G --' and double check that they match the conditions listed in **Table 2** [ [go to step #4.3](#)].

14 Run PCR



▪ When cycle is complete:

14.1 **Procedure - Gel Electrophoresis**, Steps  [go to step #15](#) -  [go to step #30.3](#) **OR**

14.2 Place the tubes in a labeled rack at  -20 °C until ready to run a gel.

14.3 **Turn off the Thermocycler.**

Procedure – Gel Electrophoresis

3h 11m 30s

15 Using a spoon, measure out  2.6 g of powdered agarose into a weigh boat.



16 Measure  260 mL of 1X **SB** buffer into a graduated cylinder.

17 Add the powdered agarose and 1X **SB** into 500mL bottle and swirl to mix.



18 Microwave the **SB** and agarose mixture for  00:01:00 with bottle lid LOOSE.

1m

19 Remove bottle from microwave (**using hot pad gloves**), tighten the lid and gently swirl.



Safety information

For Steps 19 - 27 - Use thermal hot pad gloves to protect from burns.

20 LOOSEN the lid and put in microwave for another minute.

21 Check mixture visually, if flecks are present repeat  [go to step #19](#) , loosen lid and microwave again for  00:00:30 .

30s

22 Add  26 μ L of SYBR safe to molten gel.



23 Swirl gently with a loose lid to distribute SYBR.

24 Place bottle in  55 °C water bath for  00:10:00 to temper.

10m



25 While the mixture is tempering, set up gel cast by putting two plastic stoppers into the ends of the cast, with the orange rubber side of the stopper facing out.

26 Put two (or three) 36-well combs into the first and fourth notches of the gel cast (gel cast is right side up if first notch is close to the top of the cast)

27 After the gel has tempered, slowly pour it (**using hot pad gloves**) into the gel cast.





28 If there are any noticeable bubbles, take a lighter and hold it over the bubble, which should make it pop.

29 Cover gel loosely with foil and allow it to set (up to  01:00:00) before use.

1h

Note

****SB gels will be flimsier, so be careful when removing the combs or transporting the gel****

30 Electrophoresis Conditions

30.1  10 μ L of PCR amplicon will go into each well flanked by  5 μ L of ladder.

30.2 Run at 120 volts (V) for  02:00:00 in 1x SB buffer.

2h

Note

Make sure **SB gel** goes into a gel boat with **SB buffer**.



30.3 Remove processed gel and analyze using a gel imaging system.



Appendix – Positive Control (gBlock) Sequences

31

Forward, **Reverse**

*Note: 'Reverse' denotes the 'Reverse Complement'

gBlock 3 – *su1* (344 bp) - for 4G assay and primer set P2-sul1



Escherichia coli R46 sul1 gene for sulfonamide-resistant dihydropteroate synthase Sul1, complete CDS NCBI Reference Sequence: NG_048081.1

```
CCTTCGACCCGAAGACGCGCTCGACGAGATTGTGCGGTTCTCGAGGCGCGGGTTCCGCCTTGCGACGGAGCGG
GGTCGCTGCCGACCGGCTCATCTCGATCCGGGGATGGGATTTTCTTGAGCCCGCACCAGAAACATCGCTGCAC
GTGCTGTGAACCTTCAAAAGCTGAAGTCGGCGTTGGGGCTTCCGCTATTGGTCTCGGTGTGCGGAAATCCTTCT
TGGGCGCCACCGTTGGCCTTCTGTAAAGGATCTGGGTCCAGCGAGCCTTGCGGCGGAACTTACGCGATCGGCA
ATGGCGCTGACTACGTCCGCACCCACGCGCTGGAGATCTGC
```

32 **gBlock 17 – *ermB* (424 bp)** – for 4G assay and primer set P3-*ermB*

```
ATAATAAAACAATTGAATTTAAAGAAACGATACCGTTTACGAAATTGGAACAGGTAAAGGGCATTAAACGACGAACTGGC
TAAATAAGTAAACAGGTAAAGTCTATTGAATTAGACAGTCATCTATTCAACTATCGTCAGAAAAATTAACCTGAATACTCG
TGTCACCTTAATTCACCAAGATATTCTACAGTTCAATCCCTAACAAACAGAGGTATAAAATTGTTGGGAATATTCCTTACAAT
TTAAGCACACAAATTATTAAGGAGTGGTTTTGAAAGCCGTGCGTCTGACATCTATCTGACTGTTGAAGAAGGATTCTACAAG
CGTACCTTGGATATTCACCGAACACTAGGGTTGCTCTTGCACACTCAAGTCTCGATTCAGCAATTGCTTAAGCTGCCAGCGGAA
TGCT
```

33 **gBlock 18 – *ctxm-32* (249 bp)** – for 4G assay and primer set P4-*ctxm32*

```
CGGTCAGTTCACGCTGATGGCGACGGCAACCGTCACGCTGTTGTTAGGAAGTGTGCCGCTGTATGCGCAAACGGC
GGACGTACAGCAAAAACCTTGCCGAATTAGAGCGGCAGTCGGGAGGAAGACTGGGTGTGGCATTGATTAACACAG
CAGATAATTCGCAATACTTTATCGTGTGCTGATGAGCGCTTTGCGATGTGCAGCACCAGTAAAGTGATGGCCGTGGC
CGCGGTGCTGAAGAAAAGTGAAAG
```

34 **gBlock 19 – *int1* (535 bp)** – for 4G assay and primer set P5-*int1*

```
GTGAGCGCACTCCGGCACCGCCAACCTTTCAGCACATGCGTGTAATCATCGTCGTAGAGACGTCGGAATGGCCGA
GCAGATCCTGCACGGTTCGAATGTCGTAACCGCTGCGGAGCAAGGCCGTGCGAACGAGTGCGGAGGGGTGTGC
GGTGTGGCGGGCTTCGTGATGCCTGCTTGTCTACGGCACGTTTGAAGGCGCGCTGAAAGGTCTGGTCATACATG
TGATGGCGACGCACGACACCGCTCCGTGGATCGGTGCAATGCGTGTGCTGCGCAAAAACCCAGAACACGGCCA
GGAATGCCCGGCGCGCGGATACTCCGCTCAAGGGCGTCGGGAAGCGCAACGCCGCTGCGGCCCTCGGCCTGGT
CCTTCAGCCACCATGCCGTGCACGCGACAGCTGCTCGCGCAGGCTGGGTGCCAAGCTCTCGGGTAACATCAAGG
CCCATCCTTGAGCCCTTGCCCTCCGCGACGATGATCGTGCCGTGATCGAAATCCAGATCCTTGACCCGCAGTTG
CAAACCTCACT
```



Protocol references

Castleberry, Lana. 2018. 4Gene Multiplex PCR assay troubleshooting notes.

Chen et al, 2007. Chen, J., Y. Zhongtang, F. C. Michel, Jr., T. Wittum, and M. Morrison. 2007. Development and Application of Real-Time PCR Assays for Quantification of erm Genes Conferring Resistance in Macrolides-Lincosamides-Streptogramin B in Livestock Manure and Manure Management Systems. Appl and Environ Micro. 73(14):4407-4416.

Hardwick et al., 2008 Hardwick, S.A., H.W. Stokes, S. Findlay, M. Taylor, and M.R. Gillings. 2008. Quantification of class 1 integron abundance in natural environments using real-time quantitative PCR. FEMS Microbiol. Lett. 278:207-212.

Szczepanowski et al., 2009. Szczepanowski, R., B. Linke, I. Krahn, K-H. Gartemann, T. Guetzkow, W. Eicher, A. Pühler, A. Schlueter. 2009. Detection of 140clinically relevant antibiotic-resistance genes in the plasmid metagenome of wastewater treatment plant bacteria showing reduced susceptibility to selected antibiotics. Microbiology. 155(7):2306-2319.

Thames et al., 2012 Thames, C.H., Pruden, A., James, R.E., Ray, P.P., Knowlton, K.F. Excretion of antibiotic resistance genes by dairy calves fed milk replacers with varying doses of antibiotics. Frontiers in Microbiology. (2012) 3:139 1-12.