



Nov 13, 2025

4-HNE immunofluorescence staining

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5k56jv1b/v1>

lene.devos¹, tine.wylin¹

¹Abdominal Transplantation, KU Leuven

Abdominal Transplantati...



Lene Devos

KU Leuven

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5k56jv1b/v1>

Protocol Citation: lene.devos , tine.wylin 2025. 4-HNE immunofluorescence staining . protocols.io
<https://dx.doi.org/10.17504/protocols.io.8epv5k56jv1b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working in rat liver and kidney tissue as well as pig liver and lung tissue.

Created: October 28, 2025

Last Modified: November 13, 2025

Protocol Integer ID: 231005

Keywords: 4-hydroxynonenal, lipid peroxidation, immunofluorescence, liver, kidney, lung, rat, pig, hne immunofluorescence, hydroxynonenal, immunofluorescence, hne, using ab46545

Abstract

This protocol provides an overview of the different steps involved in immunofluorescence staining of 4-hydroxynonenal (4-HNE) using ab46545. This protocol has proven to be effective in rat kidney and liver tissue as well as pig liver and lung tissue.

Materials

- Paraffin tissue section of 5 μ m
- Coplin jars (a volume of 250 ml is enough to immerse slides)
- Easy Dip™ container + staining rack
- Xylene
- Ethanol (50% - 100%)
- MiliQ (MQ)
- Citric acid (0,1 M):

	A	B	C	D
	Reagent (concentration)	Volume needed		
	Citric acid (anhydrous powder) (192,123 g/mol)	9,60 g	4,80 g	1,92 g
	MQ	500 ml	250 ml	100 ml

- Sodium citrate (0,1 M, pH 6):

	A	B	C	D
	Reagent (concentration)	Volume needed		
	Sodium citrate (powder) (258,07 g/mol)	12,91 g	6,46 g	2,58 g
	MQ	500 ml	250 ml	100 ml

Note: Start with a lower volume and adjust to pH. Afterwards, you can fill until the desired volume. Check final pH to be sure.

- Citric acid buffer:

	A	B
	Reagent (concentration)	Volume needed
	Citric acid (0,1M)	2,30 ml
	Sodium citrate (0,1M, pH 6)	10,70 ml



A	B
MQ	Add until total volume of 130 ml

- Microwave
- Slide container with wet tissues
- Tris NaCl buffer:

A	B
Reagent (concentration)	Volume needed
Tris-HCl (0,1 M, pH 7,5) (157,6 g/mol)	7,88 g
NaCl (0,15 M) (58,44 g/mol)	4,38 g
MQ	500 ml

- TNT (Tris NaCl, Tween20) buffer:

A	B
Reagent (concentration)	Volume needed
Tris NaCl buffer	50 ml
Tween20 (0,05%)	25 µl

- TNB (Tris NaCl Blocking) buffer:

A	B
Reagent (concentration)	Volume needed
Tris NaCl buffer	50 ml
Bovine Serum Albumin (0,5%)	0,25 g

- Goat serum
- Primary antibody: Rabbit Polyclonal 4-HNE antibody (1/500; Cat.# ab46545; RRID: AB_722490)



- Secondary antibody: Goat anti-Rabbit IgG (H+L) Cross-Adsorbed Secondary Antibody, Alexa Fluor™ 488 (4 µg/ml; Cat.# A-11008; RRID: AB_143165)
- DAPI (Cat.#62248)
- PBS (1x)
- Coverslips
- Slow Fade Diamond antifade (Fisher Scientific (15,441,244))
- Nail polish

Troubleshooting

Problem

Insufficient antigen retrieval

Solution

Microwave settings might differ between brands/types. Test at which wattage and time combination the antigen retrieval buffer reaches boiling point, without buffer overflowing. Keeping the lid open can help to prevent boiling over. Use this setting to preheat the buffer. Next, find a setting that helps maintain this temperature without the buffer boiling over. This setting can be used for the slides.

Before start

Prepare all buffers and reagents before start of deparaffinization



Deparaffinization

25m

1 Submerge the slides in the following reagents:

- 1.1 5' Xylene
5' Xylene
3' 100% Ethanol
3' 90% Ethanol
3' 70% Ethanol
3' 50% Ethanol
3' MQ

25m

Antigen retrieval

42m 30s

2 Fill the Easy DipTM container with citric acid buffer (±120-130 ml)

3 Preheat the citric acid buffer in the microwave for 1,5' at 750W (keep the lid of the container open).

1m 30s

4 Immerse the slides in the preheated buffer

5 Place the slides in the microwave for 2 cycles of 4,5' at 90W with 1-2 minutes in between (keep the lid of the rack elevated).

11m

6 Let the slides cool down for 30' in the buffer

30m

Incubation with primary antibody

1d 1h 5m

7 Prepare a chamber (slide container) with wet tissues on the bottom to maintain humidity

8 Wash 5' with TNT

5m

9 Dry slides carefully with tissue paper without affection the sections



10 Block for 1h with TNB + 1/50 normal goat serum

1h

11 Incubate the slides with 1st antibody (TNB + 1/50 normal goat serum + 1/500 4HNE primary antibody) overnight at 4°C

1d

Incubation with secondary antibody

1h 15m

12 Wash 3× 5' with TNT

15m

13 Dry slides carefully with tissue paper

14 Cover the chamber from light

15 Incubate the slides for 1h with 2nd antibody (TNB + 1/50 normal goat serum + 4 µg/ml 2nd antibody) at room temperature

1h

DAPI nuclear staining

35m

16 Wash 3× 5' with TNT

15m

17 Dry slides with tissue paper

18 Add DAPI (1/2000 in PBS) for 5'

5m

19 Dry slides with tissue paper

20 Add DAPI again (1/2000 in PBS) for 5'

5m



21 Wash 5' PBS at RT in glass jar or Easy Dip™ container

5m

22 Repeat wash 5' PBS at RT in glass jar or Easy Dip™ container

5m

Mounting

23 Dry slides with tissue paper

24 Add 1 drop of slow fade in the middle

25 Place a coverslip on top of the slide and remove air bubbles

26 Apply nail polish on the borders to secure the coverslip

27 Store overnight in the dark at room temperature

28 After drying/visualising, the slides can be stored in the fridge or at -20°C

Acknowledgements

We would like to thank Sai Kocherlakota, PhD for sharing his expertise in this staining.