

Apr 28, 2023

5' RACE for RNA fragments with 5' phosphate

DOI

dx.doi.org/10.17504/protocols.io.261ge3mzwi47/v1

Marta Gaglia¹

¹UW Madison



Marta Gaglia

UW Madison

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge3mzwi47/v1>

Protocol Citation: Marta Gaglia 2023. 5' RACE for RNA fragments with 5' phosphate. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ge3mzwi47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 28, 2023

Last Modified: April 28, 2023

Protocol Integer ID: 81154

Keywords: rna fragment, race for rna fragment, ligation of an rna adapter, rna adapter, rna, viral rnase, mammalian transcriptional start site, identification of fragment, fragment, phosphate, fragment end, gene, ligation, expected fragment end

Funders Acknowledgements:

NIH

Grant ID: R01AI137358

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol will allow identification of fragments of RNA that have a 5' phosphate. It relies on ligation of an RNA adapter to the 5' end of the fragment and amplification with primers that bind the adapter and gene-specific sequences ~200 nt downstream of predicted/expected fragment end. It is adapted from protocols to map mammalian transcriptional start sites, but lacks any reaction to remove the cap. We have used it to isolate fragments generated by viral RNases.

Materials

Reagents

Ambion Turbo DNase (Thermofisher – Cat#AM2239)

Recombinant RNasin RNase inhibitor (Promega – cat# N2515)

Ultrapure DNase/RNase free dH₂O (Thermofisher Cat# 10-977-015)

Phenol/chloroform with acidic pH (we use Fisher cat#BP1752I-100 without adding the additional buffer, but other should work)

3 M sodium acetate (Thermofisher – Cat# AM9740)

Glycogen (Thermofisher – Cat# AM9510)

100% Ethanol (we use VWR Cat# 89125-170, but others should be fine)

T4 RNA ligase (Ambion Cat#AM2141)

Invitrogen M-MLV reverse transcriptase (Thermofisher – Cat# 28-025-013)

Random decamers (Thermofisher – Cat#AM5722G)

dNTPs (we use Thermofisher – Cat#FERR0181, but others should work; keep RNase free) – *make a working stock of 25 mM*

Taq DNA polymerase with standard Taq buffer (New England Biolabs – Cat#M0273)

RACE RNA adapter – *make working stocks of 0.3 µg/µl*

GCUGAUGGCGAUGAAUGAACACUGCGUUUGCUGGCCUUUGAUGAAA

Primers – *make working stocks of 10 µM:*

5' RACE outer forward primer (*binds RNA adapter*):

GCTGATGGCGATGAATGAACACTG

5' RACE inner forward primer (*binds RNA adapter sequences and allows for Gibson cloning of fragments into the EcoRI site of pCDNA3.1(+)*):

GGATCCACTAGTCCAGTGTGGTGAATTCTGAACACTGCGTTTGCTGGCTTTGATG

Gene-specific outer and inner reverse primers – *design ~200 nt downstream of expected cleavage site; should allow for nested PCR.*

Equipment

Thermocycler

Water bath or heat block set at 37°C

Refrigerated microcentrifuge

Troubleshooting

1 **RNA purification:**

Extract RNA and elute/resuspend in 50 µl H₂O

2 **DNase treatment:**

To each 50 µl sample, add 10 µl of DNase reaction mix containing:

- 6 µl DNase buffer
- 1.2 µl Turbo DNase
- 1.5 µl RNasin
- 1.3 µl H₂O

(total volume = 60 µl)

Incubate for 15 min at 37°C (waterbath or heatblock).

3 **RNA extraction with phenol/chloroform:**

Add 140 µl H₂O to each tube, then 200 µl phenol/chloroform.

Vortex for 15 sec.

Spin down for 5 min at 12,000 xg at room temperature.

Transfer the top aqueous phase to a new tube with 2 µl glycogen + 20 µl 3M sodium acetate.

Add 600 µl 100% ethanol to each tube, invert a few times to mix.

Incubate at -20°C for 1h (can leave it overnight if needed).

Spin down to pellet at 16,000xg at 4°C for 20 min.

Remove supernatant and wash pellet with 1ml freshly made 70% ethanol.

Spin down again at 16,000xg at 4°C for 5 min.

Remove supernatant and spin down again for 1min.

Remove all liquid and dry pellet for a few min (not too long).

Add 12 µl water and leave at RT for 10 min to resuspend.

Measure concentrations.

At this point, can set aside 250-500 ng of DNase-treated RNA for RT-qPCR to test RNA levels and/or RNase activity.

4 **RACE adaptor ligation:**

Depending on concentrations of RNA samples, try to use as much RNA as possible for the ligation (i.e. use 10 µl of the lowest sample if one sample is very low). Between 1.5 and 6 µg of RNA works well. Dilute all samples to this amount in 10 µl H₂O in PCR tubes.

To each sample add 5 µl of ligation mix containing:

- 1.5 µl 10X RNA buffer
- 1 µl T4 RNA ligase (5 U/µl)
- 1.5 µl 5' RACE adapter (0.3 µg/µl)



- 1 µl dH₂O

(total volume= 15 µl)

Incubate samples at 25°C for 2 h in a thermocycler.

5 Reverse transcription with M-MLV RT

Mix 2 µl ligated RNA with 8 µl of primer-dNTP mix containing :

- 0.4 µl dNTP (25 mM)
- 1 µl random decamers
- 6.5 µl dH₂O

Heat mixture to 65°C for 5 min in thermocycler, then chill quickly on ice and incubate on ice for at least 1 min.

Add to the reaction 7 µl of buffer-DTT mix containing :

- 4 µl RT buffer
- 2 µl DTT (0.1 M, included in kit)
- 1 µl RNasin

Mix and incubate at 37°C for 2 min in thermocycler.

Add 1 µl of M-MLV RT to each reaction and mix.

Carry out RT reaction in thermocycler:

- 25°C for 10 min
- 37°C for 50 min
- 70°C for 15 min

6 Outer PCR with Taq polymerase:

Mix 1 µl ligated cDNA with 49 µl of PCR mix containing :

- 0.4 µl dNTP (25 mM)
- 5 µl 10x Standard Taq buffer
- 0.25 µl Taq (5 U/ µl)
- 2 µl 5' RACE outer primer (10 µM)
- 2 µl Gene-specific outer primer (10 µM)
- 39.35 µl dH₂O

Run the following program on a thermocycler:

- 94°C for 3min
- 35 cycles of :
 - 94°C for 30 sec

60°C for 30 sec
72°C for 30 sec
- 72°C for 7 min

7 **Inner PCR with Taq polymerase:**

Mix 2 µl of outer PCR reaction with 48 µl of PCR mix containing :

- 0.4 µl dNTP (25 mM)
- 5 µl 10x Standard Taq buffer
- 0.25 µl Taq (5 U/ µl)
- 2 µl 5' RACE inner primer v2(10 µM)
- 2 µl Gene-specific inner primer (10 µM)
- 38.35 µl dH₂O

Run the same PCR program as outer PCR.

8 **Fragment analysis on agarose gel**

Run inner (and if desired outer) PCR on agarose gel

We usually use a 2% agarose gel, 15 × 20 cm, and run it for 3 h at 70V

Fragment size depends on location of primers relative to cut site (remember to add adaptor/primer sequence when calculating it, since it is pretty large)