

May 29, 2024

Version 2

709.1 Staining of Dissociated Lung Cells for Cellular Indexing of Transcriptomes and Epitopes by Sequencing (CITE sequencing) V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v2>

Gautam Bandyopadhyay¹, Jeffrey Malik², Heidie Huyck³, Ravi Misra⁴, Gloria S Pryhuber⁵, John Ashton⁴

¹Department of Pediatrics, University of Rochester Medical Center;

²Genomic Research Center, University of Rochester Medical Center; ³University of Rochester; ⁴URMC;

⁵University of Rochester Medical Center

LungMap2 Consortium

URMC Pryhuber Lab



Gautam Bandyopadhyay

Department of Pediatrics, University of Rochester Medical Ce...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v2>



Protocol Citation: Gautam Bandyopadhyay, Jeffrey Malik, Heidie Huyck, Ravi Misra, Gloria S Pryhuber, John Ashton 2024. 709.1 Staining of Dissociated Lung Cells for Cellular Indexing of Transcriptomes and Epitopes by Sequencing (CITE sequencing). protocols.io <https://dx.doi.org/10.17504/protocols.io.e6nvw1jk7lmk/v2> Version created by **Gloria S Pryhuber**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 29, 2024

Last Modified: May 29, 2024

Protocol Integer ID: 100748

Keywords: dissociated lung cells for cellular indexing, staining frozen lung cell, frozen lung cells for cite, dissociated lung cell, epitopes by sequencing, sequencing, transcriptome, cellular indexing, cell

Funders Acknowledgements:

NHLBI Molecular Atlas of Lung Development Program Human Tissue Core Grant

Grant ID: U01HL14886

Abstract

This article describes step by step protocol for antibody-staining frozen lung cells for CITE sequencing.

Thawing of frozen dissociated human lung cells

- 1 Frozen cells were thawed in 37°C water bath and soon after the freezing medium melted cells were transferred into a 15 ml conical tube and the tubes were filled with staining buffer [1% BSA (w/v, Millipore) in DPBS (Biowhittaker)] to dilute the DMSO in freezing medium.
- 2 Cells were centrifuged at 800g, 4°C for 10 minutes.
- 3 After centrifugation supernatants were discarded and cells were resuspended in 1.0 ml staining buffer.

Fc receptor blocking and antibody staining

- 4 Viable cell counts were determined by trypan blue exclusion assay
- 5 Cells were then incubated in staining buffer containing 2% v/v normal mouse serum (Millipore Sigma) for 10 minutes over ice for Fc-receptor blocking.
- 6 Cells were then washed with staining buffer and incubated at a final concentration of 1×10^6 cells per 100 μ l of staining cocktail containing analyte-specific antibodies (**attached table**) for 30 min at 4°C.  Covid Supp CITE Seq Staining Anti... 14.5KB
- 7 After antibody-staining, cells were washed three times (800g, 4°C for 10 minutes) and resuspended in 100 μ l staining buffer and passed through a 100- μ m strainer (Falcon, Corning) and sent Genomic Research Center, University of Rochester Medical Center for further processing for CITE sequencing.