



Apr 26, 2016

## Algal culture harvest and RNA extraction for RNA-Seq

DOI

[dx.doi.org/10.17504/protocols.io.ewjbfcn](https://dx.doi.org/10.17504/protocols.io.ewjbfcn)

Alle Lie

Protist Research to Opti...

Caron Lab - Protistan Ec...



Alle Lie

University of Southern California

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewjbfcn>

**Protocol Citation:** Alle Lie 2016. Algal culture harvest and RNA extraction for RNA-Seq. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.ewjbfcn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** April 26, 2016

**Last Modified:** February 21, 2018

**Protocol Integer ID:** 2731

**Keywords:** rna extraction for rna, rna extraction, algal culture harvest, rna, seq, extraction

## Troubleshooting



- 1 Use sterile techniques to harvest 800 ml of culture, and divide into 4 250 ml centrifuge bottles with 200 ml each.
- 2 Spin the centrifuge bottles at 5000 g for 10 min.
- 3 Carefully decant the supernatant. Use 0.5 ml of RNALater to resuspend the pellet and transfer it to a microcentrifuge tube. Use another 0.25 ml of RNALater to wash the bottle and obtain any remaining cells. Transfer this 0.25 ml into the same microcentrifuge tube.
- 4 Depending on whether you want pseudo-replicates, you can keep the 4 harvested pellet separate or mixed.  
Store the harvested cells at or below -20oC until ready for RNA extraction
- 5 Before RNA extraction, the RNALater needs to be removed.  
Thaw samples if they are frozen. Centrifuge samples at > 10 k g for 10 min at 4oC.  
RNALater sometimes affects the buoyancy of the cells, and make them a bit harder to pellet. If you find that the cells are not pelleting, spin harder for a longer time.  
Decant supernatent with RNALater
- 6 Extract RNA using a RNeasy Plant Mini Kit (Qiagen) following the manufacturer's instructions
- 7 Remove any DNA from the extracted samples using DNase (Sigma), following manufacturer's instructions
- 8 Clean up the total RNA using a RNA Clean & Concentrator kit (Zymo), following the manufacturer's instructions  
Quantify the RNA using any preferred method, e.g. Qubit fluorometer
- 9 Store the RNA at -80oC  
It may be ideal to aliquot some RNA out so that you don't need to thaw the entire RNA sample if you need some for any reason