

Jul 29, 2015

Can Crusher-Culley Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.c3jykm>

Alexander Culley¹, Grieg Steward¹

¹Matthew Sullivan Lab, University of Arizona/Ohio State University

VERVE Net

Sullivan Lab



Bonnie Poulos

Matthew Sullivan Lab, University of Arizona, The Ohio State ...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.c3jykm>

Protocol Citation: Alexander Culley, Grieg Steward 2015. Can Crusher-Culley Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.c3jykm>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working



Created: June 18, 2015

Last Modified: November 10, 2017

Protocol Integer ID: 843

Keywords: rna virus, new genera of rna virus, rna, culley protocol culley protocol, environmental microbiology, virus, polymerase gene sequence, subtropical seawater, culley, seawater

Abstract

Culley Protocol (from Culley & Steward (2007), **New genera of RNA viruses in subtropical seawater, inferred from polymerase gene sequences**, Applied & Environmental Microbiology 73(18):5937-5944).

Guidelines

Supplies needed (per 100 mL sample):

Sample (0.2 µm filtered seawater)One 60 mL syringe
One Anotop (0.02 µm filter)
One piece of parafilm
One piece of aluminum foilOne Ziploc freezer bag
Sharpie (ultra-fine)
Can crusher



- 1 Using sterile technique, remove the plunger from the 60mL syringe.
 - 2 Attach the 0.02 μm filter (Anotop) to the end of the syringe.
 - 3 Pour 50 mL of 0.2 μm filtrate into the 60-mL syringe.
 - 4 Replace plunger and mount the syringe and filter into the can crusher.
 - 5 Push the 0.2 μm filtrate **very slowly** through the Anotop, using the can crusher and applying **gentle, steady pressure**. Watch carefully for the last drop of water to go into the Anotop, and then **stop pushing**.
- Note**
- a) Do not let any air push through the Anotop—that may cause it to burst.
 - b) Do not suck up any water into the outlet of the Anotop.
- 6 Using sterile technique, remove the Anotop from the syringe. Do not touch the filter inlet or outlet.
 - 7 Suck up the remaining 50mL of 0.2 μm filtrate for this depth into the 60-mL syringe.
 - 8 Re-attach the Anotop to the end of the syringe.
 - 9 Repeat Step 5: Push the 0.2 μm filtrate **slowly** through the Anotop, using the can crusher.
- Note**
- Be careful not to push air through the Anotop or suck any water back up into it.
- 10 Remove the Anotop from the syringe.



- 11 Write on the Anotop itself using a waterproof Sharpie: date, depth, location (or a shorthand for location if key is explained elsewhere)
- 12 Wrap the Anotop in clean parafilm, covering both the inlet and outlet.
- 13 Label a piece of aluminum foil with date, depth, and location (using waterproof Sharpie).
- 14 Wrap the parafilmed Anotop in the labeled aluminum foil.
- 15 Label a Ziploc bag with date, depths, and location.
- 16 Collect several Anotops to put into one labeled Ziploc: one location/sample site per bag.
- 17 Put the foil-wrapped Anotops into the labeled Ziploc bag.
- 18 Freeze the Ziploc in -80°C freezer (if no -80°C freezer is available, use a -20°C freezer).
- 19 Discard the syringe and filtrate (unless you need some sterile water).