

Jun 13, 2016

## CviRI Assay Conditions

DOI

<https://dx.doi.org/10.17504/protocols.io.es3begn>

DAVID DUNIGAN AND IRINA AGARKOVA

VERVE Net



Irina Agarkova

UNL

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.es3begn>

**Protocol Citation:** DAVID DUNIGAN AND IRINA AGARKOVA 2016. CviRI Assay Conditions. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.es3begn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** March 29, 2016



**Last Modified:** January 20, 2018

**Protocol Integer ID:** 2619

**Keywords:** cviri purification, nc64a chlorella, cviri assay condition, cviri assay conditions for use, chlorella, purification, nc64a, virus

## Abstract

For use in **CviRI Purification From XZ-6E Virus Infected NC64A Chlorella.**



- 1 All assays are carried out in 20.0  $\mu$ L volumes with 1  $\mu$ g of pUC19 DNA as substrate for 60 to 120 min at 25°C.

 02:00:00

- 2 The assays are electrophoresed on 1.2% agarose gels (100 mL gels with double 20 lane combs) for 1 hour at 100 volts in 1X TPE buffer.

 01:00:00

- 3 Gels are stained with 0.5  $\mu$ g/mL ethidium bromide for 30 min and visualized and photographed on a UV light box.

 00:30:00