

Dec 04, 2015

## Embedding yeast colonies for light and electron microscopy

DOI

[dx.doi.org/10.17504/protocols.io.dw97h5](https://dx.doi.org/10.17504/protocols.io.dw97h5)

Sarah Piccirillo<sup>1</sup>, Melissa G. White<sup>1</sup>, Jeffrey C. Murphy<sup>1</sup>, Douglas J. Law<sup>1</sup>, and Saul M. Honigberg<sup>1</sup>

<sup>1</sup>protocols.io

Genetics

Yeast Protocols, Tools, a...



**Saul Honigberg**

University Missouri - Kansas City

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.dw97h5>

**External link:** <http://www.genetics.org/content/184/3/707>

**Protocol Citation:** Sarah Piccirillo, Melissa G. White, Jeffrey C. Murphy, Douglas J. Law, and Saul M. Honigberg 2015. Embedding yeast colonies for light and electron microscopy. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dw97h5>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**Created:** September 30, 2015

**Last Modified:** March 10, 2018

**Protocol Integer ID:** 1729

**Keywords:** embedding yeast colony, yeast colonies for light, alkaline ph regulate pattern formation in yeast, alkaline ph regulate pattern formation, yeast, electron microscopy, pacc pathway

## Abstract

The following is a modification of a Scherz R., 2001 method for colony embedding. The protocol described here is from:

Sarah Piccirillo, *et. al.* **The Rim101p/PacC Pathway and Alkaline pH Regulate Pattern Formation in Yeast Colonies**(2010)

*Genetics*184:707–716; doi:10.1534/genetics.109.113480

Please see the **full manuscript** for additional details.

## Guidelines

To embed colonies for sectioning, we modified a previously described method (**Scherz et al. 2001**), with a major change being the substitution of Spurr's reagent as embedding medium.

### Original method:

SCHERZ, R., V. SHINDER and D. ENGELBERG, 2001**Anatomical analysis of *Saccharomyces cerevisiae* stalk-like structures reveals spatial organization and cell specialization**. *J. Bacteriol.*183: 5402–5413.






## Troubleshooting











- 1 Incubate approximately 300 colonies on agar medium for the indicated time.
- 2 Remove an isolated colony (1-2 mm in diameter) and a small amount of the underlying agar medium, and place on a microscope slide
- 3 Place several drops of 2% agar (42°C) on a microscope slide, and immediately place the colony on the agar and then place several drops of agar on top of the colony and allow to solidify.
- 4 Trim the resulting agar block with a razor blade, and place in a 3.5 ml borosilicate screw-cap vial (Fisher 03-339-21B) containing 1.5ml 2%paraformaldehyde/2%glutaraldehyde.

**Note**

All subsequent incubations and washes use 1.5 -2.0 ml and are performed in the same vial.











- 5 Fix colonies by incubating for 7 days at 4°C
- 6 Wash #1: Wash agar blocks on ice by incubating for 15 minutes with 1.5 ml of 0.15M sodium cacodylate (pH 7.2)  
 00:15:00
- 7 Wash #2: Wash agar blocks on ice by incubating for 15 minutes with 1.5 ml of 0.15M sodium cacodylate (pH 7.2)  
 00:15:00
- 8 Wash #3: Wash agar blocks on ice by incubating for 5 minutes with 1.5ml OS buffer (100 mM KH<sub>2</sub>PO<sub>4</sub>, 10 mM MgCl<sub>2</sub>, pH 6.0).  
 00:05:00
- 9 Wash #4: Wash agar blocks on ice by incubating for 5 minutes with 1.5ml OS buffer (100 mM KH<sub>2</sub>PO<sub>4</sub>, 10 mM MgCl<sub>2</sub>, pH 6.0).  
 00:05:00
- 10 If sections will be used for electron microscopy, add 1% OsO<sub>4</sub> in OS to vials to cover the agar blocks and incubate on ice in a chemical fume hood for 1 hr. Otherwise skip to step 13.  
 01:00:00



- 11 Dispose of 1% OsO<sub>4</sub> in hazardous waste.  
WashA: Wash with 1.5ml OS buffer by incubating on ice for 10 minutes.
- 12 WashB: Wash with 1.5ml OS buffer by incubating on ice for 10 minutes.
- 13 Add 1.5 ml OS and incubate overnight at 4°C.  
 18:00:00
- 14 WashA: Wash blocks with 1.5ml cold water by incubating on ice for 10 minutes.
- 15 WashB: Wash blocks with 1.5ml cold water by incubating on ice for 10 minutes.
- 16 Wash #1: Add 1.5ml cold 25% ethanol and incubate on ice for 10 minutes.  
 00:10:00
- 17 Wash #2: Add 1.5ml cold 50% ethanol and incubate on ice for 10 minutes.  
 00:10:00
- 18 Wash #3: Add 1.5ml cold 75% ethanol and incubate on ice for 10 minutes.  
 00:10:00
- 19 Wash #4: Add 1.5ml cold 95% ethanol and incubate on ice for 10 minutes.  
 00:10:00
- 20 Wash #5: Add 1.5ml cold 100% ethanol and incubate on ice for 10 minutes.  
 00:10:00
- 21 Wash #6: Add 1.5ml cold 100% ethanol and incubate on ice for 10 minutes.  
 00:10:00
- 22 Remove ethanol and resuspend in 1.5ml cold 100% ethanol. Leave the blocks overnight at 4°C in 100% ethanol.  
 18:00:00

## Spurr's treatment



- 23 Make Spurr's reagent by stirring slowly under chemical fume hood 5 grams ERL4221, 4 grams DER736 and 13 grams NSA for 20 minutes(Electron Microscopy Sciences).
- 24 Add 0.15 grams DMAE(Electron Microscopy Sciences) and stir for 20 minutes. De-gas for 1-2 hrs.  
 00:15:00
- 25 Wash blocks with 1.5ml unopened room temp 100% ethanol by incubating at room temp for 10 minutes. Repeat 4 more times.  
 00:30:00
- 26 Treatment#1: Remove ethanol and add 1.5ml of a 2:1 ratio of 100% ethanol:Spurr's reagent.  
 00:15:00
- 27 Treatment#1: Rotate vial for 15 minutes on wheel at room temperature.  
 00:30:00
- 28 Treatment#1: Allow to stand for 30 minutes at room temperature.  
 00:15:00
- 29 Treatment#2: Remove 100%ethanol:Spurr's reagent and add 1.5ml of a 2:1 ratio of 100% ethanol:Spurr's reagent. Rotate vial for 15 minutes on wheel at room temperature.  
 00:30:00
- 30 Treatment#2: Allow to stand for 30 minutes at room temperature.  
 00:15:00
- 31 Treatment #3: Remove 100% ethanol:Spurr's reagent and add 1.5ml of a 2:1 ratio of 100%ethanol: Spurr's reagent. Rotate vial for 15 minutes on wheel at room temperature.  
 00:30:00
- 32 Treatment #3: Allow to stand for 30 minutes at room temperature.  
 00:15:00
- 33 Treatment #4: Remove 100% ethanol: Spurr's reagent and add 1.5ml of a 1:1 ratio of 100%ethanol:Spurr's reagent. Rotate vials for 15 minutes on wheel at room temperature.  
 00:30:00
- 34 Treatment #4: Allow to stand for 30 minutes at room temperature.



04:00:00

- 35 Treatment#5: Remove 100% ethanol:Spurr's reagent and add 1.5ml of a 1:1 ratio of 100%ethanol: Spurr's reagent. Rotate vials for 15 minutes on wheel at room temperature.

18:00:00

- 36 Treatment#5: Allow to stand for 30 minutes at room temperature.

18:00:00

- 37 Remove 100% ethanol:Spurr's reagent and add 1.5ml Spurr's reagent to vial. Incubate for 4hrs at room temperature.

04:00:00

- 38 Replace with 1.5ml Spurr's reagent and rotate the vial overnight on the wheel at room temperature.

18:00:00

- 39 Replace with 1.5ml Spurr's reagent and rotate the vial until late afternoon on the wheel at room temperature.

- 40 Remove the Spurr's reagent and replace with 1.5ml freshly made Spurr's reagent. Rotate the vial overnight on the wheel at room temperature.

00:00:15

- 41 The next day replace with 1.5ml Spurr's reagent and rotate until the following day on the wheel at room temperature.

- 42 The next day replace with 1.5ml Spurr's reagent and rotate until the following day on the wheel at room temperature.

- 43 Place each agar block in a mold with 0.2ml Spurr's reagent and incubate at 60°C for four hours.

- 44 Top off the molds with Spurr's reagent and incubate at 60°C for 3 days.

- 45 Collect sections(0.5u) from the central region of the colony in a drop of dH2O on a glass slide. Dry slide on a 52°C heat block.

- 46 Stain sections with 1% toluidine blue,1%Sodium Borate for 5-15 seconds.



- 47 Wash slide under a stream of dH<sub>2</sub>O. Dry on heat block. Cover in Permount(Fisher SP15-100). Examine by light microscopy.