



Jun 07, 2015

Embryoid bodies generation

DOI

dx.doi.org/10.17504/protocols.io.cy7xzm

Vivian Liu



Vivian Liu

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.cy7xzm>

Protocol Citation: Vivian Liu 2015. Embryoid bodies generation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.cy7xzm>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 05, 2015

Last Modified: November 09, 2017



Protocol Integer ID: 767

Keywords: embryoid bodies generation

Guidelines

MATERIALS

Medium: NO LIF ES medium

Dish: BD Petri Dish.

Trypsin

4% PFA

PBS

Troubleshooting



Seed

- 1 Add Trypsin 2ml and gently shake plate
- 2 Add 8ml ES medium without LIF to inactivate trypsin
- 3 Centrifuge at 500rpm for 10min
 00:10:00
- 4 Add 8ml of EB medium to each 10cm petri-dish while centrifuging
- 5 Gently resuspend pellet in 5ml EB medium by 5ml pipette, up and down for 5 times
- 6 Take 100ul cell suspension and dilute into 1ml. Count.
- 7 Take appropriate number of cells and dilute into 1×10^6 cells/ml
- 8 Gently re-suspend by 5ml pipette, up and down for 4 times
- 9 Seed 2ml of cells/dish
- 10 Save 2×10^6 cells for RNA extraction

Changing Medium

- 11 Transfer medium containing Ebs to 50ml tubes

Note

Medium should be changed every other day



12 Carefully add 5ml medium to dishes immediately and put back in incubator

Note

Try not to disturb the attached Ebs.

13 Sink cells in tubes by gravity for 5min

 00:05:00

Note

It may need **longer** @ first time as EB is small at this time point.

14 Discard supernatant

15 Add 5ml medium to pellets

16 Mix samples by inverting a few times gently

17 Distribute samples to different dishes

Note

If Ebs settle down during distributing, invert tube gently for a few times again.

Outgrowth of Ebs

18 Prepare gelatin coated chamber slides (2 well)

19 Add 1.5ml of medium into each well

20 Transfer 1 dish of day4 Ebs into 15ml tube

21 Wait for Ebs sink to the bottom of tube by gravity for 5min



00:05:00

- 22 Discard supernatant
- 23 Add 10ml EB medium into pellets
- 24 Seed 0.5ml of EB to each well
- 25 Change medium everyday: suck supernatant, replace with fresh EB medium

Immunostaining

- 26 Take 2 wells for each group on day6, day8 and day12
- 27 Suck the medium carefully
- 28 Fix with 4% PFA for 30min @RT. (in hood)
 00:30:00
- 29 Wash #1 with 2ml of PBS
 00:05:00
- 30 Wash #2 with 2ml of PBS
 00:05:00
- 31 Wash #3 with 2ml of PBS
 00:05:00
- 32 Store samples in PBS, @ 4°C