

Feb 23, 2016

Glyoxal/borate and RNase T1-based detection of inosine residues in RNA

DOI

dx.doi.org/10.17504/protocols.io.ejjbckn

Matus Valach¹

¹Université de Montréal, Montreal, Quebec, Canada



Matus Valach

Université de Montréal, Montreal, Quebec, Canada

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ejjbckn>

Protocol Citation: Matus Valach 2016. Glyoxal/borate and RNase T1-based detection of inosine residues in RNA. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ejjbckn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: February 16, 2016



Last Modified: March 24, 2018

Protocol Integer ID: 2379

Keywords: inosine residue by rnase t1, inosine residues in rna, rna adenosine, rnase t1, detection of inosine residue, inosine residue, mitochondrial transcript, rna, mitochondrial transcripts of diplomonid protist, important modification of rna, inosine, many metazoan nuclear transcript, several types trna, diplomonid protist, recent adaptation of cattenoz et al, digestion at g, glyoxal

Abstract

Adenosine-to-inosine (A-to-I) deamination is a functionally important modification of RNA that occurs in many metazoan nuclear transcripts, in several types tRNAs, and in mitochondrial transcripts of diplomonid protists. Several conceptually different methods exist to identify inosine residues in RNA, including the approach that exploits the cleavage after a guanosine or inosine residue by RNase T1. The digestion at Gs is blocked after an incubation in a glyoxal/borate solution and leads predominantly to cleavage 3' to Is. This protocol blends the procedure originally conceived by Morse and Bass (DOI: 10.1021/bi9709607) and its more recent adaptation of Cattenoz *et al.* (DOI: 10.1261/rna.036202.112).

Guidelines

- The primer extension assay to map inosine residues achieves single-nucleotide resolution, but requires more material. The RT-PCR assay is more sensitive, but for single-nucleotide resolution, it needs to be coupled to 5' adapter-ligation and sequencing. For a practical primer extension with sufficiently intense bands for cleavage intermediate detection, the target RNA in a sample should amount to >0.1 pmol (>10 ng of a 300 nt long transcript), optimally >0.5 pmol.
- The described procedure and quantities work well for a transcript, which represents 1-2% of the RNA sample. For the primer extension assay, if more RNA sample is required to achieve >0.1 pmol of the target, scale up all components, so that guanosine residues in the sample RNA are evenly labeled with glyoxal-borate adducts.
- If a purified or synthetic RNA is being analyzed, complement the sample to the final amount of 2.5 µg RNA with non-target RNA (e.g. total yeast tRNAs) to avoid over-digestion with RNase T1. Alternatively, for an RT-PCR assay, use 1 µg of poly(A) RNA as a substrate.



Materials

MATERIALS

☒ DMSO Merck MilliporeSigma (Sigma-Aldrich) Catalog #472301

☒ Ethanol

☒ Glyoxal (40%) Merck MilliporeSigma (Sigma-Aldrich) Catalog #50649

☒ Sodium Phosphate monobasic

☒ Boric acid Fisher Scientific Catalog #BP1681

☒ Sodium Hydroxide Fisher Scientific Catalog #BP359500

☒ Tris(hydroxymethyl)aminomethane Thermo Fisher Scientific Catalog #327360010

☒ RNase T1 Thermo Scientific Catalog #EN0541

☒ Hydrochloric Acid

Troubleshooting

Safety warnings

- ❗ DMSO is a very efficient solvent and easily enters through skin. Glyoxal can cause oral, dermal, and respiratory difficulties. Home-made Trizol substitute contains phenol, guanidine thiocyanate, and ammonium thiocyanate, which are all toxic. Work in gloves and use the fume hood.

Before start

- Prepare deionized glyoxal (see the protocol: [Preparation of deionized glyoxal](#)).
- Prepare 100 mM sodium phosphate buffer, pH 7.0.
- Prepare 1 M sodium borate buffer, pH 7.5.
- Prepare Tris-borate buffer (10 mM Tris-HCl, pH 7.8; 1 M sodium borate, pH 7.5).

- 1 Mix the following components (100 μ l):

	Component	Amount [μl]	Final concentration
	2.5 μ g RNA (in water)	34	25 ng/ μ l
	100 mM Sodium Phosphate, pH 7.0	10	10 mM
	DMSO	50	50%
	40% Deionized Glyoxal	6	2.4%

- 2 Incubate for 45 min at 37 °C.
- 3 Dilute the reaction by adding an equal volume (100 μ l) of 1 M sodium borate, pH 7.5.
- 4 Immediately precipitate with five volumes of 100% ethanol (relative to the diluted mixture; 1000 μ l).
- 5 Spin at >15,000 g for 30 minutes at 4 °C.
- 6 Discard the supernatant.
 - *A white pellet of 10-20 μ l is expected to have formed.*
- 7 Add 1 volume of 70% (aq.) ethanol to wash the pellet.
- 8 Directly proceed to a spin at >15,000 g for 5 minutes at 4 °C.
- 9 Discard the supernatant.
- 10 Let the pellet dry at room temperature, so that no traces of liquid remain (30-60 minutes). Alternatively, dry under vacuum for 10 minutes (without heating).
- 11 Solubilize the pellet in 40 μ l of Tris-borate buffer (10 mM Tris-HCl, pH 7.8; 1 M sodium borate, pH 7.5). The final volume of the resuspended sample should be 47-48 μ l.

- *After adding the buffer, wait for 5-10 minutes before resuspending the pellet. This makes the pellet less sticky and prevents sample loss.*
- 12 Add RNase T1 to a desired final concentration. If required, dilute the RNase in the Tris-borate buffer. Use the Tris-borate buffer to complete the volume to 50 µl.
 - *Optimize the amount using the enzyme concentration gradient by monitoring the progress of the reaction by a Northern blot or RT-PCR. If the inosine content of the investigated site is close to 100%, 1000 U of RNase T1 will cut essentially all molecules of a transcript present at ~0.5 pmol or less in 3 minutes at 37 °C.*
 - 13 Incubate for 3 minutes at 37 °C.
 - 14 Dilute the reaction with an equal volume of ice-cold RNase-free water and place on ice.
 - **Pause point:** *Proceed directly to the next step or snap-freeze the reactions in liquid nitrogen and store at -80 °C.*
 - 15 Extract the RNA, e.g. using the home-made Trizol substitute (see the corresponding protocol: **RNA extraction using the 'home-made' Trizol substitute**).
 - 16 Solubilize the pelleted RNA in 50 µl of 100 mM sodium phosphate pH 7.0.
 - 17 Add an equal volume (50 µl) of 100% DMSO and mix well.
 - 18 De-glyoxalate the RNA by incubating it for 3 hours at 65 °C (with gentle mixing every 15-30 minutes).
 - 19 Dilute the reaction with an equal volume of ice-cold RNase-free water (100 µl).
 - 20 Extract the RNA, e.g. using the home-made Trizol substitute (see the corresponding protocol: **RNA extraction using the 'home-made' Trizol substitute**).
 - 21 Resuspend the RNA in RNase-free water. Store at -80 °C or proceed directly to downstream analysis (Northern blot hybridization, primer extension, RT-PCR, RNA-Seq).