



# MojoSort™ Isolation Kits No Wash Protocol

DOI

[dx.doi.org/10.17504/protocols.io.e2wbgfe](https://dx.doi.org/10.17504/protocols.io.e2wbgfe)



Kelsey Miller<sup>1</sup>

<sup>1</sup>BioLegend

BioLegend

Tech. support email: [tech@biolegend.com](mailto:tech@biolegend.com)



Sam Li

BioLegend

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)



DOI: <https://dx.doi.org/10.17504/protocols.io.e2wbgfe>

External link:

[http://www.biolegend.com/media\\_assets/support\\_protocol/MojoSort\\_Isolation\\_Kits\\_No\\_Wash\\_Protocol\\_V03\\_10192015.pdf](http://www.biolegend.com/media_assets/support_protocol/MojoSort_Isolation_Kits_No_Wash_Protocol_V03_10192015.pdf)

**Protocol Citation:** Kelsey Miller . MojoSort™ Isolation Kits No Wash Protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e2wbgfe>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Created:** June 03, 2016

**Last Modified:** August 29, 2019

**Protocol Integer ID:** 2870

**Keywords:** sample with the biotin antibody cocktail, sample with the biotin, biotin antibody cocktail, wash protocol target cell, biotin, untouched cell, incubation with magnetic streptavidin nanobead, antibody cocktail, target cell, functional assay, wash protocol, magnetic streptavidin nanobead, antibody, mojosort, cell, streptavidin nanobead, cells of interest, nanobead

## Abstract

Target cells are depleted by incubating your sample with the biotin antibody cocktail followed by incubation with magnetic Streptavidin Nanobeads (Cat. No. 480015/480016). The magnetically labeled fraction is retained by the use of a magnetic separator. The untouched cells are collected. These are your cells of interest; do not discard the liquid. Some of the downstream applications include functional assays, gene expression, phenotypic characterization, etc.


## Guidelines

### Reagents and instruments required:

- MojoSort™ Buffer (5X) (Cat. No. 480017)
- MojoSort™ Magnet (Cat. No. 480019) or compatible magnetic separation system
- Adjustable pipettes
- 70 µm filters (one per sample)
- 5 mL (12 × 75 mm) polystyrene tubes
- Reagents for sample preparation
- Reagents and instruments (Flow cytometer) to determine yield and purity

### Protocol:

This protocol has been optimized to remove washing steps after antibody cocktail and nanobeads incubations, resulting in a shorter and more convenient protocol. This procedure is optimized for the isolation of  $10^7$  to  $2 \times 10^8$  cells per tube. If working with fewer than  $10^7$  cells, keep volumes as indicated for  $10^7$  cells. For best results, optimize the conditions to your specific cell number and tissue. Prepare fresh MojoSort™ Buffer solution by diluting the 5X concentrate with sterile distilled water.



## Troubleshooting



- 1 Prepare cells from your tissue of interest without lysing erythrocytes
- 2 In the final wash of your sample preparation, resuspend the cells in MojoSort™ Buffer by adding up to 4 mL in a 5 mL (12 × 75 mm) polystyrene tube.  
*Note: Keep MojoSort™ Buffer on ice throughout the procedure*
- 3 Filter the cells with a 70 µm cell strainer, centrifuge at 300 × g for 5 minutes, and resuspend in an appropriate volume of MojoSort™ Buffer. Count and adjust the cell concentration to  $1 \times 10^8$  cells/mL.

00:05:00

- 4 Aliquot 100 µL of cell suspension ( $10^7$  cells) into a new tube. Add 10 µL of the Biotin-Antibody Cocktail, mix well and incubate on ice for 15 minutes. Scale up the volume accordingly if separating more cells. For example, add 100 µL for  $1 \times 10^8$  cells. When working with less than  $10^7$  cells, use indicated volumes for  $10^7$  cells.

**Optional:** Take an aliquot before adding the cocktail to monitor purity and yield.

00:15:00

- 5 Resuspend the beads by vortexing, maximum speed, 5 touches. Without washing, add 10 µL of Streptavidin Nanobeads. Mix well and incubate on ice for 15 minutes. Scale up the volume accordingly if separating more cells. For example, add 100 µL for  $1 \times 10^8$  cells. When working with less than  $10^7$  cells, use indicated volumes for  $10^7$  cells

00:15:00

- 6 Add 3 mL of MojoSort™ Buffer.

*Note: If you observe aggregates, filter the suspension. To maximize yield, you can disrupt the aggregates by pipetting the solution up and down.*

- 7 Place the tube in the magnet for 5 minutes

00:05:00

- 8 Pour out and collect the liquid. These are your cells of interest; **DO NOT DISCARD.**

- 9 If needed, add 3 mL of MojoSort™ Buffer and repeat steps 7 and 8 with the magnetically labeled fraction up to two times, and then pool the unlabeled fractions.

*Note: Repeating the magnetic separation increases the yield, without a strong impact on the purity. The yield will typically increase about 8 – 10% with a second separation, and about 2 – 5% with a third separation. The purity may decrease 1 – 2% with each separation. Optional: Take a small aliquot before placing the tube in the magnet to*

*monitor purity and yield. Keep unused cells to be used as control or other applications if needed.*