



Mar 15, 2016

One-step growth curves for Cellulophaga phages

One-step growth curves for Cellulophaga phages

1



dx.doi.org/10.17504/protocols.io.eqfbdtm

VERVE Net Sullivan Lab



Bonnie Poulos

Matthew Sullivan Lab, University of Arizona, The Ohio State ...

One-step growth curves are used to make determinations about the life cycle of a virus on a particular host. By following a virus infection during one life cycle phase of host a growth curve can be constructed and the burst size can be calculated.

DOI

dx.doi.org/10.17504/protocols.io.eqfbdtm

Bonnie Poulos: One-step growth curves for Cellulophaga phages. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eqfbdtm>



Added points for clarification.

One-step growth curves for Cellulophaga phages, Verve Team

_____ protocol ,

Mar 15, 2016

Mar 28, 2018

2535

Note: One-step growth experiment instructions are for MOI=0.1.

Before performing a one-step experiment, you need to determine when the host is growing exponentially. You must also determine the titer of the phage lysate that will be used for the experiment. These are described in the protocol.

Host Growth Curve

- 1 Inoculate a new culture; ie, pick a colony into a new flask containing 10 ml of MLB media

Cellulophaga should grow at room temperature on the benchtop.

- 2 Immediately after the transfer, take a 'time 0' growth reading (T0)

Take the T0 growth reading at 595nm using 200ul of culture, in duplicate. Make sure the microtiter plate you are using is clean inside and out, with no scratches or spots on its surface, as it will interfere with the light reading. Subtract the media blank reading from the culture reading to arrive at the OD595 for the culture.



Cellulophaga growth reading
by Verve Team,
University of Arizona

PREVIEW

RUN



- 2.1 Pipet 200 μ l of MLB media into wells A1 and A2 of a white microtiter plate

This is your 'blank'.

- 2.2 Pipet 200 μ l of sample (the new culture you just inoculated) into wells B1 and B2 of the same plate

Ensure that there are no bubbles in the wells, as they will affect your readings. Pipet away any bubbles.

- 2.3 Read the plate on the plate reader

- 3 Continue taking readings as performed in step 2 periodically

Graph the results as you go! It is best to infect the host in med-exponential (log linear) phase, when $OD \cong 0.02$.

You can start with longer intervals (1–2 hours) until you start to see growth, then shorter intervals (15–30 minutes) until the growth starts to level off. If it's taking a while, you can go back to reading at longer intervals.

Phage Lysate Titer

4 Do a plaque assay to determine the PFU/ml of the lysate you plan to use

5 Calculate the volume needed for 10^7 phages

If 10^7 phages are contained in less than $1\mu\text{l}$, you will need to dilute the lysate prior to performing the growth curve experiment.

One-Step Growth Experiment

6 Determine the concentration of your culture at the time you start the infection

Use a correlation of readings from the plate reader and cell counts (CFU, DAPI, or FCM counts) to estimate this.

7 Calculate the volume of host culture needed for 10^8 cells

8 Pipet this amount into a 1.5 ml tube

9 Add 10^7 phages to the tube and start your timer for 15 minutes to allow the phages to adsorb to the host cells


🕒 00:15:00

10 After 15 minutes, dilute the infection 1:1000 in MLB media in a 250 ml flask.

If you have 50mL of MLB in the flask, add 50ul of host for a 1:1000 dilution.

11 Take a sample immediately after dilution – this is ‘time 0’

12 Steps for centrifuged sample:

 Centrifuged Sample Steps
by Verve Team,
University of Arizona

[PREVIEW](#) [RUN](#)

12.1 Pipet 100 µl from the flask into 900 µl of MSM in a 15 ml tube (you are diluting your sample 10x: 10^{-1})

Once you know how many phages to expect, you know what dilutions of your early samples to plate to get good counts For example, if the T0 expected concentration is 10^4 , there should be 100 plaques if you plate 100 µl of a 10^{-1} dilution.

Note that if you are using a different MOI, you will need to calculate the expected number of phage at T0 to guide you in what dilution to plate This will depend on the total volume of the initial infection (ie, the volume of cells plus phages) So the concentration at T0 should be total phage added/volume of infection, divided by 1000 (for the 1:1000 dilution) Convert this to phages per ml.

12.2 Vortex briefly

12.3 Centrifuge at 5 min at 1000 rpm

 00:05:00

12.4 Very carefully remove the tube (do not disturb the pellet!) and plate 100 µl

13 Steps for samples that are not centrifuged:



Non-Centrifuged Sample Steps
by Verve Team,
University of Arizona

PREVIEW

RUN



13.1 Pipet 100 μ l from the flask into 900 μ l of MSM in a 1.5 ml tube

13.2 Vortex briefly

13.3 Plate 100 μ l

14 Continue sampling in this way for 8 hours

08:00:00

At later time points, more dilutions will need to be plated. On the first trial, be generous with what you plate (ie, plate 10^{-2} , 10^{-3} , 10^{-4} , 10^{-5}) and use the results as a guide for what you should plate in repeat experiments.

T1 and T2 can be at 1 and 2 hours, respectively, but then switch to every 30 minutes for the duration of the experiment.

15 Store the filtered samples at 4°C

16 The next day, count the plaques on all plates that have a countable number of them

17 Decide which dilution gives the best count at each time point for the next time you do this same phage-host pair

Depending on the size of the plaques, a "good" count will be somewhere between 10 and a few hundred.

18 The next day, count any new plaques that have appeared and add these to your original count

19 Count again on the third day

20 Calculate PFU/ml at each time point for both the centrifuged (free phage only) and not centrifuged (total phage) samples

21 Graph the results

22 Calculate burst size



Calculating burst size
by Verve Team,
University of Arizona

PREVIEW

RUN



22.1 Take the FREE phage average of the time points on the plateau before the burst (A)

22.2 Take the FREE phage average of the time points on the plateau after the burst (B)

22.3 Subtract A from B; This is the total burst or new phages released (C)

22.4 Divide C by the number of infecting phage (TOTAL phages at T0 minus FREE at T0); This is the burst size

