

Jan 30, 2015

Version 1

PCR with Phusion® Polymerase (M0530) V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.cdjs4m>



New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com



New England Biolabs

New England Biolabs

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.cdjs4m>

External link: <https://www.neb.com/protocols/1/01/01/pcr-protocol-m0530>

Protocol Citation: New England Biolabs: PCR with Phusion® Polymerase (M0530). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.cdjs4m>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

Created: July 11, 2014

Last Modified: March 28, 2018

Protocol Integer ID: 139

Keywords: Phusion High Fidelity DNA Polymerase, phusion dna polymerase, successful PCR using HiFi DNA polymerase, fidelity dna polymerase, pcr protocol, pcr, m0530, dna, protocol

Abstract

This is the PCR protocol for Phusion® High-Fidelity DNA Polymerase (M0530)

Guidelines

OVERVIEW

The following guidelines are provided to ensure successful PCR using Phusion® DNA Polymerase. These guidelines cover routine PCR. Amplification of templates with high GC content, high secondary structure, low template concentrations or long amplicons may require further optimization.

PROTOCOL

Reaction Setup:

We recommend assembling all reaction components on ice and quickly transferring the reactions to a thermocycler preheated to the denaturation temperature (98°C). All components should be mixed and centrifuged prior to use. It is important to add Phusion DNA Polymerase last in order to prevent any primer degradation caused by the 3' → 5' exonuclease activity. Phusion DNA Polymerase may be diluted in 1X HF or GC Buffer just prior to use in order to reduce pipetting errors. Please note that protocols with Phusion DNA Polymerase may differ from protocols with other standard polymerases. As such, conditions recommended below should be used for optimal performance.

Component	20 µl Reaction	50 µl Reaction	Final Concentration
Nuclease-free water	to 20 µl	to 50 µl	
5X Phusion HF or GC Buffer	4 µl	10 µl	1X
10 mM dNTPs	0.4 µl	1 µl	200 µM
10 µM Forward Primer	1 µl	2.5 µl	0.5 µM
10 µM Reverse Primer	1 µl	2.5 µl	0.5 µM
Template DNA	variable	variable	< 250 ng
DMSO (optional)	(0.6 µl)	(1.5 µl)	3%
Phusion DNA Polymerase	0.2 µl	0.5 µl	1.0 units/50 µl PCR

Notes: Gently mix the reaction. Collect all liquid to the bottom of the tube by a quick spin if necessary. Overlay the sample with mineral oil if using a PCR machine without a heated lid.

Transfer PCR tubes from ice to a PCR machine with the block preheated to 98°C and begin thermocycling:

Thermocycling conditions for a routine PCR:

STEP	TEMP	TIME
------	------	------

Initial Denaturation	98°C	30 seconds
25-35 Cycles	98°C45-72°C72°C	5-10 seconds10-30 seconds15-30 seconds per kb
Final Extension	72°C	5-10 minutes
Hold	4-10°C	

General Guidelines:

1. Template:

Use of high quality, purified DNA templates greatly enhances the success of PCR. Recommended amounts of DNA template for a 50 µl reaction are as follows:

DNA	Amount
genomic	50 ng–250 ng
plasmid or viral	1 pg–10 ng

If the template DNA is obtained from a cDNA synthesis reaction, the volume added should be less than 10% of the total reaction volume.

2. Primers:

Oligonucleotide primers are generally 20–40 nucleotides in length and ideally have a GC content of 40–60%. Computer programs such as **Primer3** can be used to design or analyze primers. The final concentration of each primer in a reaction using Phusion DNA Polymerase may be 0.2–1 µM, while 0.5 µM is recommended.

3. Mg++ and additives:

Mg++ is critical to achieve optimal activity with Phusion DNA Polymerase. The final Mg++ concentration in 1X Phusion HF and GC Buffer is 1.5 mM. Excessive Mg++ can prevent full denaturation of DNA as well as cause non-specific binding of primers. The optimal Mg++ concentration is affected by dNTP concentration, the template being used and supplements that are added to the reaction. This can also be affected by the presence of chelators (e.g. EDTA). Mg++ can be optimized in 0.5 mM increments using the MgCl₂ provided.

Amplification of difficult targets, such as those with GC-rich sequences or secondary structure, may be improved by the presence of additives such as DMSO (included). A final concentration of 3% DMSO is recommended, although concentration can be optimized in 2% increments. It is important to note that if a high concentration of DMSO is used, the annealing temperature must be lowered as it decreases the primer T_m (2). Phusion DNA polymerase is also compatible with other additives such as formamide or glycerol.

4. Deoxynucleotides:

The final concentration of dNTPs is typically 200 μ M of each deoxynucleotide. Phusion cannot incorporate dUTP.

5. Phusion DNA Polymerase Concentration:

We generally recommend using Phusion DNA Polymerase at a concentration of 20 units/ml (1.0 units/50 μ l reaction). However, the optimal concentration of Phusion DNA Polymerase may vary from 10–40 units/ml (0.5–2 units/50 μ l reaction) depending on amplicon length and difficulty. Do not exceed 2 units/50 μ l reaction, especially for amplicons longer than 5 kb.

6. Buffers:

5X Phusion HF Buffer and 5X Phusion GC Buffer are provided with the enzyme. HF buffer is recommended as the default buffer for high-fidelity amplification. For difficult templates, such as GC-rich templates or those with secondary structure, GC buffer can improve reaction performance. GC buffer should be used in experiments where HF buffer does not work. Detergent-free reaction buffers are also available for applications that do not tolerate detergents (e.g. microarray, DHPLC).

7. Denaturation:

An initial denaturation of 30 seconds at 98°C is sufficient for most amplicons from pure DNA templates. Longer denaturation times can be used (up to 3 minutes) for templates that require it.

During thermocycling, the denaturation step should be kept to a minimum. Typically, a 5–10 second denaturation at 98°C is recommended for most templates.

8. Annealing:

Annealing temperatures required for use with Phusion tend to be higher than with other PCR polymerases. **The NEB Tm calculator should be used to determine the annealing temperature when using Phusion.** Typically, primers greater than 20 nucleotides in length anneal for 10–30 seconds at 3°C above the T_m of the lower T_m primer. If the primer length is less than 20 nucleotides, an annealing temperature equivalent to the T_m of the lower primer should be used. A temperature gradient can also be used to optimize the annealing temperature for each primer pair. For two-step cycling, the gradient can be set as high as the extension temperature.

For high T_m primer pairs, two-step cycling without a separate annealing step can be used.

9. Extension:

The recommended extension temperature is 72°C. Extension times are dependent on amplicon length and complexity. Generally, an extension time of 15 seconds per kb can be used. For complex amplicons, such as genomic DNA, an extension time of 30 seconds per kb is recommended. Extension time can be increased to 40 seconds per kb for cDNA templates, if necessary.

10. Cycle number:

Generally, 25–35 cycles yields sufficient product.

11. 2-step PCR:

When primers with annealing temperatures $\geq 72^\circ\text{C}$ are used, a 2-step thermocycling protocol is recommended.

Thermocycling conditions for a routine 2-step PCR:

	STEP	TEMP	TIME
	Initial Denaturation	98°C	30 seconds
	25-35 Cycles	98°C 72°C	5-10 seconds 15-30 seconds per kb
	Final Extension	72°C	5-10 minutes
	Hold	4-10°C	

12. PCR product:

The PCR products generated using Phusion DNA Polymerase have blunt ends; if cloning is the next step, then blunt-end cloning is recommended. If TA-cloning is preferred, then DNA should be purified prior to A-addition, as Phusion DNA Polymerase will degrade any overhangs generated.

Addition of an untemplated -dA can be done with Taq DNA Polymerase ([NEB #M0267](#)) or Klenow exo- ([NEB #M0212](#)).

References:

1. Chester, N. and Marshak, D.R. (1993). Analytical Biochemistry. 209, 284-290.

Materials

MATERIALS

 Phusion High-Fidelity DNA Polymerase - 100 units **New England Biolabs Catalog #M0530S**

Safety warnings

 Please note that protocols with Phusion DNA Polymerase may differ from protocols with other standard polymerases. As such, conditions recommended below should be used for optimal performance.

Before start

Annealing temperatures should be determined using the [NEB Annealing Temp Calculator](#).



1 Set up the following reaction on ice.

Component	20 μ l Reaction	50 μ l Reaction	Final Concentration
Nuclease-free water	to 20 μ l	to 50 μ l	
5X Phusion HF or GC Buffer	4 μ l	10 μ l	1X
10 mM dNTPs	0.4 μ l	1 μ l	200 μ M
10 μ M Forward Primer	1 μ l	2.5 μ l	0.5 μ M
10 μ M Reverse Primer	1 μ l	2.5 μ l	0.5 μ M
Template DNA	variable	variable	< 250 ng
DMSO (optional)	(0.6 μ l)	(1.5 μ l)	3%
Phusion DNA Polymerase	0.2 μ l	0.5 μ l	1.0 units/50 μ l PCR

Protocol

NAME



Mixture for M0530 Phusion PCR

CREATED BY

New England Biolabs

[Preview](#)

1.1 Nuclease-free water

1.2 5X Phusion HF or GC Buffer

Note

GC buffer should be used in experiments where HF buffer does not work. Detergent-free reaction buffers are also available for applications that do not tolerate detergents (e.g. microarray, DHPLC).

**Note**

5X Phusion HF Buffer and 5X Phusion GC Buffer are provided with the enzyme. HF buffer is recommended as the default buffer for high-fidelity amplification. For difficult templates, such as GC-rich templates or those with secondary structure, GC buffer can improve reaction performance.

1.3 10 mM dNTPs**Note**

Phusion cannot incorporate dUTP.



Deoxynucleotide Solution Mix - 8 umol of each **New England Biolabs Catalog #N0447S**

1.4 10 μ M Forward Primer**1.5 10 μ M Reverse Primer****1.6 Template DNA****1.7 DMSO (optional)****Note**

It is important to note that if a high concentration of DMSO is used, the annealing temperature must be lowered as it decreases the primer T_m (2).

Note

Amplification of difficult targets, such as those with GC-rich sequences or secondary structure, may be improved by the presence of additives such as DMSO (included). A final concentration of 3% DMSO is recommended, although concentration can be optimized in 2% increments.

1.8 Phusion DNA Polymerase**2 Gently mix the reaction.**



- 3 Collect all liquid to the bottom of the tube by a quick spin if necessary and overlay the sample with mineral oil if using a PCR machine without a heated lid.
- 4 Quickly transfer PCR tubes from ice to a PCR machine with the block preheated to 98°C and begin thermocycling.