



🔒 Quick Protocol for Monarch® PCR & DNA Cleanup Kit (5 ?g) (NEB #T1030) V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.ejxbcpn>



New England Biolabs¹

¹New England Biolabs



Isabel Gautreau

New England Biolabs

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: <https://dx.doi.org/10.17504/protocols.io.ejxbcpn>

External link: <https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-pcr-dna-cleanup-kit-5-g-t1030>

Protocol Citation: New England Biolabs . Quick Protocol for Monarch® PCR & DNA Cleanup Kit (5 ?g) (NEB #T1030). protocols.io <https://dx.doi.org/10.17504/protocols.io.ejxbcpn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: February 18, 2016

Last Modified: May 14, 2021

Protocol Integer ID: 2391



Keywords: Monarch PCR Cleanup Kit, Monarch DNA cleanup kit, dna cleanup kit, pcr, t1030, neb, dna, quick protocol, full protocol

Abstract

This is the "quick" version of Monarch® PCR & DNA Cleanup Kit (5 µg) Protocol (NEB #T1030). For the full protocol, please click [here](#).

Guidelines

For detailed protocol and more information, visit www.neb.com/T1030

The full protocol is available [here](#).

The video protocol is available [here](#).

Materials

MATERIALS

 Monarch® PCR & DNA Cleanup Kit (5 µg) **New England Biolabs Catalog #T1030**

Before start

Add 4 volumes of ethanol ($\geq 95\%$) to one volume of DNA Wash Buffer.

- For 50-prep kit, add 20 ml of ethanol to 5 ml of Monarch DNA Wash Buffer
- For 250-prep, kit add 100 ml of ethanol to 25 ml of Monarch DNA Wash Buffer

All centrifugation steps should be carried out at 16,000 x g (~13,000 RPM).



- 1 **Dilute sample with DNA Cleanup Binding Buffer according to the table below. Mix well by pipetting up and down or flicking the tube. Do not vortex.** A sample volume of 20–100 μl is recommended. For smaller samples, TE can be used to adjust the volume. For diluted samples larger than 800 μl , load a portion of the sample, proceed with step 2, and then repeat as necessary.

	Sample Type	Ratio of Binding Buffer: Sample	Example
	dsDNA > 2 kb (plasmids, gDNA)	2:1	200 μl : 100 μl
	dsDNA < 2 kb (some amplicons, fragments)	5:1	500 μl : 100 μl
	ssDNA (cDNA, M13)	7:1	700 μl : 100 μl

Note

A sample volume of 20–100 μl is recommended. For smaller samples, TE can be used to adjust the volume. For diluted samples larger than 800 μl , load a portion of the sample, proceed with step 2, and then repeat as necessary.

- 2 **Insert column into collection tube and load sample onto column. Spin for 1 minute at 16,000 x g, then discard flow-through.**

 00:01:00

- 3 **Re-insert column into collection tube. Add 200 μl DNA Wash Buffer (with ethanol added) and spin for 1 minute at 16,000 x g.** Discarding flow-through is optional.

 00:01:00

- 4 **Repeat Step 3.** (Step 3: Re-insert column into collection tube. Add 200 μl DNA Wash Buffer and spin for 1 minute at 16,000 x g. Discarding flow-through is optional).

 00:01:00

- 5 **Transfer column to a clean 1.5 ml microfuge tube.** Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for 1 minute.

Note

Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for 1 minute.

- 6 **Add $\geq 6 \mu\text{l}$ of DNA Elution Buffer to the center of the matrix. Wait for 1 minute, then spin for 1 minute at 16,000 x g to elute the DNA.**

 00:02:00



Note

Typical elution volumes are 6–20 μ l. Nuclease-free water (pH 7–8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated. For larger size DNA (≥ 10 kb), heating the elution buffer to 50°C prior to use can improve yield.