

Aug 03, 2016 Version 1

RNA Extraction from Drosophila Tissues using TRIzol Reagent V.1

DOI

dx.doi.org/10.17504/protocols.io.fgrbjv6

Elizabeth Allen¹

¹University of Massachusetts Medical School



Elizabeth Allen

University of Massachusetts Medical School

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.fgrbjv6

Protocol Citation: Elizabeth Allen 2016. RNA Extraction from Drosophila Tissues using TRIzol Reagent. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.fgrbjv6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: August 03, 2016



Last Modified: March 20, 2018

Protocol Integer ID: 3313

Keywords: rna extraction from drosophila tissue, rna extraction, using trizol reagent, trizol reagent, invitrogen life technologies trizol manual, trizol reagent this protocol, drosophila tissue, rna, trizol, drosophila

Abstract

This protocol is adapted from the Invitrogen Life Technologies Trizol manual.

Materials

MATERIALS

✕ Isopropanol

✕ Chloroform

✕ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

✕ Bio Plas Disposable Homogenization Pestles **Capitol Scientific Catalog #BPI-4040-PB**

✕ Microcentrifuge Tubes

✕ Temperature-regulated centrifuge

✕ Ultrapure Distilled, Nuclease Free Water

✕ Filter Tips

✕ Dry Ice

✕ 75% Ethanol

STEP MATERIALS

✕ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

✕ Chloroform

✕ Microcentrifuge Tubes

✕ Isopropanol

✕ 75% Ethanol

✕ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

✕ Chloroform

✕ Microcentrifuge Tubes

✕ Isopropanol

✕ 75% Ethanol

Protocol materials

✕ Chloroform

✕ Microcentrifuge Tubes

✕ Isopropanol

✕ 75% Ethanol

✕ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

✕ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

✕ Dry Ice

✕ Isopropanol

✕ Chloroform

✕ Bio Plas Disposable Homogenization Pestles **Capitol Scientific Catalog #BPI-4040-PB**

✕ Filter Tips

✕ Isopropanol

✕ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

✕ Chloroform

✕ Microcentrifuge Tubes

✕ 75% Ethanol

✕ Microcentrifuge Tubes

✕ Temperature-regulated centrifuge

✕ Ultrapure Distilled, Nuclease Free Water

✕ 75% Ethanol

✕ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

✕ Chloroform

✕ Microcentrifuge Tubes

✕ Isopropanol

✕ 75% Ethanol

Safety warnings

❗ TRIzol Reagent and Chloroform are toxic (inhalation, contact, and ingestion). Always use a fume hood, and wear protective clothing, eyewear, and gloves.



Before start

1. Prepare an **RNase-free** working area, wipe down barrels of micropipettes, use **filter tips** and RNase-free microcentrifuge tubes, and always wear **gloves**.
2. **Snap-freeze tissue** as it is collected by adding tissue to a microfuge tube on dry ice.
3. This protocol calls for **50-100 mg of tissue** to be homogenized in a microcentrifuge tube in **1 mL of TRIzol Reagent**. *Use caution here - homogenization with a motorized homogenizer may result in overflow.



- 1 Add 1 mL TRIzol Reagent to 50-100 mg of frozen *Drosophila* tissue in a 1.5 mL microcentrifuge tube, and homogenize immediately with a disposable plastic pestle.

 1 mL

 TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

- 2 **Centrifuge** the sample at 12,000 x *g* for 10 minutes at 4°C.

*Pellet contains ECM, polysaccharides, and high molecular weight DNA; **supernatant contains the RNA**. In high fat samples, a layer of fat collects above the supernatant.

 00:10:00

- 3 Remove and discard the fatty layer.

- 4 Transfer the cleared supernatant to a new tube.

- 5 Incubate the sample for 5 minutes at room temperature.

- 6 Add 0.2 mL of chloroform, and cap the tube securely.

 Chloroform

- 7 **Shake** the tube vigorously **by hand** for **15 seconds**.

 00:00:15

- 8 **Incubate** for 2-3 minutes at room temperature.

 00:03:00

- 9 **Centrifuge** the sample at 12,000 x *g* for **15** minutes at **4°C**.

*The mixture separates into a lower red phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. **RNA remains in colorless aqueous phase** (~50% of the total volume).

 00:15:00

- 10 Remove the aqueous phase of the sample by angling the tube at 45° and pipetting the solution out.

Place the aqueous phase into a **new tube**.



*Avoid drawing any of the interphase or organic layer into the pipette.

*The interphase and organic phenol-chloroform phases can be saved for DNA or protein isolation if desired (saved overnight at 4°C); however, the protocols for these procedures will not be discussed here. Please refer to the manufacturer's TRIzol Reagent manual.

 Microcentrifuge Tubes

- 11 Add **0.5 mL** of 100% isopropanol to the aqueous phase.

 Isopropanol

- 12 Incubate sample at room temperature for 10 minutes.

 00:10:00

- 13 **Centrifuge** at 12,000 x *g* for **10 minutes** at 4°C.

*RNA is often visible prior to centrifugation, and forms a gel-like pellet on the side and bottom of tube.

 00:10:00

- 14 Remove all supernatant from the tube, leaving the RNA pellet.

- 15 **Wash** the pellet with **1 mL** 75% ethanol.

*RNA can be stored in 75% ethanol at least 1 year at -20°C, or at least 1 week at 4°C.

 1 mL

 75% Ethanol

- 16
1. Briefly **vortex** the sample.
 2. **Centrifuge** the tube at 7,500 x *g* for 5 minutes at 4°C.
 3. Discard the wash.

 00:05:00

- 17 Vacuum or air **dry** the **RNA pellet** for 5-10 minutes.

*Do not dry the pellet by vacuum centrifuge.

*Do not allow the RNA to dry completely.

 00:10:00

- 18 **Resuspend** the RNA pellet in RNase-free water by passing the solution up and down several times through a pipette tip.



19 **Incubate** in a water bath or heat block set at 55-60°C for **10-15 minutes**.

 00:15:00

20 Proceed to downstream applications, such as DNase treatment or cDNA synthesis, or store at -70°C