

Jan 21, 2016

Screening Recombinant Clones by PCR

DOI

<https://dx.doi.org/10.17504/protocols.io.c48yzv>

Matthew Sullivan¹

¹Matthew Sullivan Lab

VERVE Net

Sullivan Lab



Verve Team

University of Arizona

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.c48yzv>

Protocol Citation: Matthew Sullivan 2016. Screening Recombinant Clones by PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.c48yzv>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working



Created: June 19, 2015

Last Modified: November 13, 2017

Protocol Integer ID: 896

Keywords: recombinant clones by pcr, screening recombinant clone, pcr

Guidelines

This protocol describes conditions using either ABI reagents or PGC Scientific Reagents. The reactions have been tested on the gradient Eppendorf, the Gene Techne and the ABI 9700 thermal cyclers and all machines use the same cycling conditions.

When testing plasmid DNA preps, use 0.5 µL of DNA as the template in place of the bacterial colony.

Recommended vector primers:

Primer	Sequence (5' to 3')	Fragment Size
pUC 18 Forward	TGT AAA ACG ACG GCC AGT	109 bp ¹²
pUC 18 Reverse	TCA CAC AGG AAA CAG CTA TGA C	

Thermal Cycling Protocol:

95°C	5 minutes	1 cycle
95°C 55°C 70°C	30 seconds 30 seconds 30 seconds*	35 cycles
70°C	7 minutes	1 cycle
4-10°C		Hold cycle

*increase to 1 minute if insert is great than 1 Kb

Reagent catalog numbers

ABI Reagents:

- dNTP mix N808-0007
- AmpliTaq Gold N808-0247

PGC Scientific Reagents

- dNTP mix 62-6113-40
- Taq DNA polymerase (w/10x Buffer containing MgCl₂) 62-6086-02

In pBlueScript these primers bind at nt 599 and nt 833 = 234 bp amplicon.

In pGEM-t Easy the primers bind at nt 2975 and nt 197 = 242 bp amplicon.



- 1 Transfer the transformed clones to a stock LB-Amp or other appropriate plated media.
- 2 Grow for 24 hr at 37°C.

 24:00:00

- 3 Prepare forward and reverse primers for the vector used for transformation at 100 ng/μl in DDH₂O. See table in guidelines for recommended primers to use with pUC18 vector.

Note

When choosing primers, find ones that flank the multiple cloning region of your vector

- 4 Label thin-walled tubes.

Note

Include a vector control and a known plasmid control when first starting out. The vector control will yield the smallest fragment in comparison to your transformed clones and the plasmid control will yield a fragment that is the size of the vector control plus the size of the insert.

- 5 Prepare 1 ml of Master Mix:

ABI Reagents			PGC Scientific Reagents			Final Concentration
DDH ₂ O	690 μl		DDH ₂ O	833.5 μl		-----
10x Buffer II	100 μl		10x Buffer w/MgCl ₂	100 μl		1x Buffer
dNTP mix	80 μl		dNTP mix	20 μl		200 μM each
25mM MgCl ₂	80 μl					2mM
Primer F	20 μl		Primer F	20 μl		0.33 μM
Primer R	20 μl		Primer R	20 μl		0.33 μM
AmpliTaq Gold	10 μl		Taq Polymerase	6.5 μl		1U/20 μl rxn.

**Note**

This protocol describes conditions using either ABi reagents or PGC Scientifics reagents.

6 Vortex and spin Master Mix and dispense 20 µl per thin-walled tube.

7 Cap all tubes.

8 Set a 10 µl pipetter to 1 µl.

Note

You do not need to use barrier tips.

9 Place transformant plate on template and touch the pipette tip to the appropriate colony.

10 Pipette up and down in the appropriate PCR tube.

11 Re-cap tube, change tip and repeat steps 9-10 until all colonies have been added.

Note

The colony is your DNA template.

Note

When testing plasmid DNA preps, use 0.5 µl of DNA as the template in place of the bacterial colony.

12 Perform PCR using the thermal cycling protocol:

	95°C	5 minutes	1 cycle
	95°C 55°C 70°C	30 seconds 30 seconds 30 seconds*	35 cycles
	70°C	7 minutes	1 cycle



	4-10°C		Hold cycle
--	--------	--	------------

*increase to 1 minute if insert is great than 1 Kb

13 Remove tubes at end of thermal cycling.

14 Add 2 μ l 10x gel loading buffer to each tube.

 2 μ L

15 Spin to bring contents to bottom of tube and load 10 μ l into lane of a 1% agarose gel in 0.5x TBE buffer.

Note

Be sure to have a DNA ladder in the first and last lane of each gel so that an accurate assessment of the size of the inserts can be determined.