

Aug 10, 2016

Wash Buffer (10mM Tris-HCl pH 7.5, 80% EtOH)

DOI

dx.doi.org/10.17504/protocols.io.fjmbkk6

Sarah Hessen-Schmidt¹

¹University of Southern California



Sarah Hessen-Schmidt

University of Southern California

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.fjmbkk6

Protocol Citation: Sarah Hessen-Schmidt 2016. Wash Buffer (10mM Tris-HCl pH 7.5, 80% EtOH). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.fjmbkk6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: August 10, 2016

Last Modified: March 28, 2018

Protocol Integer ID: 3405

Materials

STEP MATERIALS

 Ethyl alcohol, Pure 200 proof, for molecular biology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**

 Ethyl alcohol, Pure 200 proof, for molecular biology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**



Protocol materials



Ethyl alcohol, Pure 200 proof, for molecular biology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**



Ethyl alcohol, Pure 200 proof, for molecular biology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**



Ethyl alcohol, Pure 200 proof, for molecular biology **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**



- 1 In an autoclaved bottle add 2.5ml of 1M Tris-HCl pH 7.5 solution

Protocol



NAME

1M Tris-HCl

CREATED BY

Sarah Hessen-Schmidt

PREVIEW

- 1.1 Dissolve 30.275g Trizma base in 200ml MilliQ

Trizma Base **Catalog #93362**

- 1.2 Adjust to pH 7.5 with HCl (~35ml)

- 1.3 Adjust total volume to 250ml with MilliQ

- 1.4 Autoclave

- 2 Add 200ml 100% EtOH

Ethyl alcohol, Pure 200 proof, for molecular biology **Merck MilliporeSigma**
(Sigma-Aldrich) **Catalog #E7023**

- 3 Adjust volume to 250ml by adding 47.5ml MilliQ