


Feb 13, 2024

A/Ci method using Dynamic Assimilation Technique

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5xpj4g1b/v1>



Mauricio Tejera-Nieves¹, Luke M. Gregory^{1,2}, Aaron Saathoff³, Berkley J. Walker^{1,2}

¹MSU-DOE Plant Research Laboratory; ²Plant Biology Department, Michigan State University;

³Veterinary Medicine Research and Development, Zoetis Inc.



Mauricio Tejera-Nieves

MSU-DOE Plant Research Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5xpj4g1b/v1>

Protocol Citation: Mauricio Tejera-Nieves, Luke M. Gregory, Aaron Saathoff, Berkley J. Walker 2024. A/Ci method using Dynamic Assimilation Technique. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5xpj4g1b/v1>

Manuscript citation:

Mauricio Tejera-Nieves, Do Young Seong, Lucas Reist, Berkley J. Walker*. The Dynamic Assimilation Technique measures photosynthetic CO₂ response curves with similar fidelity as steady-state approaches in half the time. Submitted to *Journal of Experimental Botany* (MS ID#: JEXBOT/2023/311662)

*Corresponding author: berkley@msu.edu

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 12, 2024

Last Modified: February 14, 2024

Protocol Integer ID: 93453

Keywords: Photosynthesis, rubisco carboxylation limitation, gas exchange methods, RuBP regeneration limitation, A/Ci response, Triose phosphate limitation, sink limitation, climate change, dynamic assimilation technique the net co₂ assimilation, net co₂ assimilation, response to intercellular co₂ concentration, fundamental measurement in photosynthesis, intercellular co₂ concentration, plant physiology research, consistent across all co₂ protocol, photosynthesis, co₂ protocol, using dynamic assimilation technique, arabidopsi, oxygenase carboxylation limitation, estimations of biochemical limitation, arabidopsis thaliana

Funders Acknowledgements:

U.S. Department of Energy, Office of Science, Office of Biological and Environmental Research

Grant ID: DE-SC0018409

U.S. Department of Energy, Office of Science, Basic Energy Sciences

Grant ID: DE- FG02-91ER20021

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Abstract

The net CO₂ assimilation (A) response to intercellular CO₂ concentration (C_i) is a fundamental measurement in photosynthesis and plant physiology research. Here, we optimized and compile a non-steady state Dynamic Assimilation Technique (DAT) protocol for A/C_i measurements. We tested this protocol in Tobacco (*Nicotina tabacum*), arabidopsis (*Arabidopsis thaliana*; Col-0), soybean (*Glycine max* 'AG21XF1', Asgrow, New Haven, CT, USA), potato (*Solanum tuberosum*; 'PA01N32-1'), apple (*Malus domestica* 'Jonagold Decoster PP#80' on 'EMLA 106' rootstock) and the extremophile plant *Rhyza Stricta*, and compare the parameter estimations against conventional steady state A/C_i measurements. The DAT protocol reduced the measurement time by almost half without compromising estimations accuracy or precision. Estimations of biochemical limitations to photosynthesis were very consistent across all CO₂ protocols, with slight differences in ribulose 1·5- biphosphate carboxylase/oxygenase carboxylation limitation.

Materials

The A/Ci measurement method using the dynamic assimilation technique requires a LI-6800 (Li-cor BioSciences, Nebraska, USA) equipped with either the Small Light Source (6800-02) or Multiphase Flash Fluorometer (6800-01A) chamber. In case of using the flash fluorometer chamber, all fluorescence features must be turned off.

Equipment	
LI-6800	NAME
Portable Photosynthesis System	TYPE
LI-COR BioSciences	BRAND
LI-6860	SKU
https://www.licor.com/env/products/photosynthesis/LI-6800/?pk_campaign=18941898409&pk_kwd=licor%20co2&pk_source=google&pk_medium=cpc&pk_content=635270912316&gad_source=1&gclid=CjwKCAiA75itBhA6EiwAkho9e9fB2v2ddF937QTG-ipuETzk3NJhDIW7ULWEwp3o7paLJjeojIIPJRo	LINK

Equipment	
Small Light Source (6800-02)	NAME
LI-6800 chamber head and light source	TYPE
LICOR BIOSCIENCES	BRAND
LI-6850	SKU
https://www.licor.com/env/products/photosynthesis/LI-6800/small-light-source	LINK

Equipment

Multiphase Flash Fluorometer (6800-01A)	NAME
LI-6800 chamber head, light source and flash fluorometer	TYPE
LI-COR BioSciences	BRAND
6800-01A	SKU
https://www.licor.com/env/products/photosynthesis/LI-6800/fluorometer	LINK



Prepare LI-6800 for measurements

1 Check instrument chemicals and run warm-up test

- 1.1 Follow LICOR [Start up checklist](#) before clamping the LI-6800 into the first leaf.
- 1.2 Run [warm-up test](#) on an empty and closed chamber
Start Up > Warmup Tests > Select Warmup Test from dropdown menu > Start

30m

Prepare LI-6800 for Dynamic measurements

- 2 **Enable Dynamic Computations:** Start Up > Chamber Setup > Check 'Add Dynamic Equations'
- 3 **Disable Fluorometer functions:** Start Up> Chamber Setup > under 'Fluorometer function as' check 'Light Source only'

- 4 The dynamic assimilation technique relies on a continuous ramping of the CO₂ concentration during the A/Ci measurement. To avoid estimation errors associated with apparent differences in the leaf sample and reference measurement channels, two main steps must be considered. Both must be performed with an empty chamber.

4.1 Acquire CO₂ and H₂O range match.

Constants > Range Match > select 'CO₂ Range Match' from dropdown menu > Start > Tap 'Continue' on the pop-up window. Repeat and select 'H₂O Range Match'.

If Range Match does not appear in the Constant Tab, double check 'Add Dynamic Equations' under Chamber Setup in the Start Up tab is checked.

Look into [Acquiring range match data](#) in the Licor Operating Instructions website for further details.

When acquiring CO₂ and H₂O range match data, it is crucial for the consistency of the measurements that ramp up values overlap reasonably close with the ramp down values.

If not the Flow_s/Flow_r parameter needs to be adjusted to compensate;

Flow_s/Flow_r should be lowered if the "coming back" (Ramp down; ▽ symbols) values are higher than the "going up" (Ramp up; Δ symbols) values, and vice-versa. Usually small changes in the Flow_s/Flow_r parameter lead to larger changes in the trajectory of the ramp up and down values, we suggest using +/- 0.02 increments.

15m

CO₂ and H₂O range matches should be done at least once before collecting new data for CO₂ response curves. If measurements are collected where there are large changes in ambient air temperature, both range matches should be done every 2 - 3 hours, or at least once in the morning and once in the afternoon.

4.2 Run CO₂ dynamic tuning

5m

Constants > Dynamic > Utilities/Test (bottom right corner of the screen) > Select 'CO₂ Test Current' from dropdown menu > Start > Tap 'Continue' on the pop-up window
Look into The dynamic tests in the Licor Operating Instructions website for further details.

Upon completion of the test, the user looks for minimum variation in A_{dyn}, and no correlation with ramp direction.

'Before' and 'After' statistical summary appears on screen next to the graph. If the standard deviation of the "After" is lower than the "Before" stats, the user should add this result by tapping the 'Yes' button. Otherwise, the user may tap 'No' and discard this

Clamp onto your first leaf

- 5 Set environmental conditions in the leaf chamber to mimic the ambient conditions the leaf was experiencing at the time of measurement. Please follow LICOR Clamping onto your first leaf instructions for a complete description.
- 6 **Open log file:** Log Setup > Open a Log File > Select folder > New File > Type file name > Done

Run Dynamic Assimilation CO₂ Response Curve

- 7 Once measurements are stable over 2 - 4 minutes, start the Dynamic Assimilation CO₂ Response Curve.
Open Dynamic Assimilation Background Program : Programs > BP builder > In the Dynamic folder select 'DAT_CO₂_Continuous.py' > Start BP
- 8 **Set Program Options.** When the Background Program is started a dialog box opens with the following options:
 1. *Starting CO₂*: Initial CO₂ concentration of the ramping. We recommend using a Low-to-High protocol to avoid oscillations in assimilation. Suggested value: 50 µmol mol⁻¹
 2. *Pre ramp wait*: Length of equilibration period at the initial CO₂ concentration before the ramping starts. Suggested value: 3 min
 3. *Ending CO₂*: Final CO₂ concentration of the ramping. Suggested value 1500 - 2000 µmol mol⁻¹



4. *Ramp rate*: Sets the rate at which CO₂ changes in the leaf chamber during the program. It affects how long it takes to complete the response curve. Suggested value 100 - 200 $\mu\text{mol mol}^{-1} \text{min}^{-1}$ (Reduce rate if assimilation oscillates after it has passed the initial linear portion)
 5. *Logging interval*: Frequency of data collection during the protocol. Suggested value = 1
 6. *When done, go to*: Final CO₂ concentration after the protocol has been completed.
Suggested value: 420 $\mu\text{mol m}^{-2} \text{s}^{-1}$
- Tap 'Continue' bottom to start the Dynamic Assimilation CO₂ Response Curve.

Close log file

- 9 Upon completion of last Dynamic Assimilation CO₂ Response Curve, close log file:
Log Setup > Logging to *Filename* > Close Log

Follow LICOR [Transferring files to a computer](#) for instructions on how to download data from the instrument.

Additional Notes

- 10 **Range CO₂ and H₂O matching:**
 - The Acquire program has 'Normal' and 'Quick' modes. The 'Quick' mode runs at a faster rate, and tends to exaggerate the difference between ramping up and down values. If the Flow_s/Flow_r parameter is optimized using the 'Quick' mode, the overlap will be even better using the 'Normal' mode.
 - The overlap for the "going up" and "coming back" values do not need to be perfect, but they should be reasonably close together. Please refer to Figure 1 for a good overlapping example.




Figure 1: Response of the difference between CO₂ concentration measured in the sample and reference channel (y-axis) to increasing (Δ symbols) and decreasing (∇ symbols) CO₂ concentrations (x-axis). Users look for changes in increasing and decreasing concentrations to be as overlapped as possible.

11 **Dynamic Assimilation Background Program considerations:**

- Suggested values may vary from plant to plant and in different environments, we recommend running a couple of test runs to optimize the program options.
- The *Starting CO₂* could be as low as 5 $\mu\text{mol mol}^{-1}$, avoid 0 $\mu\text{mol mol}^{-1}$
- Extend *Pre ramp wait* if assimilation shows a initial decrease after ramping has started
- Use a 1-second averaging window when logging data at 1 Hz as the higher data density helped overcome the additional noise. Log Setup > Logging Options > Check 'Use additional averaging time' > enter '1' on the box to the right
- *Logging interval* could be increased to 5 (one observation every 5 seconds) if data look noisy. Alternatively, moving averages or filtering could be implemented after the measurement.
- Use ambient or growing CO₂ concentration for *When done, go to*
- Matching options are ignored while running a Dynamic Assimilation CO₂ Response Curve.

12 **Steady State:**

For reliable data it is crucial the leaf biology is at, or very near, steady-state conditions in the leaf chamber before initiating measurements. Usually leaves have reached steady-state once, over a span of 2 - 4 minutes, net CO₂ shows no positive or negative trend and oscillations are within 2 - 5% of the mean.

13 **Environmental Conditions:**



Light: If ambient and chamber light levels are not adequately matched, a period of time to allow for light adaption will be needed (especially in cases where light levels are higher in the chamber than preceding ambient conditions).

Water vapor & temperature: In some cases, using a set VPD value could occasionally introduce artifacts during the CO₂ ramp. If this is the case, set reference H₂O to a level that resulted in 45-50% relative humidity in an empty chamber; leaf transpiration would add enough additional water vapor so that leaf VPD remained fairly stable throughout the DAT curve.