



Feb 10, 2025

A β ELISA

 [Nature Cell Biology](#)

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj9e6plk5/v1>

Ching-Chieh Chou¹, Judith Frydman¹

¹Stanford University



Ching-Chieh Chou

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj9e6plk5/v1>

External link: <https://doi.org/10.1038/s41556-025-01623-y>

Protocol Citation: Ching-Chieh Chou, Judith Frydman 2025. A β ELISA. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj9e6plk5/v1>



Manuscript citation:

Chou, CC., Vest, R., Prado, M.A. *et al.* Proteostasis and lysosomal repair deficits in transdifferentiated neurons of Alzheimer's disease. *Nat Cell Biol* **27**, 619–632 (2025). <https://doi.org/10.1038/s41556-025-01623-y>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 10, 2025

Last Modified: February 10, 2025

Protocol Integer ID: 119915

Keywords: ASAPCRN, a β elisa measure a β 42 levels in human neuron, a β elisa measure a β 42 level, human neuron, neuron

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Aligning Science Across Parkinson's

Grant ID: ASAP-024268

Abstract

Measure A β 42 levels in human neurons.

- 1 Human A β 42 ELISA Kit (Thermo Fisher Scientific).
- 2 Human tNeurons are first trypsinized, washed with ice-cold PBS and pelleted by centrifugation for 5 min at 1,000 X G.
- 3 Cells are lysed with RIPA buffer (25 mM Tris-HCl pH 7.5, 150 mM NaCl, 1% NP-40, 0.5% sodium deoxycholate) supplemented with protease inhibitors (Roche).
- 4 Collect the supernatants and protein concentrations are determined by BCA assay (Thermo Fisher Scientific).
- 5 Human A β 42 ELISA Kit is used to detect and quantify the A β 42 levels in total tNeuron lysates.
- 6 Prepare standards of 1000, 500, 250, 125, 62.5, 31.25, 15.63, and 0 pg/mL human A β 42 in Standard Diluent Buffer.
- 7 50 μ L of standards, control medium, cell lysates are added to each well of a 96-well plate.
- 8 Add 50 μ L anti-human A β 42 antibody per well and incubate 3 hr in the dark at room temperature with shaking.
- 9 Remove the solution and rinse wells 4 times with 1X Wash Buffer.
- 10 Add 100 μ L Anti-Rabbit IgG HRP into each well and incubate in the dark for 30 min at room temperature.
- 11 Remove the solution and rinse wells 4 times with 1X Wash Buffer.
- 12 Add 100 μ L Stabilized Chromogen to each well and incubate in the dark for 30 min at room temperature.



- 13 Add 100 μ L Stop Solution to each well. The plate is analyzed within 2 hr according to manufacturer's protocol and A β 42 values are normalized to total protein concentration of lysates.