

Apr 20, 2024

Version 2

A General Guide To Generate Different Humanized Mouse Models V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyp4zwlx9/v2>



Mohsen Khosravi-Maharlooei¹, Markus Holzl¹, Austin Chen¹, Megan Sykes¹

¹Columbia Center for Translational Immunology, Columbia University, New York

Human Islet Research N...



Sandy Beshir

City of Hope

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyp4zwlx9/v2>

Protocol Citation: Mohsen Khosravi-Maharlooei, Markus Holzl, Austin Chen, Megan Sykes 2024. A General Guide To Generate Different Humanized Mouse Models. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2lyp4zwlx9/v2> Version created by **Sandy Beshir**

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 20, 2024

Last Modified: April 20, 2024

Protocol Integer ID: 98522

Keywords: native thymus of nsg mice, human fetal thymus preparation, protocols on nsg mouse thymectomy, humanized mice, human fetal thymus preparation for more detail, type of humanized mouse, different humanized mouse model, nsg mouse thymectomy, mouse models from nsg mice, different humanized mouse models this protocol detail, fetal thymus preparation, humanized mouse, human fetal thymus, humanized mouse model, human thymus, human thymus piece, nsg mice, human cd34, isolation from human bone marrow, native thymus, cell isolation from fetal liver, cell isolation, mice, nsg mouse, thymectomy procedure, human bone marrow, mouse

Funders Acknowledgements:

NIH

Grant ID: UC4 DK104207

Abstract

This protocol details how to create humanized mouse models from NSG mice. Four different variants of humanized mice can be generated based on whether or not the native thymus is retained or if a human thymus piece is transplanted. Which type of humanized mouse is desired depends on the goals for the experiment. See our protocols on NSG mouse thymectomy, human CD34+ cell isolation, and human fetal thymus preparation for more details on some of the steps in this protocol ("[Thymectomy procedure to remove native thymus of NSG mice](#)", "[Human CD34+ cell isolation from fetal liver, and fetal thymus preparation](#)", "[CD34+ isolation from human bone marrow](#)").

Note

Corresponding Authors

Mohsen Khosravi-Maharlooei

Email: mkm2182@cumc.columbia.edu

Austin Chen

Email: ac4274@cumc.columbia.edu

Tel: 425-283-6900

- 1 NOD-scid common gamma chain knockout (NOD.Cg-Prkdcscid Il2rgtm1Wjl/SzJ) (NSG) and NSG-(K^b D^b)^{null} (IA)^{null} (NSG MHC KO) mice are obtained from the Jackson Laboratory or bred in-house and housed in a specific pathogen-free microisolator environment. Adult mice are used at 6–10 wk of age.
- 2 Discarded human fetal thymus and liver tissues of gestational age 17–20 wk are obtained without identifiers from Advanced Bioscience Resources (ABR, Alameda, CA). Informed consent from women choosing to donate fetal remains for research are obtained by ABR.
- 3 Fetal thymus fragments are cryopreserved in 10% dimethyl sulfoxide (DMSO) and 90% human AB serum (Atlanta Biologicals).
- 4 Fetal livers (FL) are cut into small pieces and incubated in  5 mL of  0.1 mg/mL liberase (Liberase™ Research Grade, 0540119001, Roche) for 20 minutes in a  37 °C CO₂ incubator. Every 5–7 minutes, the samples are vortexed and pipetted up and down to facilitate the digestion process. Digested cells are washed and CD34⁺ cells (referred to as HSCs) are positively sorted with magnetic cell sorting (MACS) (Miltenyi Biotec). From each FL, we harvest 10–30 million CD34⁺ cells with a purity of 70–90%. CD34⁺ HSCs are also cryopreserved in 10% DMSO and 90% human AB serum.
- 5 Four types of HIS mice can be generated with human FL-derived CD34⁺ HSCs: HIS mice with no thymus, with a native mouse thymus, with a human thymus graft, or with both human and mouse thymi. To generate HIS mice with no thymus, NSG mice are thymectomized and 2 weeks later, sublethally irradiated (100cGy) and injected i.v. with 1–2×10⁵ human FL-derived CD34⁺ HSCs. HIS mice with a native mouse thymus are generated similarly, except they were not thymectomized. Similar to the first group, HIS mice with a human thymus graft are thymectomized. Two weeks later, the mice are injected i.v. with HSCs and transplanted with autologous (to human HSCs) human fetal thymus fragments measuring about 1 mm³ under the kidney capsule. HIS mice with both human and mouse thymi are generated similar to the HIS mice with a human thymus, except they were not thymectomized.
- 6 To ensure that pre-existing T cells in the transplanted donor human thymus are not able to persist, we freeze and thaw the thymus tissues and also pipett them up and down several times before transplantation to release as many thymocytes as possible. To further deplete passenger thymocytes that might migrate to the periphery, an anti-human CD2 antibody is injected to the mice in 2 weekly doses (400µg/mouse, intraperitoneally).
- 7 HIS mice humanized at the neonatal stage are generated by injecting FL HSCs into neonates 1–2 days after birth (intra-liver, 10⁵ HSCs in 50µl PBS injected using a 28-gauge

needle). Neonatal mice receive half of the irradiation dose (50cGy) that the adult mice receive (100cGy).

- 8 To generate Personalized Immune (PI) HIS mice with a human thymus graft, thymectomized NSG mice are sublethally irradiated (100cGy) and injected i.v. with 2×10^5 bone marrow (BM)-derived CD34⁺ HSCs from an adult donor. Two weeks after thymectomy, the mice are transplanted with human fetal thymus fragments (matched to the BM HSCs for a class I and a class II HLA allele) measuring about 1 mm³ under the kidney capsule. Similar to generation of mice with fetal HSCs, measures are applied to deplete the pre-existing thymus graft-derived T cells.