



Jul 19, 2024

A method for RNA extraction, cDNA synthesis and and quantitative PCR from NIH-3T3 fibroblasts

 Proceedings of the National Academy of Sciences of the United States of America

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzjo84lx1/v1>

Sreeja V Nair^{1,2}, Suzanne Pfeffer^{1,2}

¹Department of Biochemistry, Stanford University School of Medicine; ²Aligning Science Across Parkinson's



Suzanne R Pfeffer

Stanford University School of Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzjo84lx1/v1>

External link: <https://doi.org/10.1073/pnas.2504774122>

Protocol Citation: Sreeja V Nair, Suzanne Pfeffer 2024. A method for RNA extraction, cDNA synthesis and and quantitative PCR from NIH-3T3 fibroblasts. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzjo84lx1/v1>

**Manuscript citation:**

Nair SV, Jaimon E, Adhikari A, Nikoloff J, Pfeffer SR Lysosomal glucocerebrosidase is needed for ciliary Hedgehog signaling: A convergent pathway contributing to Parkinson's disease. Proceedings of the National Academy of Sciences of the United States of America 122(31). doi: [10.1073/pnas.2504774122](https://doi.org/10.1073/pnas.2504774122)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 19, 2024

Last Modified: July 19, 2024

Protocol Integer ID: 103750

Keywords: ASAPCRN, mouse embryonic fibroblast, method for rna extraction, rna extraction, protocol for rna extraction, rna, cDNA synthesis

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: 000463

Abstract

Here we describe a protocol for RNA extraction, cDNA synthesis and quantitative PCR from NIH-3T3 fibroblasts. This method has also been used with mouse embryonic fibroblasts.

Materials

1. NIH-3T3 fibroblasts (Thermo Fisher Scientific #R76107)
2. Chloroform (Fisher Scientific #C298)
3. Trizol reagent (Invitrogen #15596026)
4. Isopropanol (Sigma #I9516)
5. 80% EtOH in RNase-free water
6. High-Capacity cDNA Reverse Transcription Kit (Applied Biosystems #4368814)
7. RNase OUT recombinant (Invitrogen # 100000840)
8. PowerUp SYBR Green Master Mix (Applied Biosystems #A25742)
9. Nanodrop (NanoDrop One, ThermoFisher)
10. ViiA 7 Real-Time PCR System (Applied Biosystems)



RNA extraction, cDNA synthesis and quantitative PCR from NIH-3T3 fibroblasts

1 RNA extraction

1.1 Plate 0.3×10^6 cells NIH-3T3 fibroblasts or MEF cells in individual wells of a 6-well plate.

1.2 Add 0.5 mL Trizol reagent to each well and incubate for 00:05:00 at Room temperature . Scrape the cells with a 1mL pipette tip during the 5 min incubation, to help cell lysis.

5m

1.3 Transfer the Trizol lysate to a fresh 1.5 mL tube. NOTE: Trizol lysate can be stored at -80 °C (Optional STOP step)

1.4 Add 100 μ l chloroform to the Trizol lysate and vortex for 00:00:30 .

30s

1.5 Centrifuge the samples at 14000 x g for 00:15:00 at 4 °C to separate aqueous, interphase and organic phases.

15m

1.6 Transfer the top layer aqueous phase to a fresh 1.5 mL tube and add an equal volume of isopropanol; gently tap the tubes to mix. NOTE: Use a P200 pipette set at 50 μ l to collect the aqueous phase. This helps to collect aqueous phase without disturbing the organic phase.

1.7 Store the tube at -20 °C (Optional STOP step) or leave the tubes at Room temperature for 00:05:00 and centrifuge at 13000 x g for 00:10:00 .

15m

1.8 Discard the supernatant using a 1 mL pipette add 0.5 mL 80 % (v/v) ethanol to the tube

1.9 Centrifuge at 13000 x g for 00:10:00 at Room temperature .

10m

- 1.10 Repeat steps 8-9
- 1.11 Remove the supernatant using a 1 mL pipette. Decant the remaining volume onto a Kimwipe tissue. A RNA pellet should be visible at the bottom of the tube. Care should be taken to avoid losing the RNA pellet while removing the supernatant. Allow the RNA pellet to air dry.
- 1.12 Resuspend the pellet in  50 µL RNase-free water and check the quality and concentration of RNA by Nanodrop. Purity of RNA is determined by the ratio of absorbances at 260 nm to 280 nm. A ratio >2 is considered ideal.
- 1.13 Store the RNA at  -80 °C or proceed to cDNA synthesis as described below.

cDNA synthesis using High-Capacity cDNA Reverse Transcription Kit

- 1.14
 1. Prepare a 2X master mix on ice as per the manufacturer's instructions: An example 2X master mix for 4 samples is shown below. Preparing a master mix for 4.5 samples provides adequate volume in case of pipetting error.

Component	Volume (µl) per sample	Volume (µl) for 4.5 samples
10 X RT buffer	2	9
10 X RT Random Primers	2	9
100 mM dNTP mix	0.8	3.6
RNase OUT (40 U/µl)	0.20	0.9
Reverse transcriptase (50 U/µl)	1	4
Water	4	18.5

Example Master Mix calculation

- 1.15
 - 2) Briefly spin the tube and add  10 µL of the Master Mix to labeled PCR tubes for each sample (4 tubes here) on  On ice .



- 1.16 3) Add  500 ng RNA (calculate the volume required for 500 ng from the concentration of RNA extracted) to each tube and bring the total reaction volume to  20 μ L with sterile water  On ice .
- 1.17 4) Tap the tubes once or twice to mix and set up the cDNA synthesis reaction in a PCR machine as follows: 25°C-10 min, 37°C-120 min, 85°C-5 min, 4°C-hold. Store the samples at  -20 °C or proceed to qPCR reactions on ice as described below.
- 1.18 5) Store the samples at  -20 °C or proceed to qPCR reactions on ice as described below.

Quantitative PCR

- 1.19 Dilute the cDNA 20 fold in sterile water.
- 1.20 Prepare a 2X master mix for each RNA to be analyzed. Here, as an example, we are analyzing Gli1 levels in each sample with GAPDH as an internal control. Calculations are shown below:
- 1.21

	Gli1	GAPDH
2X Master mix	90 μ l	90 μ l
10 μ M Forward Primer	3.6 μ l	3.6 μ l
10 μ M Reverse Primer	3.6 μ l	3.6 μ l
Water	28.8 μ l	28.8 μ l
Total	126 μ l	126 μ l

Example: For four sample determinations in quadruplicate: $4 \times 4 = 16$ + 2 extra (in case of pipetting errors) = 18; 18 samples X 10 μ l (reaction volume per well of a 384 well plate)= 180 μ l total volume; for 18 samples, 126 μ l is the total volume **before adding cDNA**.



- 1.22 Briefly spin the tubes, dispense 7 μ l of the master mix and 3 μ l cDNA to each well of a 384 well plate as assigned below:

	Sample 1	Sample 2	Sample 3	Sample 4
Gli1	B2, B3, B4, B5	B6, B7, B8, B9	B10, B11, B12, B13	B14, B15, B16, B17
GAPDH	C2, C3, C4, C5	C6, C7, C8, C9	C10, C11, C12, C13	C14, C15, C16, C17
NTC	D2, D3, D4, D5			

- 1.23 Seal the plate, centrifuge the plate briefly to collect the solutions to the bottom of the well. Set the PCR reaction using the fast cycling mode in a Vii7A real time PCR system:

Enzyme activation at 95°C, 2 min (1 cycle)

40 cycles of denaturation at 95°C, 1 sec and anneal/extend/acquire at 60°C , 20 sec

Dissociation curve conditions (melt curve stage): 95°C-15 sec, 60°C-1 min, and 95°C-15 sec.