

Jul 11, 2023

A new metabarcoding approach to survey diversity at the species level of Arcellinida (Amoebozoa: Tubulinea)

 [Molecular Ecology Resources](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm2389g3p/v1>

R. González-Miguéns¹

¹Real_Jardín_Botánico_CSIC



R. González-Miguéns

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm2389g3p/v1>

External link: <https://doi.org/10.1111/1755-0998.13771>

Protocol Citation: R. González-Miguéns 2023. A new metabarcoding approach to survey diversity at the species level of Arcellinida (Amoebozoa: Tubulinea). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm2389g3p/v1>

**Manuscript citation:**

González-Miguéns, R., Cano, E., Guillén-Oterino, A., Quesada, A., Lahr, D. J. G., Tenorio-Rodríguez, D., de Salvador-Velasco, D., Velázquez, D., Carrasco-Braganza, M. I., Patterson, R. T., Lara, E., & Singer, D. (2023). A needle in a haystack: A new metabarcoding approach to survey diversity at the species level of Arcellinida (Amoebozoa: Tubulinea). *Molecular Ecology Resources*, 23, 1034–1049. <https://doi.org/10.1111/1755-0998.13771>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 22, 2022

Last Modified: July 11, 2023

Protocol Integer ID: 73112

Keywords: diversity at the species level, species level, new metabarcoding approach, arcellinida, amoebozoa, pcr protocol, diversity, pcr protocol for the paper

Abstract

PCR protocol for the paper:" A needle in a haystack: a new metabarcoding approach to survey diversity at the species level of Arcellinida (Amoebozoa: Tubulinea)"



Polymerase Chain Reaction protocol

1

First PCR

Primers:

Primer forward: LCO 1490 (5' GGTCAACAAATCATAAAGATATTGG 3')

Primer reverse: HCO 2198 (5' TAACTTCAGGGTGACCAAAAAATCA 3')

Reagents for the PCR:

	A	B
	Polymerase	6 µL
	BSA	1 µL
	Primer forward	1 µL (10 µmol)
	Primer reverse	1 µL (10 µmol)
	Destilated water	1 µL
	Sample (eDNA)	1 µL

Total reaction mix per sample= 10 µL


Thermocycler program:

A	B	C	D
	Temperature	Time	cycles
Initial denaturation	96 °C	5 min	1
Denaturation	94 °C	15 s	40
Annealing	40 °C	15 s	
Extension	72 °C	90 s	
Final extension	72 °C	10 min	1

Denaturation temperatures and times may vary depending on the polymerase used. This protocol was tested with the MyTaq™ Red Mix™ polymerase.

Gel electrophoresis:

Analyze the results of your PCR reaction via gel electrophoresis. The resultant fragment was 692 bp long:



example of four PCR products with one negative control on the right. HyperLadder™ 1kb as the molecular weight marker.

2

Second PCR

Primers:

Primer forward: LCO 1490 (5' GGTCAACAAATCATAAAGATATTGG 3')

Primer reverse: ArCOIR (5' CCACYNGAATGWGCTARAATACC 3')

Reagents for the PCR:

	A	B
	Polymerase	12 µL



	A	B
	Primer forward	1 μ L (10 μ mol)
	Primer reverse	1 μ L (10 μ mol)
	Destilated water	6 μ L
	Sample (first PCR product)	1 μ L

Total reaction mix per sample= 20 μ L

Thermocycler program:

	A	B	C	D
		Temperatur e	Time	cycles
	Initial denaturatio n	96 °C	5 min	1
	Denaturatio n	94 °C	15 s	40
	Annealing	55 °C	15 s	
	Extension	72 °C	90 s	
	Final extension	72 °C	10 min	1

Denaturation temperatures and times may vary depending on the polymerase used. This protocol was tested with the MyTaq™ Red Mix™ polymerase.

Gel electrophoresis:

Analyze the results of your PCR reaction via gel electrophoresis. The resultant fragment was 407 bp long