


Oct 28, 2024

A Portable Device for Lab-Free, Versatile Nucleic Acid Extraction - Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxyj4wl8j/v1>



Anthony J Politza¹, Tianyi Liu¹, Aneesh Kshirsagar¹, Ming Dong¹, Md. Ahasan Ahamed¹, Weihua Guan¹

¹The Pennsylvania State University



Anthony J Politza

The Pennsylvania State University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxyj4wl8j/v1>

Protocol Citation: Anthony J Politza, Tianyi Liu, Aneesh Kshirsagar, Ming Dong, Md. Ahasan Ahamed, Weihua Guan 2024. A Portable Device for Lab-Free, Versatile Nucleic Acid Extraction - Protocol . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kxygxyj4wl8j/v1>

Manuscript citation:

Politza, AJ et. al. Development and Validation of a Portable Device for Lab-Free Versatile Nucleic Acid Extraction. *Biotechniques*, 2024.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

As of Oct. 2024, this protocol has been tested and validated for use with viral RNA from plasma or saliva, as well as for miRNA from saliva.

Created: October 14, 2024

Last Modified: October 28, 2024

Protocol Integer ID: 109832

Keywords: Sample preparation, Nucleic Acid Testing, Portable Sample Prep, Battery-powered, Field-deployable, Versatile RNA Extraction, easier extraction of viral rna, processing qiagen qiaamp rna, qiagen qiaamp rna, versatile nucleic acid extraction, viral rna, μ l of nuclease, rna, extraction concentration, extraction, easier extraction, extraction protocol, limiting extraction concentration, nuclease, portable centrifuge, benchtop extraction

Funders Acknowledgements:

National Science Foundation

Grant ID: 2319913

National Science Foundation

Grant ID: 2045169

National Institutes of Health

Grant ID: R33AI147419

National Institutes of Health

Grant ID: R33HD105610

USDA

Grant ID: NIFA 2022-11225

Abstract

To enable faster, easier extraction of viral RNA outside of traditional laboratories, we use a custom, automated, and portable centrifuge for processing Qiagen QIAamp RNA extractions.

Viral RNA is eluted into approx. 60 μ l of nuclease free water that can then be analyzed by downstream testing (PCR, gel electro, Qubit). Our results showed very strong agreement with benchtop extractions and probit analysis demonstrated a limiting extraction concentration of 0.75 cp/ μ l at 95% confidence.

Guidelines

This has now been tested on a variety of patient samples including HIV in plasma, SARS-CoV-2 in saliva, SARS-CoV-2 in VTM (from nasopharyngeal swab), and miRNA in saliva (as of October 15, 2024).



Materials

RNA EXTRACTION

- Qiagen - QIAamp Viral RNA Mini Kit

PRIMERS

- Primers were purchased from IDT.

RT-PCR

- Applied Biosystems - Fast Virus 1-Step Master Mix

Protocol materials



Applied Biosystems™ TaqMan™ Fast Virus 1-Step Master Mix for qPCR **Applied Biosystems (ThermoFisher Scientific) Catalog #44-444-32**



QIAamp Viral RNA Mini Kit **Qiagen Catalog #52904**

Safety warnings

 Please follow standard health and safety guidelines when working with viral samples.



Sample Collection

- 1 Label a  1.5 mL tube and add each  50 µl Sample .

Note

The following reagent volumes and steps are assuming a 50 µl starting sample.

Note

If you need to process more than 50 µl, increase the lysis buffer and ethanol proportionally to your new sample volume (550 µl : 50 µl).

A maximum of 750 µl will fit into the top of the mini-column, therefore if you increase sample volume you will need to add additional repetitions of the RNA/DNA binding steps.

****Please see the manufacturer's manual for additional details.****

- 2 Gather  QIAamp Viral RNA Mini Kit **Qiagen Catalog #52904**

RNA/DNA Lysis

- 3 Add  550 µL Lysis buffer to the  50 µl Sample

10m

Mix well by vortexing or pipetting.

Incubate for  00:10:00 at  Room temperature

Expected result

 600 µL DNA Lysate

Power On Portable Centrifuge



4 Turn on the power switch located on the rear of the device.

5

Expected result

A bright red LED should illuminate on the front of the device, indicating the device is ready to be used.

Note

Our device is pre-programmed to manage the spin down steps and automatically adjust the timings depending on the appropriate step. Example: Step 1. 60s Step 2. 60s etc.

At this time the use of other devices is not recommended or tested with this protocol. We imagine mini-centrifuges operating on external battery power could work as an alternative if they maintain higher than 1750g while spinning. (Can check with equation listed below)

$$RPM = \sqrt{\frac{RCF}{1.118 \times 10^{-5} \times r}}$$

RNA/DNA Binding

6 Add  550 µL Ethanol (95%) to the  600 µL DNA Lysate

Expected result

 1150 µL DNA Lysate

7 Mix well and pipette  750 µL DNA Lysate / 1150µl into the top of the QIAamp Mini Spin Column.



- 8 Place the entire QIAamp Mini Spin Column & tube into the portable centrifuge.
- 9 Close the lid of the device, check the red LED is illuminated, and push the front button beside the LED.

Note

If the red LED is not illuminated check that Step 4 was completed, and that the device was charged overnight.

- 10 Wait  00:01:00 for the liquid to spin down.

1m

- 11 Remove the QIAamp Mini Spin Column & discard the waste tube.
- 12 Insert the QIAamp Mini Spin Column into a fresh collection tube.
- 13 Repeat Steps 7 → 12.

RNA/DNA Wash 1

- 14 Add  500 µL Buffer AW1 to the top of QIAamp Mini Spin Column.
- 15 Place the entire QIAamp Mini Spin Column & tube into the portable centrifuge.
- 16 Close the lid of the device, check the red LED is illuminated, and push the front button beside the LED.
- 17 Wait  00:01:00 for the liquid to spin down.



18 Remove the QIAamp Mini Spin Column & discard the waste tube.

19 Insert the QIAamp Mini Spin Column into a fresh collection tube.

RNA/DNA Wash 2

20 Add  500 µL Buffer AW2 to the top of QIAamp Mini Spin Column.

21 Place the entire QIAamp Mini Spin Column & tube into the portable centrifuge.

22 Close the lid of the device, check the red LED is illuminated, and push the front button beside the LED.

23 Wait  00:01:00 for the liquid to spin down.

24 Remove the QIAamp Mini Spin Column & discard the waste tube.

25 Insert the QIAamp Mini Spin Column into a fresh collection tube.

RNA/DNA Drying

26 Add  500 µL Ethanol (95%) to the top of QIAamp Mini Spin Column.

27 Place the entire QIAamp Mini Spin Column & tube into the portable centrifuge.

28 Close the lid of the device, check the red LED is illuminated, and push the front button beside the LED.



29 Wait  00:03:00 for the liquid to spin down.

3m

30 Remove the QIAamp Mini Spin Column & discard the waste tube.

31 Insert the Mini Column into a fresh  1.5 mL tube and label accordingly.

RNA/DNA Elution

32 Add  80 μ L Nuclease Free Water to the top of QIAamp Mini Spin Column. Be sure to drop liquid into the center of the column membrane.

Note

In this protocol we used Nuclease Free Water to elute our samples for highest compatibility with downstream PCR. For other applications replace with Buffer AVE.

Note

In this protocol we used 80 μ l Nuclease Free Water to elute our samples. For other applications this volume could be adjusted. For highest yield a two-stage elution with 1/2 and 2/2 elution volume can be used, but was not tested in this protocol. Please see the Qiagen Viral RNA Mini Kit manual for more details when designing for your specific application.

33 Incubate  00:01:00 at  Room temperature

1m

34 Place the entire QIAamp Mini Spin Column & tube into the portable centrifuge.

35 Close the lid of the device, check the red LED is illuminated, and push the front button beside the LED.

36 Wait  00:01:00 for the liquid to spin down.

37 Collect the labeled  1.5 mL tube.

Expected result

The  1.5 mL tube now contains ~50-60 µl of RNA sample in Nuclease Free Water. Please store appropriately (or move immediately to downstream testing).

38 Discard the QIAamp Mini Spin Column.

RT-PCR

5m 53s

39 Thaw at  Room temperature and place on ice:

1. PCR Master Mix
2. PCR primers
3. Water
4. RNA Sample

The PCR Mix used in this protocol:

 Applied Biosystems™ TaqMan™ Fast Virus 1-Step Master Mix for qPCR **Applied Biosystems (ThermoFisher Scientific) Catalog #44-444-32**

Primers were purchased from Integrated DNA Technologies (IDT).

40 Setup the following reaction:

Component (Stock Concentration)	Volume
PCR Master Mix (4X)	 6.25 µL
Forward primer (10 µM)	 1.50 µL
Reverse primer (10 µM)	 1.50 µL
RNA Sample (variable)	 10 µL
Nuclease Free Water	 5.50 µL
Total	 25 µL



41 Set-up the following program on the thermal cycler:

5m 53s

Step	Temperature	Time	Cycles
Reverse Transcription	55 °C	00:05:00	1
Heat Activation	95 °C	00:00:20	1
Denaturation	95 °C	00:00:03	45-60
Annealing & Extension	60 °C	00:00:30	45-60
Hold	4 °C		

Protocol references

<https://www.qiagen.com/us/products/diagnostics-and-clinical-research/sample-processing/qiaamp-viral-rna-kits>

<https://www.idtdna.com/pages/products/custom-dna-rna/dna-oligos/custom-dna-oligos>

<https://doi.org/10.1128/jcm.41.10.4531-4536.2003>

<https://www.idtdna.com/pages/landing/coronavirus-research-reagents/cdc-assays>