



Apr 11, 2020

A protocol for massively parallel diagnosis and genome sequencing of SARS-CoV-2

DOI

<https://dx.doi.org/10.17504/protocols.io.betrjem6>

Leigh Monahan¹, Kay Anantanawat¹, Joyce To¹, Aaron Darling¹


¹University of Technology Sydney

Coronavirus Method De...



Leigh G Monahan

University of Technology Sydney



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.betrjem6>

Protocol Citation: Leigh Monahan, Kay Anantanawat, Joyce To, Aaron Darling 2020. A protocol for massively parallel diagnosis and genome sequencing of SARS-CoV-2. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.betrjem6>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

This protocol is currently in the early stages of testing. Feedback and suggestions are welcomed.

Created: April 08, 2020

Last Modified: April 11, 2020

Protocol Integer ID: 35409

Keywords: strategy for rna capture, genome sequencing of sar, cov genome, rna capture, viral rna, dna probe, rna extraction, diagnostic testing at an unprecedented scale, patient swab samples via dna, genome sequencing, stranded dna probe, sequencing data, many thousand diagnostic test, rna hybridisation, rna, sequencing protocol, sample reverse transcription, genome, whole genome, patient swab sample, parallel diagnosis, need for sar, reverse transcription, diagnostic testing, barcoded probe set, diagnostic test, stranded dna, entire sar, transcription

Abstract

Managing the current COVID-19 pandemic requires diagnostic testing at an unprecedented scale. However, the crisis has revealed severe deficiencies in our capacity to perform such testing. Here we build on the work of others (1-3) to develop a protocol that not only enables many thousand diagnostic tests to be run in parallel, but also provides near-whole genome sequencing data to facilitate phylogenetic analysis and contact tracing. One of the key features of our protocol is a magnetic bead-based strategy for RNA capture that may circumvent the need for SARS-CoV-2 RNA extraction, currently one of the major bottlenecks in both reagent supply and hands-on sample processing time. This approach also eliminates the requirement for per-sample reverse transcription, significantly reducing per-sample costs.

The major steps in our wet lab workflow can be summarised as follows. First, we generate a collection of bead-bound, single-stranded DNA probes tiling the entire SARS-CoV-2 genome. Multiple uniquely barcoded probe sets are prepared, each of which are used to capture viral RNA directly from patient swab samples via DNA/RNA hybridisation. Samples are combined together at this stage, enabling cDNA synthesis to be performed in a single pooled reaction. Finally, multiplex PCR is used to generate a library of overlapping amplicons ready for Illumina sequencing.

A schematic overview of the workflow is attached below, along with a more detailed figure depicting the various stages of library preparation.

References

1. https://docs.google.com/document/d/1kP2w_uTMSep2UxTCOnUhh1TMCjWvHEY0sUUpkJHPYV4/preview
2. <https://www.protocols.io/view/ncov-2019-sequencing-protocol-bbmuik6w/abstract>
3. <https://www.biorxiv.org/content/10.1101/2020.03.20.001008v1.full.pdf>

Attachments



[Workflow Overview.pdf...](#)

73KB



[Library Prep Schemat...](#)

79KB



Materials

Reagents

	ITEM	SUPPLIER	CATALOGUE NUMBER
		General consumables, chemicals and equipment	
	Ultra Pure 1M Tris-HCl pH 7.5	Invitrogen	15567027
	EDTA	Sigma-Aldrich	EDS-100G
	NaCl	Sigma-Aldrich	S3014-500G
	Ultra Pure DNase/RNase-free distilled water (500 mL)	Invitrogen	10977015
	Ethanol,	Sigma-	E7023-



	200 proof for molecular biology (500 mL)	Aldrich	500 mL
	DNA LoBind tubes, 1.5 mL (250 tubes)	Eppendorf	30108051
	0.2 mL MAX YMu m Recovery Thin Wall PCR Tubes	Fisher Biotec	PCR-02-L-C
	20X TE buffer - RNAse-free	Invitrogen	T11493
	Sodium citrate	Sigma-Aldrich	C8532-100G
	NaOH	Sigma-Aldrich	S8045-500G
	PCR plate - Twintec LoBind DNA (Epp	Sigma-Aldrich	EP0030129512-25EA

	end orf)		
	LiCl	Sig ma- Aldri ch	L96 50- 100 G
	LiDS	Sig ma- Aldri ch	L978 1-5G
	DTT	Sig ma- Aldri ch	D97 79- 250 MG
	PCR foil	Epp end orf	3012 779 0
		Step 1: Prep arati on of barc ode d bea ds for RNA capt ure	
	Dyn abea ds MyO ne Stre ptavi din C1	Invitr oge n	650 01
	Cust om DNA oligo s	IDT	-
	2X B&W buff er	See Buff er list	-
	1X SSC	See Buff	-



	buff er	er list	
	Prim eST AR GXL Poly mer ase	TaK aRa	R05 0A
		Step 2: Sam ple colle ctio n and prep arati on	
	Univ ersal Tran spor t Medi um (UT M)	COP AN	306 C
	Lysis /bin ding buff er	See Buff er list	-
	Qubi t HS RNA Assa y Kit	Invitr oge n	Q32 855
	Qubi t tube s	Invitr oge n	Q32 856
		Step 3: RNA capt ure	
	Was hing	See Buff	-



	buffer A	er list	
		Step 4: cDNA A synthesis	
	SuperScript™ IV First-Strand Synthesis System	Invitrogen	18091200
	Washing buffer B	See Buffer List	-
		Step 5: Library amplification	
	Q5 Hot Start HiFi DNA Polymerase	NEB	M0493S
	Custom DNA oligos	IDT	-
	Qubit dsDNA HS	Invitrogen	Q33230



	Assay Kit		
	PrimeSTAR GXL Polymerase	TaKaRa	R050A
	SPRI select beads	Beckman Coulter	B23318
		Step 6: Library quantitation and sequencing	
	Bioanalyzer HS dsDNA	Agilent	5067-4626
	Qubit dsDNA HS	Invitrogen	Q33230

Buffer list

	BUFFER NAME	COMPONENTS

	2X B&W buff er	10 mM Tris- HCl (pH 7.5)
		1 mM EDT A
		2 M NaCl
	1X SSC buff er	0.15 M NaCl
		0.01 5 M sodi um citra te
		Adju st pH to 7.0 with NaO H
	TE buff er	10 mM Tris- HCl (pH 7.5)
		1 mM EDT A
	Lysis /bin ding buff er	100 mM Tris- HCl (pH 7.5)
		500 mM

		LiCl
		10 mM EDTA, pH 8
		1% LiDS
		5 mM dithiothreitol (DTT)
	Washing buffer A	10 mM Tris-HCl (pH 7.5)
		150 mM LiCl
		1 mM EDTA
		0.1% LiDS
	Washing buffer B	10 mM Tris-HCl (pH 7.5)
		150 mM LiCl

**Individual primers [version 1]**

	NAM E	SEQ UEN CE
	[Biot in]- barc ode d bea d teth er_1	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC AGA TGT GTA TAA GAG ACA GGA ACG TAA TCG GCT CGG GTC GTC TG
	[Biot in]- barc ode d bea d teth er_2	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC AGA TGT GTA TAA GAG ACA GAT TGA



		TTA CGG GCT CGG GTC GTC TG
	[Biot in]- barc ode d bea d teth er_3	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC AGA TGT GTA TAA GAG ACA GTC CGT CAT CTG GCT CGG GTC GTC TG
	[Biot in]- barc ode d bea d teth er_4	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC AGA TGT GTA TAA GAG ACA GAC TAG TTG CGG



		GCT CGG GTC GTC TG
	[Biot in]- barc ode d bea d teth er_5	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC AGA TGT GTA TAA GAG ACA GAC CAT ATA AGG GCT CGG GTC GTC TG
	[Biot in]- barc ode d bea d teth er_6	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC AGA TGT GTA TAA GAG ACA GAT CTC GAT GAG GCT CGG



		GTC GTC TG
	[Biot in]- barc ode d bea d teth er_7	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC AGA TGT GTA TAA GAG ACA GAC TTA TAG TTG GCT CGG GTC GTC TG
	[Biot in]- bea d teth er	/5bi osg/ TCT AGA GCC ACC AGC GGC ATA GTA ATC GTC GGC AGC GTC
	Barc ode d linke r_1	CAG ACG ACC CGA GCC TAC CAT GCT TCT GTC TCT TAT



		ACA CAT CTG ACG CTG CCG ACG A
	Barcode linker_2	CAG ACG ACC CGA GCC AAT GGA GCC GCT GTC TCT TAT ACA CAT CTG ACG CTG CCG ACG A
	Barcode linker_3	CAG ACG ACC CGA GCC CTA ATC CGG TCT GTC TCT TAT ACA CAT CTG ACG CTG CCG ACG A
	Barcode linker_4	CAG ACG ACC CGA GCC TGC CTA ACT CCT GTC TCT TAT



		ACA CAT CTG ACG CTG CCG ACG A
	Barcode linker_5	CAG ACG ACC CGA GCC AAT AGT AGG CCT GTC TCT TAT ACA CAT CTG ACG CTG CCG ACG A
	Barcode linker_6	CAG ACG ACC CGA GCC GCA ACT TGG ACT GTC TCT TAT ACA CAT CTG ACG CTG CCG ACG A
	Barcode linker_7	CAG ACG ACC CGA GCC GAG ATC GTC ACT GTC TCT TAT



		ACA CAT CTG ACG CTG CCG ACG A
	Illumina P5 index primer_1	AAT GAT ACG GCG ACC ACC GAG ATC TAC ACT CGC TGA ATC GTC GGC AGC GTC
	Illumina P5 index primer_2	AAT GAT ACG GCG ACC ACC GAG ATC TAC ACC TCC AGG ATC GTC GGC AGC GTC
	Illumina P7 index primer_1	CAA GCA GAA GAC GGC ATA CGA GAT AAC CGC GGG TGA CTG GAG TTC AGA CGT GTG



		CTC TTC CGA TCT
	Illumina P7 index primer_2	CAA GCA GAA GAC GGC ATA CGA GAT GGT TAT AAG TGA CTG GAG TTC AGA CGT GTG CTC TTC CGA TCT
	Illumina P5 flow cell primer	AAT GAT ACG GCG ACC ACC GAG ATC TAC AC

SARS-CoV-2 probe pool 1 [version 1]

	NAM E	SEQ UEN CE
	ncov _cap _1_R1 GHT	CAG CAC ATC TAG GTT TCG TCC GCA GAC



		GAC CCG AGC C
	ncov _cap _3_R IGH T	CCT CAT GGT CAT GTT ATG GTT GAG CCA GAC GAC CCG AGC C
	ncov _cap _5_R IGH T	GTT ACC CGT GAA CTC ATG CGT GCA GAC GAC CCG AGC C
	ncov _cap _7_RI GHT	GCA GAC ACC TTT TGA AAT TAA ATT GGC CAG ACG ACC CGA GCC
	ncov _cap _9_R IGH T	GGT GAA ACT TCA TGG CAG ACG GCA GAC GAC CCG AGC C



	ncov _cap _11_R IGH T	GTA AGG GTG GTC GCA CTA TTG CCA GAC GAC CCG AGC C
	ncov _cap _13_ RIG HT	CTG CTT CCA CAA GTG CTT TTG TCA GAC GAC CCG AGC C
	ncov _cap _15_ RIG HT	CGT GTT TTA CAG AAG GCC GCT ACA GAC GAC CCG AGC C
	ncov _cap _17_R IGH T	AGA GAC GGT TGG GAA ATT GTT AAA TTT CAG ACG ACC CGA GCC
	ncov _cap _19_ RIG HT	AGA AGA AAC TGG CCT ACT



		CAT GCC AGA CGA CCC GAG CC
	ncov _cap _21_ RIG HT	TGG TAA CAA ACA ATA CCT TCA CAC TCA CAG ACG ACC CGA GCC
	ncov _cap _23_ RIG HT	ACC ACT GGG CAT TGA TTT AGA TGA CAG ACG ACC CGA GCC
	ncov _cap _25_ RIG HT	AGA GCA AGA AGA AGA TTG GTT AGA TGA CAG ACG ACC CGA GCC
	ncov _cap _27_ RIG HT	GTG GTT GTT AAT GCA GCC AAT GTC AGA CGA CCC



		GAG CC
	ncov _cap _29_ RIG HT	CAG CAC GAA GTT CTA CTT GCA CCA GAC GAC CCG AGC C
	ncov _cap _31_ RIG HT	ACC TTC AGT TGA ACA GAG AAA ACA AGC AGA CGA CCC GAG CC
	ncov _cap _33_ RIG HT	AAA GGC TGG TGG CAC TAC TGA ACA GAC GAC CCG AGC C
	ncov _cap _35_ RIG HT	CAG AAG AAA CAC GCA AAT TAA TGC CCA GAC GAC CCG AGC C
	ncov _cap	AAG TGC



	37 RIG HT	CAG CTA CAG TTT CTG TCA GAC GAC CCG AGC C
	ncov _cap _39_ RIG HT	GTT ATC ACC TTT GAC AAT CTT AAG ACA CTC AGA CGA CCC GAG CC
	ncov _cap _41_ RIG HT	TGA TCC TAG TTT TCT GGG TAG GTA CAC AGA CGA CCC GAG CC
	ncov _cap _43_ RIG HT	GAA GCT GCT AAC TTT TGT GCA CTC AGA CGA CCC GAG CC
	ncov _cap _45_ RIG HT	CGT GTG GTA AAC AAG CTA



		CAA AAT ATC TCA GAC GAC CCG AGC C
	ncov _cap _47_ RIG HT	CAA AGG TCC TAT TAC GGA TGT TTT CTA CCA GAC GAC CCG AGC C
	ncov _cap _49_ RIG HT	ACA TTT TTC CCT GAC TTA AAT GGT GAT GCA GAC GAC CCG AGC C
	ncov _cap _51_ RIG HT	TGA AGA AGT AGT GGA AAA TCC TAC CAT CAG ACG ACC CGA GCC
	ncov _cap _53_ RIG HT	TGC TGT TAA TAG TGT CCC



		TTG GGC AGA CGA CCC GAG CC
	ncov _cap _55_ RIG HT	CTC AAC CGC TGC TTT AGG TGT TCA GAC GAC CCG AGC C
	ncov _cap _57_ RIG HT	TGG TTT TTG GCA TATA TTC TTT TCA CTA GGC AGA CGA CCC GAG CC
	ncov _cap _59_ RIG HT	TCA ACT TGT ATG ATG TGT TAC AAA CGT AAC AGA CGA CCC GAG CC
	ncov _cap _61_ RIG HT	ACA GTG AAG AAT GGT TCC ATC CAT CAG



		ACG ACC CGA GCC
	ncov _cap _63_ RIG HT	TGT TAC TAG ATC AGG CAT TAG TGT CTG CAG ACG ACC CGA GCC
	ncov _cap _65_ RIG HT	AGG GTT TGT TGA TTC AGA TGT AGA AAC CAG ACG ACC CGA GCC
	ncov _cap _67_ RIG HT	ACG AAA ACA AAT ACG TAG TGC TGC CAG ACG ACC CGA GCC
	ncov _cap _69_ RIG HT	TGT CAT GTC TAA ACA TAC TGA CTT TTC AAG CAG ACG ACC



		CGA GCC
	ncov _cap _71_ RIG HT	ATT ACG CAC AAC TAA TGG TGA CTT TTC AGA CGA CCC GAG CC
	ncov _cap _73_ RIG HT	ACG TTA TGT GCT CAT GGA TGG CCA GAC GAC CCG AGC C
	ncov _cap _75_ RIG HT	GGT GCT TTG GAC ATA TCA GCA TCT CAG ACG ACC CGA GCC
	ncov _cap _77_ RIG HT	AGC ACA TAT TCA GTG GAT GGT TAT GTT CAG ACG ACC CGA GCC
	ncov _cap	TCT AAA



	79 RIG HT	GTT GCG TAG TGA TGT GCT ACA GAC GAC CCG AGC C
	ncov _cap _81_ RIG HT	CCC ATC TGG TAA AGT TGA GGG TTG CAG ACG ACC CGA GCC
	ncov _cap _83_ RIG HT	GCT TAA GGT TGA TAC AGC CAA TCC TCA GAC GAC CCG AGC C
	ncov _cap _85_ RIG HT	TCT TTT TGT TAC ATG CAC CAT ATG GAA TTC AGA CGA CCC GAG CC
	ncov _cap _87_ RIG HT	CCT CTA ACA CAA GAC



		CAT GTT GAC ACA GAC GAC CCG AGC C
	ncov _cap _89_ RIG HT	CAC ACC ACT GGT TGT TAC TCA CAC AGA CGA CCC GAG CC
	ncov _cap _91_ RIG HT	TGC TAG TTG GGT GAT GCG TAT TCA GAC GAC CCG AGC C
	ncov _cap _93_ RIG HT	TCT CTG TTA CTT CTA ACT ACT CAG GTG TAC AGA CGA CCC GAG CC
	ncov _cap _95_ RIG HT	ACT GAC TCT TGG TGT TTA TGA TTA CTT



		AGT CAG ACG ACC CGA GCC
	ncov _cap _97_ RIG HT	GGG CTC AAT GTG TCC AGT TAC ACA GAC GAC CCG AGC C
	ncov _cap _99_ RIG HT	TGC TAA TGG TGA TTC TGA AGT TGT TCT CAG ACG ACC CGA GCC
	ncov _cap _101_ RIG HT	GCA AGA GAT GGT TGT GTT CCC TCA GAC GAC CCG AGC C
	ncov _cap _103_ RIG HT	GCA TGG CCT CTT ATT GTA ACA GCC AGA CGA CCC GAG CC



	ncov _cap _105_ _RIG HT	ATG GAA CTG GTA CTA TCT ATA CAG AAC TGC AGA CGA CCC GAG CC
	ncov _cap _107_ _RIG HT	ACA CAC ACT GGT ACT GGT CAG GCA GAC GAC CCG AGC C
	ncov _cap _109_ _RIG HT	AAA CAC AGT CTG TAC CGT CTG CCA GAC GAC CCG AGC C
	ncov _cap _111_ _RIG HT	AGG ACG AAG ATG ACA ATT TAA TTG ATT CTT CAG ACG ACC CGA GCC
	ncov _cap _113_ _RIG HT	GCA TTT TGA



	RIG HT	TGA AGG TAA TTG TGA CAC ACA GAC GAC CCG AGC C
	ncov _cap _115_ RIG HT	GGT ATG ATT TCG GTG ATT TCA TAC AAA CCC AGA CGA CCC GAG CC
	ncov _cap _117_ RIG HT	TGT TAA CTG TTT GGA TGA CAG ATG CCA GAC GAC CCG AGC C
	ncov _cap _119_ RIG HT	TGT ATG CTG CTG ACC CTG CTA TCA GAC GAC CCG AGC C
	ncov _cap _121_ RIG HT	CCA ACA ATG TGT GAT



		ATC AGA CAA CTA CTC AGA CGA CCC GAG CC
	ncov _cap _123 _RIG HT	CGT AAT GTC ATC CCT ACT ATA ACT CAA ATG CAG ACG ACC CGA GCC
	ncov _cap _125 _RIG HT	TTA TGG CCT CAC TTG TTC TTG CTC AGA CGA CCC GAG CC
	ncov _cap _127_ RIG HT	AAT GCA CTT TTA TCT ACT GAT GGT AAC AAA CAG ACG ACC CGA GCC
	ncov _cap _129 _RIG HT	GGC TAG CAT AAA GAA CTT



		TAA GTC AGT TCA GAC GAC CCG AGC C
	ncov _cap _131_ RIG HT	GGT ACA CTT ATG ATT GAA CGG TTC GCA GAC GAC CCG AGC C
	ncov _cap _133_ _RIG HT	ACA CAC CGC ATA CAG TCT TAC AGC AGA CGA CCC GAG CC
	ncov _cap _135_ _RIG HT	TGT GCT AAT GGA CAA GTT TTT GGT TCA GAC GAC CCG AGC C
	ncov _cap _137_ RIG HT	TGC TGT CTG ACA GAG AAT TAC ATC TTC



		AGA CGA CCC GAG CC
	ncov _cap _139 _RIG HT	GCC ATT AAG TGC ACC TAC ACT AGT CAG ACG ACC CGA GCC
	ncov _cap _141_ RIG HT	CGC TGT TGA TGC ACT ATG TGA GAC AGA CGA CCC GAG CC
	ncov _cap _143 _RIG HT	CAA TGC CAG ATT ACG TGC TAA GCC AGA CGA CCC GAG CC
	ncov _cap _145 _RIG HT	AGG GTG TTA TCA CGC ATG ATG TCA GAC GAC CCG AGC C



	ncov _cap _147_ RIG HT	CAG CTC ACT CTT GTA ATG TAA ACA GAT TCA GAC GAC CCG AGC C
	ncov _cap _149_ _RIG HT	TCA AAA CTG AAG GTT TAT GTG TTG ACA CAG ACG ACC CGA GCC
	ncov _cap _151_ RIG HT	TAC CTT TAC AGC TAG GTT TTT CTA CAG GCA GAC GAC CCG AGC C
	ncov _cap _153_ _RIG HT	GGC TTT GAG TTG ACA TCT ATG AAG TAT TTT CAG ACG ACC CGA GCC



	ncov _cap _155_ _RIG HT	AGC TAG TTG TGA TGC AAT CAT GAC TCA GAC GAC CCG AGC C
	ncov _cap _157_ _RIG HT	GTA GAA TGG AAG TTC TAT GAT GCA CAG CAG ACG ACC CGA GCC
	ncov _cap _159_ _RIG HT	CAT TCC ACA CAC CAG CTT TTG ATA ACA GAC GAC CCG AGC C
	ncov _cap _161_ _RIG HT	AGC TTG TGG GTT TAC AAA CAA TTT GAC AGA CGA CCC GAG CC
	ncov _cap _163	AAT GTA GCA



	_RIG HT	TTT GAG CTT TGG GCC AGA CGA CCC GAG CC
	ncov _cap _165 _RIG HT	AGA AAT GCC CGT AAT GGT GTT CTC AGA CGA CCC GAG CC
	ncov _cap _167_ RIG HT	AGC TAT GGA TGA ATT CAT TGA ACG GTA CAG ACG ACC CGA GCC
	ncov _cap _169 _RIG HT	CCC AAG ATT TAT CTG TAG TTT CTA AGG TTG CAG ACG ACC CGA GCC
	ncov _cap _171_ RIG HT	GTC GCA AAA TAT ACT CAA CTG



		TGT CAC AGA CGA CCC GAG CC
	ncov _cap _173 _RIG HT	GAT CTC ATT ATT AGT GAT ATG TAC GAC CCT CAG ACG ACC CGA GCC
	ncov _cap _175 _RIG HT	AAT TGG ATG TAA TTA TCT TGG CAA ACC ACA GAC GAC CCG AGC C
	ncov _cap _177_ RIG HT	TGA TGT TCT TGT TAA CAA CTA AAC GAA CAC AGA CGA CCC GAG CC
	ncov _cap _179 _RIG HT	TGG GAC CAA TGG TAC TAA GAG



		GTT CAG ACG ACC CGA GCC
	ncov _cap _181_ RIG HT	CAA AAG TTG GAT GGA AAG TGA GTT CAC AGA CGA CCC GAG CC
	ncov _cap _183_ RIG HT	ACT AGG TTT CAA ACT TTA CTT GCT TTA CAC AGA CGA CCC GAG CC
	ncov _cap _185_ RIG HT	AGA GTC CAA CCA ACA GAA TCT ATT GTC AGA CGA CCC GAG CC
	ncov _cap _187_ RIG HT	GCA GAT TCA TTT GTA ATT AGA GGT GAT GAC



		AGA CGA CCC GAG CC
	ncov _cap _189 _RIG HT	ATC TAT CAG GCC GGT AGC ACA CCA GAC GAC CCG AGC C
	ncov _cap _191_ RIG HT	ACA AAA AGT TTC TGC CTT TCC AAC ACA GAC GAC CCG AGC C
	ncov _cap _193 _RIG HT	CAA ACA CGT GCA GGC TGT TTA ACA GAC GAC CCG AGC C
	ncov _cap _195 _RIG HT	TCT ACC AGT GTC TAT GAC CAA GAC ACA GAC GAC CCG AGC C



	ncov _cap _197_ RIG HT	CCA TCA AAA CCA AGC AAG AGG TCA CAG ACG ACC CGA GCC
	ncov _cap _199_ RIG HT	TCT GGT TGG ACC TTT GGT GCA GCA GAC GAC CCG AGC C
	ncov _cap _201_ RIG HT	ACG CTT GTT AAA CAA CTT AGC TCC CAG ACG ACC CGA GCC
	ncov _cap _203_ RIG HT	TGG AAA GGG CTA TCA TCT TAT GTC CTC AGA CGA CCC GAG CC
	ncov _cap _205_ RIG HT	GTC TGG TAA CTG TGA TGT



		TGT AAT AGG ACA GAC GAC CCG AGC C
	ncov _cap _207 _RIG HT	CTC ATC GAT CTC CAA GAA CTT GGA CAG ACG ACC CGA GCC
	ncov _cap _209 _RIG HT	GGA ACT GTA ACT TTG AAG CAA GGT GCA GAC GAC CCG AGC C
	ncov _cap _211 _RIG HT	AAC AGT TTA CTC ACA CCT TTT GCT CCA GAC GAC CCG AGC C
	ncov _cap _213 _RIG HT	TCA TTA CTT CAG GTG ATG GCA CAA CAG



		ACG ACC CGA GCC
	ncov _cap _215 _RIG HT	CCA ATT TAT GAT GAA CCG ACG ACG CAG ACG ACC CGA GCC
	ncov _cap _217 _RIG HT	AAA ACC TTC TTT TTA CGT TTA CTC TCG TCA GAC GAC CCG AGC C
	ncov _cap _219 _RIG HT	TGT CTT CTA CAA TTT GCC TAT GCC ACA GAC GAC CCG AGC C
	ncov _cap _221 _RIG HT	CTT CTC AAC GTG CCA CTC CAT GCA GAC GAC CCG AGC C



	ncov _cap _223 _RIG HT	GTC GCT ACA GGA TTG GCA ACT ATC AGA CGA CCC GAG CC
	ncov _cap _225 _RIG HT	TGA AGA GCA ACC AAT GGA GAT TGA CAG ACG ACC CGA GCC
	ncov _cap _227 _RIG HT	TGT TCA TCA GAC AAG AGG AAG TTC ACA GAC GAC CCG AGC C
	ncov _cap _229 _RIG HT	TTT TCT TAG GAA TCA TCA CAA CTG TAG CCA GAC GAC CCG AGC C
	ncov _cap _231	TGC CAG GAA CCT



	_RIG HT	AAA TTG GGT AGC AGA CGA CCC GAG CC
	ncov _cap _233 _RIG HT	GGC CCC AAG GTT TAC CCA ATA ATC AGA CGA CCC GAG CC
	ncov _cap _235 _RIG HT	TGG ACT TCC CTA TGG TGC TAA CAA CAG ACG ACC CGA GCC
	ncov _cap _237 _RIG HT	TTC AAC TCC AGG CAG CAG TAG GCA GAC GAC CCG AGC C
	ncov _cap _239 _RIG HT	GGT CCA GAA CAA ACC CAA GGA ACA GAC GAC CCG



		AGC C
	ncov _cap _241 _RIG HT	CGC ATA CAA AAC ATT CCC ACC AAC CAG ACG ACC CGA GCC
	ncov _cap _243 _RIG HT	CGT TTT CGC TTT TCC GTT TAC GAC AGA CGA CCC GAG CC
	ncov _cap _245 _RIG HT	TGG AAG AGC CCT AAT GTG TAA AAT TAA TTT CAG ACG ACC CGA GCC
	RNA se_P _bea d_re v	GTG GAG ACA GCC GCT CAC AGA CGA CCC GAG CC

SARS-CoV-2 probe pool 2 [version 1]



	NAM E	SEQ UEN CE
	ncov _cap _2_R IGH T	TCT TAA AGA TGG CAC TTG TGG CTT CAG ACG ACC CGA GCC
	ncov _cap _4_R IGH T	AGG AGC TGG TGG CCA TAG TTA CCA GAC GAC CCG AGC C
	ncov _cap _6_R IGH T	GCT TCA TGC ACT TTG TCC GAA CCA GAC GAC CCG AGC C
	ncov _cap _8_R IGH T	TGA TGG CTT TAT GGG TAG AAT TCG ACA GAC GAC CCG AGC C



	ncov _cap _10_ RIG HT	GCT GTT GTT AAA ATT TAT TGT CCA GCA TCA GAC GAC CCG AGC C
	ncov _cap _12_ RIG HT	TGG AGA AGG TTC CGA AGG TCT TCA GAC GAC CCG AGC C
	ncov _cap _14_ RIG HT	GGA ATA TTG GTG AAC AGA AAT CAA TAC TGA CAG ACG ACC CGA GCC
	ncov _cap _16_ RIG HT	AGG TGG TGT TGT TCA GTT GAC TCA GAC GAC CCG AGC C
	ncov _cap _18_ RIG HT	TTT GGC TTT



	RIG HT	GTG TGC TGA CTC TCA GAC GAC CCG AGC C
	ncov _cap _20_ RIG HT	TTA GAA CAA CCT ACT AGT GAA GCT GTT CAG ACG ACC CGA GCC
	ncov _cap _22_ RIG HT	GCC TAT ACA GTT GAA CTC GGT ACA GCA GAC GAC CCG AGC C
	ncov _cap _24_ RIG HT	TGT GAA GAA GAA GAG TTT GAG CCA CAG ACG ACC CGA GCC
	ncov _cap _26_ RIG HT	GGA ACT TAC ACC AGT TGT TCA GAC



		TCA GAC GAC CCG AGC C
	ncov _cap _28_ RIG HT	CTT AAA GTG GGT GGT AGT TGT GTT TTC AGA CGA CCC GAG CC
	ncov _cap _30_ RIG HT	CAC AAA TGT CTA CTT AGC TGT CTT TGA TCA GAC GAC CCG AGC C
	ncov _cap _32_ RIG HT	CAA TCT TCA TCC AGA TTC TGC CAC TCA GAC GAC CCG AGC C
	ncov _cap _34_ RIG HT	CAG TGC TTA AAA AGT GTA AAA GTG CCC AGA



		CGA CCC GAG CC
	ncov _cap _36_ RIG HT	GTA AAA CAA CTG TAG CGT CAC TTA TCA CAG ACG ACC CGA GCC
	ncov _cap _38_ RIG HT	GTT CCT ATA AAG ATT GGT CCT ATT CTG GAC AGA CGA CCC GAG CC
	ncov _cap _40_ RIG HT	ACT TAT TTG GAT GGA GCT GAT GTT ACT AAC AGA CGA CCC GAG CC
	ncov _cap _42_ RIG HT	CTT CTA TTA AAT GGG CAG ATA ACA ACT GTC AGA



		CGA CCC GAG CC
	ncov _cap _44_ RIG HT	AAG AGT CTT GAA CGT GGT GTG TCA GAC GAC CCG AGC C
	ncov _cap _46_ RIG HT	TGG TAC ATT TAC TTG TGC TAG TGA GTA CAG ACG ACC CGA GCC
	ncov _cap _48_ RIG HT	AAC CAA CCA TAT CCA AAC GCA AGC AGA CGA CCC GAG CC
	ncov _cap _50_ RIG HT	CCT GGT GTA TAC GTT GTC TTT GGA CAG ACG ACC CGA GCC



	ncov _cap _52_ RIG HT	CAC ACA GAT CTA ATG GCT GCT TAT GTC AGA CGA CCC GAG CC
	ncov _cap _54_ RIG HT	GCA AAG AAT ACT GTT AAG AGT GTC GGT ACA GAC GAC CCG AGC C
	ncov _cap _56_ RIG HT	GTA GTG TTT GTC TTA GTG GTT TAG ATT CTT CAG ACG ACC CGA GCC
	ncov _cap _58_ RIG HT	ATG GCC CCG ATT TCA GCT ATG GCA GAC GAC CCG AGC C
	ncov _cap	CTG TGC



	60 RIG HT	TGG TAG TAC ATT TAT TAG TGA TGC AGA CGA CCC GAG CC
	ncov _cap _62_ RIG HT	AGC TAA TAA CAC TAA AGG TTC ATT GCC TAC AGA CGA CCC GAG CC
	ncov _cap _64_ RIG HT	ACT CAA AAC ACT AGT TGC AAC TGC CAG ACG ACC CGA GCC
	ncov _cap _66_ RIG HT	GCT TGT ATT GAC TGT AGT GCG CCA GAC GAC CCG AGC C
	ncov _cap _68_ RIG HT	AGC ACT TAA GGG TGG



		TAA AAT TGT TAA TCA GAC GAC CCG AGC C
	ncov _cap _70_ RIG HT	GCT GAT TTT GAC ACA TGG TTT AGC CCA GAC GAC CCG AGC C
	ncov _cap _72_ RIG HT	GTG TTT TGG CTG CTG AAT GTA CAC AGA CGA CCC GAG CC
	ncov _cap _74_ RIG HT	GCT GGT GTT TGT GTA TCT ACT AGT GGC AGA CGA CCC GAG CC
	ncov _cap _76_ RIG HT	GAG CTT TTG GTG AAT ACA GTC ATG



		TAG CAG ACG ACC CGA GCC
	ncov _cap _78_ RIG HT	TCC ACA AAG CAT TTC TAT TGG TTC TTT CAG ACG ACC CGA GCC
	ncov _cap _80_ RIG HT	GCT TGT TGT CAT CTC GCA AAG GCA GAC GAC CCG AGC C
	ncov _cap _82_ RIG HT	TGA AGA CAT GCT TAA CCC TAA TTA TGA AGA CAG ACG ACC CGA GCC
	ncov _cap _84_ RIG HT	TGC TAT GAG GCC CAA TTT CAC TCA GAC GAC CCG



		AGC C
	ncov _cap _86_ RIG HT	TTT TAG CTT GGT TGT ACG CTG CCA GAC GAC CCG AGC C
	ncov _cap _88_ RIG HT	CCA TAT TGG GTA GTG CTT TAT TAG AAG ATG CAG ACG ACC CGA GCC
	ncov _cap _90_ RIG HT	TGT CTG CTT TTG CAA TGA TGT TTG TCA GAC GAC CCG AGC C
	ncov _cap _92_ RIG HT	ACT GTG TAT GAT GAT GGT GCT AGG CAG ACG ACC CGA GCC



	ncov _cap _94_ RIG HT	TGC CCT ATT TTC TTC ATA ACT GGT AAT ACA CAG ACG ACC CGA GCC
	ncov _cap _96_ RIG HT	GGC AAA CCT TGT ATC AAA GTA GCC CAG ACG ACC CGA GCC
	ncov _cap _98_ RIG HT	GCT TTG TGA AGA AAT GCT GGA CAC AGA CGA CCC GAG CC
	ncov _cap _100_ RIG HT	TGA CCC AAA TGT ATA AAC AGG CTA GAC AGA CGA CCC GAG CC
	ncov _cap _102	AAA TAC GTG TGA



	_RIG HT	TGG TAC AAC ATT TAC TTC AGA CGA CCC GAG CC
	ncov _cap _104 _RIG HT	ACA AAC TGC TTG CAC TGA TGA CAC AGA CGA CCC GAG CC
	ncov _cap _106 _RIG HT	CTT GGT AGT TTA GCT GCC ACA GTC AGA CGA CCC GAG CC
	ncov _cap _108 _RIG HT	TTT GGT GGT GCA TCG TGT TGT CCA GAC GAC CCG AGC C
	ncov _cap _110_ RIG HT	GTA TAC AGG GCT TTT GAC ATC TAC AAT



		GCA GAC GAC CCG AGC C
	ncov _cap _112_ RIG HT	AAT AGA CGG TGA CAT GGT ACC ACC AGA CGA CCC GAG CC
	ncov _cap _114_ RIG HT	CGC CAA CTT AGG TGA ACG TGT ACA GAC GAC CCG AGC C
	ncov _cap _116_ RIG HT	TGT TGA CAC TGA CTT AAC AAA GCC TTC AGA CGA CCC GAG CC
	ncov _cap _118_ RIG HT	GTG AGA AAA ATA TTT GTT GAT GGT GTT CCC AGA CGA CCC



		GAG CC
	ncov _cap _120 _RIG HT	GCT TTT CAA ACT GTC AAA CCC GGC AGA CGA CCC GAG CC
	ncov _cap _122 _RIG HT	CCT AGA CAA ATC AGC TGG TTT TCC ACA GAC GAC CCG AGC C
	ncov _cap _124 _RIG HT	TTG GAA CAA GCA AAT TCT ATG GTG GCA GAC GAC CCG AGC C
	ncov _cap _126 _RIG HT	TCA CTA TAT GTT AAA CCA GGT GGA ACC CAG ACG ACC CGA GCC



	ncov _cap _128 _RIG HT	AGA CTT TGT GAA TGA GTT TTA CGC ATC AGA CGA CCC GAG CC
	ncov _cap _130 _RIG HT	CTC TCA ACA TAC AAT GCT AGT TAA ACA GGC AGA CGA CCC GAG CC
	ncov _cap _132 _RIG HT	AAG CTA CAT GAT GAG TTA ACA GGA CAC AGA CGA CCC GAG CC
	ncov _cap _134 _RIG HT	TCT TGT CTG TTA ATC CGT ATG TTT GCC AGA CGA CCC GAG CC
	ncov _cap	TTT AGC



	_136 _RIG HT	TAA CAC CTG TAC TGA AAG ACT CAG ACG ACC CGA GCC
	ncov _cap _138 _RIG HT	CAA ATA GGA GAG TAC ACC TTT GAA AAA GGC AGA CGA CCC GAG CC
	ncov _cap _140 _RIG HT	AAG GTT GGT ATG CAA AAG TAT TCT ACA CCA GAC GAC CCG AGC C
	ncov _cap _142 _RIG HT	GTA CTG TAA ATG CAT TGC CTG AGA CCA GAC GAC CCG AGC C
	ncov _cap _144	TAG GTC CAG



	_RIG HT	ACA TGT TCC TCG GCA GAC GAC CCG AGC C
	ncov _cap _146 _RIG HT	GCC TCA AAG ATT TTG GGA CTA CCA CAG ACG ACC CGA GCC
	ncov _cap _148 _RIG HT	CGT AGG AAT GTG GCA ACT TTA CAA GCA GAC GAC CCG AGC C
	ncov _cap _150 _RIG HT	CCC GCG AAG AAG CTA TAA GAC ACA GAC GAC CCG AGC C
	ncov _cap _152 _RIG HT	CCT CAT ACC ACT TAT GTA CAA AGG ACT



		CAG ACG ACC CGA GCC
	ncov _cap _154 _RIG HT	GCA TCA TTC TAT TGG ATT TGA TTA CGT CTC AGA CGA CCC GAG CC
	ncov _cap _156 _RIG HT	CTT GTA GAA AGG TTC AAC ACA TGG TTC AGA CGA CCC GAG CC
	ncov _cap _158 _RIG HT	TGC AAT GTC GAT AGA TAT CCT GCT AAC AGA CGA CCC GAG CC
	ncov _cap _160 _RIG HT	CCA CTA AAG TCT GCT ACG TGT ATA ACA CAG ACG



		ACC CGA GCC
	ncov _cap _162 _RIG HT	TGA TGG ACA ACA GGG TGA AGT ACC AGA CGA CCC GAG CC
	ncov _cap _164 _RIG HT	CTC CAG CAC ATAT ATC TAC TAT TGG TGT CAG ACG ACC CGA GCC
	ncov _cap _166 _RIG HT	AGA AGC CGT AAA AAC ACA GTT CAA TTC AGA CGA CCC GAG CC
	ncov _cap _168 _RIG HT	TGG TTT ACA TCT ACT GAT TGG ACT AGC CAG ACG ACC CGA GCC



	ncov _cap _170 _RIG HT	AAT CTA GTC AAG CGT GGC AAC CCA GAC GAC CCG AGC C
	ncov _cap _172 _RIG HT	TTG CAC CAG GTA CAG CTG TTT TCA GAC GAC CCG AGC C
	ncov _cap _174 _RIG HT	TGG AGG TTC CGT GGC TAT AAA GAT CAG ACG ACC CGA GCC
	ncov _cap _176 _RIG HT	CAA ATC CAA TTC AGT TGT CTT CCT ATT CTC AGA CGA CCC GAG CC



	ncov _cap _178 _RIG HT	CAG TGT GTT AAT CTT ACA ACC AGA ACT CCA GAC GAC CCG AGC C
	ncov _cap _180 _RIG HT	TTC GAA GAC CCA GTC CCT ACT TCA GAC GAC CCG AGC C
	ncov _cap _182 _RIG HT	GCC TAT TAA TTT AGT GCG TGA TCT CCC AGA CGA CCC GAG CC
	ncov _cap _184 _RIG HT	TGA AAA TGG AAC CAT TAC AGA TGC TGC AGA CGA CCC GAG CC
	ncov _cap _186	ACT GTG TTG



	_RIG HT	CTG ATT ATT CTG TCC TCA GAC GAC CCG AGC C
	ncov _cap _188 _RIG HT	CGT TAT AGC TTG GAA TTC TAA CAA TCT TGA CAG ACG ACC CGA GCC
	ncov _cap _190 _RIG HT	TTC TTT TGA ACT TCT ACA TGC ACC AGC AGA CGA CCC GAG CC
	ncov _cap _192 _RIG HT	AAC AAA TAC TTC TAA CCA GGT TGC TGC AGA CGA CCC GAG CC
	ncov _cap _194 _RIG HT	ACG TAG TGT AGC TAG



		TCA ATC CAT CAC AGA CGA CCC GAG CC
	ncov _cap _196 _RIG HT	ACT GGA ATA GCT GTT GAA CAA GAC ACA GAC GAC CCG AGC C
	ncov _cap _198 _RIG HT	ATT GCT GCT AGA GAC CTC ATT TGT CAG ACG ACC CGA GCC
	ncov _cap _200 _RIG HT	GAA CCA AAA ATT GAT TGC CAA CCA ATC AGA CGA CCC GAG CC
	ncov _cap _202 _RIG HT	AAG TTT GCA GAC ATA TGT GAC TCA ACC



		AGA CGA CCC GAG CC
	ncov _cap _204 _RIG HT	TGG AAA AGC ACA CTT TCC TCG TCA GAC GAC CCG AGC C
	ncov _cap _206 _RIG HT	TTG ATT TAG GTG ACA TCT CTG GCA CAG ACG ACC CGA GCC
	ncov _cap _208 _RIG HT	GCT GTA TGA CCA GTT GCT GTA GTC AGA CGA CCC GAG CC
	ncov _cap _210 _RIG HT	CTT ATT GTT GGC GTT GCA CTT CTT CAG ACG ACC CGA GCC



	ncov _cap _212 _RIG HT	GCC GTT CCA AAA ACC CAT TAC TTC AGA CGA CCC GAG CC
	ncov _cap _214 _RIG HT	TGT ATT ACA CAG TTA CTT CAC TTC AGA CTC AGA CGA CCC GAG CC
	ncov _cap _216 _RIG HT	TTC TTG CTT TCG TGG TAT TCT TGC CAG ACG ACC CGA GCC
	ncov _cap _218 _RIG HT	CAT GGC AGA TTC CAA CGG TAC TCA GAC GAC CCG AGC C



	ncov _cap _220 _RIG HT	GGA TCA CCG GTG GAA TTG CTA TCA GAC GAC CCG AGC C
	ncov _cap _222 _RIG HT	ACG CTG TGA CAT CAA GGA CCT GCA GAC GAC CCG AGC C
	ncov _cap _224 _RIG HT	GTT GAC TTT CAG GTT ACT ATA GCA GAG ACA GAC GAC CCG AGC C
	ncov _cap _226 _RIG HT	ATT CAC CAT TTC ATC CTC TAG CTG ATC AGA CGA CCC GAG CC
	ncov _cap _228	ATT TGT GCT TTT



	_RIG HT	TAG CCT TTC TGC TAC AGA CGA CCC GAG CC
	ncov _cap _230 _RIG HT	CGT GTC CTA TTC ACT TCT ATT CTA AAT GGC AGA CGA CCC GAG CC
	ncov _cap _232 _RIG HT	CCC AAA ATC AGC GAA ATG CAC CCA GAC GAC CCG AGC C
	ncov _cap _234 _RIG HT	GCA GTC CAG ATG ACC AAA TTG GCC AGA CGA CCC GAG CC
	ncov _cap _236 _RIG HT	CCT CAA GGA ACA ACA TTG CCA AAC



		AGA CGA CCC GAG CC
	ncov _cap _238 _RIG HT	ACA ACA ACA AGG CCA AAC TGT CCA GAC GAC CCG AGC C
	ncov _cap _240 _RIG HT	ATT GGC ATG GAA GTC ACA CCT TCA GAC GAC CCG AGC C
	ncov _cap _242 _RIG HT	TTC TTC CTG CTG CAG ATT TGG ACA GAC GAC CCG AGC C
	ncov _cap _244 _RIG HT	AGT GTG TAA CAT TAG GGA GGA CTT GCA GAC GAC CCG AGC C



	RNA se_P _bea d_re v	GTG GAG ACA GCC GCT CAC AGA CGA CCC GAG CC
--	----------------------------------	--

SARS-CoV-2 primer pool 1 [version 1]

	NAM E	SEQ UEN CE
	ncov _cap _1_L EFT	ACG TGT GCT CTT CCG ATC TTC TCT TGT AGA TCT GTT CTC TAA ACG AAC
	ncov _cap _3_L EFT	ACG TGT GCT CTT CCG ATC TCT GTT TTA CAG GTT CGC GAC GT
	ncov _cap _5_L EFT	ACG TGT GCT CTT CCG ATC TCC TCA TGT GGG CGA



		AAT ACC AG
	ncov _cap _7_L EFT	ACG TGT GCT CTT CCG ATC TAT GGC TAC CCT CTT GAG TGC AT
	ncov _cap _9_L EFT	ACG TGT GCT CTT CCG ATC TGG GGA ATG TCC AAA TTT TGT ATT TCC C
	ncov _cap _11_L EFT	ACG TGT GCT CTT CCG ATC TCC ACT TGC GAA TTT TGT GGC AC
	ncov _cap _13_ LEF T	ACG TGT GCT CTT CCG ATC TCT ATT GGG TTC CAC GTG



		CTA GC
	ncov _cap _15_ LEF T	ACG TGT GCT CTT CCG ATC TTC AAA CAA ATT GTT GAA TCC TGT GGT
	ncov _cap _17_L EFT	ACG TGT GCT CTT CCG ATC TTG TTC ACA TCT GAT TTG GCT ACT AAC A
	ncov _cap _19_ LEF T	ACG TGT GCT CTT CCG ATC TCC TGT GCA AAG GAA ATT AAG GAG AGT
	ncov _cap _21_ LEF T	ACG TGT GCT CTT CCG ATC TTC CCA CAG AAG TGT



		TAA CAG AGG A
	ncov _cap _23_ LEF T	ACG TGT GCT CTT CCG ATC TAG TGC AAG GTT ACA AGA GTG TGA
	ncov _cap _25_ LEF T	ACG TGT GCT CTT CCG ATC TGG CTT CAC ATA TGT ATT GTT CTT TCT ACC
	ncov _cap _27_ LEF T	ACG TGT GCT CTT CCG ATC TAG GAC AAT CAG ACA ACT ACT ATT CAA ACA
	ncov _cap _29_ LEF T	ACG TGT GCT CTT CCG ATC TGG CTA CTA



		ACA ATG CCA TGC AAG
	ncov _cap _31_ LEF T	ACG TGT GCT CTT CCG ATC TTA AGA GTT TGT GTA GAT ACT GTT CGC A
	ncov _cap _33_ LEF T	ACG TGT GCT CTT CCG ATC TTC TGG AAG AAA CTA AGT TCC TCA CAG
	ncov _cap _35_ LEF T	ACG TGT GCT CTT CCG ATC TTT ATAT AAC CAC TTA CCC GGG TCA GG
	ncov _cap _37_ LEF T	ACG TGT GCT CTT CCG ATC TAC AAG



		AGG GTG TGG TTG ATT ATG G
	ncov _cap _39_ LEF T	ACG TGT GCT CTT CCG ATC TAC TTC TTC TTC TAA AAC ACC TGA AGA ACA
	ncov _cap _41_ LEF T	ACG TGT GCT CTT CCG ATC TCG CAA GTT GTG GAC ATG TCA ATG
	ncov _cap _43_ LEF T	ACG TGT GCT CTT CCG ATC TTG TCA GCA TTA AAT CAC ACT AAA AAG TGG
	ncov _cap _45_ LEF T	ACG TGT GCT CTT CCG ATC



		TAA TGA GTT ACT TGT TTC AAC ATG CCA
	ncov _cap _47_ LEF T	ACG TGT GCT CTT CCG ATC TAC CTT TTG TTA TGA TGT CAG CAC CA
	ncov _cap _49_ LEF T	ACG TGT GCT CTT CCG ATC TAG AAA GAC AAT TCT TAT TTC ACA GAG CAA
	ncov _cap _51_ LEF T	ACG TGT GCT CTT CCG ATC TTG TTT GGC ATG TTA ACA ATG CAA CT
	ncov _cap _53_ LEF T	ACG TGT GCT CTT



	LEF T	CCG ATC TAC CGA AGT TGT AGG AGA CAT TAT ACT TAA
	ncov _cap _55_ LEF T	ACG TGT GCT CTT CCG ATC TAA TTC TAG AAT TAA AGC ATC TAT GCC GAC
	ncov _cap _57_ LEF T	ACG TGT GCT CTT CCG ATC TTT GAA CTC TAC TAA TGT CAC TAT TGC AAC C
	ncov _cap _59_ LEF T	ACG TGT GCT CTT CCG ATC TGC AAT TGT TTT TCA GCT ATT TTG



		CAG T
	ncov _cap _61_ LEF T	ACG TGT GCT CTT CCG ATC TTG TCT ATG CTA ATG GAG GTA AAG GC
	ncov _cap _63_ LEF T	ACG TGT GCT CTT CCG ATC TAA GAC TTA TGA AAG ACA TTC TCT CTC TCA
	ncov _cap _65_ LEF T	ACG TGT GCT CTT CCG ATC TTA GTG CGG AAG TTG CAG TTA AAA TG
	ncov _cap _67_ LEF T	ACG TGT GCT CTT CCG ATC TGC GAT AGT TGT AAT



		AAC TAT ATG CTC ACC
	ncov _cap _69_ LEF T	ACG TGT GCT CTT CCG ATC TAC TTA CCT TTT AAG TTG ACA TGT GCA A
	ncov _cap _71_ LEF T	ACG TGT GCT CTT CCG ATC TTC ACT CGT GAC ATA GCA TCT ACA GA
	ncov _cap _73_ LEF T	ACG TGT GCT CTT CCG ATC TAA CAT CTG TTA CAC ACC ATC AAA ACT T
	ncov _cap _75_ LEF T	ACG TGT GCT CTT CCG ATC TTT



		CTG AGT ACT GTA GGC ACG GC
	ncov _cap _77_ LEF T	ACG TGT GCT CTT CCG ATC TAT GAG GTT TAG AAG AGC TTT TGG TGA
	ncov _cap _79_ LEF T	ACG TGT GCT CTT CCG ATC TTG TTC ACA CCT TTA GTA CCT TTC TGG
	ncov _cap _81_ LEF T	ACG TGT GCT CTT CCG ATC TAG TAC AAG TAT TTT AGT GGA GCA ATG GA
	ncov _cap _83_ LEF T	ACG TGT GCT CTT CCG ATC



		TAA CGG TCT TTG GCT TGA TGA CG
	ncov _cap _85_ LEF T	ACG TGT GCT CTT CCG ATC TTA AGT TTG TTC GCA TTC AAC CAG GA
	ncov _cap _87_ LEF T	ACG TGT GCT CTT CCG ATC TTG GAC CTT TTG TTG ACA GGC AA
	ncov _cap _89_ LEF T	ACG TGT GCT CTT CCG ATC TGG AAT TGC CGT TTT AGA TAT GTG TGC
	ncov _cap _91_ LEF T	ACG TGT GCT CTT CCG ATC TGT



		CCA GAG TAC TCA ATG GTC TTT GT
	ncov _cap _93_ LEF T	ACG TGT GCT CTT CCG ATC TAC TGT GTT ATG TAT GCA TCA GCT GT
	ncov _cap _95_ LEF T	ACG TGT GCT CTT CCG ATC TAC AAC TGT CAT GTT TTT GGC CAG
	ncov _cap _97_ LEF T	ACG TGT GCT CTT CCG ATC TAC TAC TCC CAC CCA AGA ATA GCA
	ncov _cap _99_ LEF T	ACG TGT GCT CTT CCG ATC TTG GTT



		TCA CTA CTT TCT GTT TTG CTT
	ncov _cap _101_ LEF T	ACG TGT GCT CTT CCG ATC TAA TTT GAC CGT GAT GCA GCC AT
	ncov _cap _103_ LEF T	ACG TGT GCT CTT CCG ATC TAC CTC TTA CAA CAG CAG CCA AAC
	ncov _cap _105_ LEF T	ACG TGT GCT CTT CCG ATC TGC TTA GTC CTG TTG CAC TAC GAC
	ncov _cap _107_ LEF T	ACG TGT GCT CTT CCG ATC TGG TAT GGT ACT



		TGG TAG TTT AGC TGC
	ncov _cap _109_ LEF T	ACG TGT GCT CTT CCG ATC TTC AGG CAA TAA CAG TTA CAC CGG
	ncov _cap _111_ LEF T	ACG TGT GCT CTT CCG ATC TAG TCA GCT GAT GCA CAA TCG TT
	ncov _cap _113_ LEF T	ACG TGT GCT CTT CCG ATC TAC ACT TTC TCT AAC TAC CAA CAT GAA GA
	ncov _cap _115_ LEF T	ACG TGT GCT CTT CCG ATC TAA GGA CTG GTA



		TGA TTT TGT AGA AAA CCC
	ncov _cap _117_ LEF T	ACG TGT GCT CTT CCG ATC TTG TTA ATG CCT ATA TTA ACC TTG ACC AGG
	ncov _cap _119_ LEF T	ACG TGT GCT CTT CCG ATC TTG TTT TAT TCT CTA CAG TGT TCC CAC C
	ncov _cap _121_ LEF T	ACG TGT GCT CTT CCG ATC TTG TTG CTT TTC AAA CTG TCA AAC CC
	ncov _cap _123_ LEF T	ACG TGT GCT CTT CCG ATC



		TGA TAA GTA CTT TGA TTG TTA CGA TGG TGG
	ncov _cap _125_ _LEF T	ACG TGT GCT CTT CCG ATC TTC AAA AAT TAT TGA AAT CAA TAG CCG CCA
	ncov _cap _127_ _LEF T	ACG TGT GCT CTT CCG ATC TTG TAG CTT GTC ACA CCG TTT CT
	ncov _cap _129_ _LEF T	ACG TGT GCT CTT CCG ATC TTG TCC GCA ATT TAC AAC ACA GAC T
	ncov _cap _131_ T	ACG TGT GCT CTT



	LEF T	CCG ATC TAT GTC TGA AGC AAA ATG TTG GAC TG
	ncov _cap _133_ LEF T	ACG TGT GCT CTT CCG ATC TAT GCT TAC CCA CTT ACT AAA CAT CCT
	ncov _cap _135_ LEF T	ACG TGT GCT CTT CCG ATC TGT TGT AAA TGC TGT TAC GAC CAT GT
	ncov _cap _137_ LEF T	ACG TGT GCT CTT CCG ATC TAG CGA TAA TGT TAC TGA CTT TAA TGC AA
	ncov _cap	ACG TGT



	_139 _LEF T	GCT CTT CCG ATC TCC ACT TAA CCG AAA TTA TGT CTT TAC TGG
	ncov _cap _141_ _LEF T	ACG TGT GCT CTT CCG ATC TAC CCA ACA CTC AAT ATC TCA GAT GAG
	ncov _cap _143_ _LEF T	ACG TGT GCT CTT CCG ATC TAA TGT AGT AGA ATT ATA CCT GCA CGT GC
	ncov _cap _145_ _LEF T	ACG TGT GCT CTT CCG ATC TAC CAG AAT ATT TCA ATT CAG TGT



		GTA GAC
	ncov _cap _147_ LEF T	ACG TGT GCT CTT CCG ATC TGG TAA GAG AAT TCC TTA CAC GTA ACC C
	ncov _cap _149_ LEF T	ACG TGT GCT CTT CCG ATC TCT GAT AGA GAC CTT TAT GAC AAG TTG CA
	ncov _cap _151_ LEF T	ACG TGT GCT CTT CCG ATC TAG GAC ATG ACC TAT AGA AGA CTC ATC T
	ncov _cap _153_ LEF T	ACG TGT GCT CTT CCG ATC TTT TTC CAG



		AGT TAG TGC TAA ACC ACC
	ncov _cap _155_ LEF T	ACG TGT GCT CTT CCG ATC TTA GAC GTG CCA CAT GCT TTT CC
	ncov _cap _157_ LEF T	ACG TGT GCT CTT CCG ATC TTG TTA AGC GTG TTG ACT GGA CT
	ncov _cap _159_ LEF T	ACG TGT GCT CTT CCG ATC TCT ATT CTT ATG CCA CAC ATT CTG ACA AA
	ncov _cap _161_ LEF T	ACG TGT GCT CTT CCG ATC TTC TGA CAG



		TCC ATG TGA GTC TCA
	ncov _cap _163 _LEF T	ACG TGT GCT CTT CCG ATC TTG GAA CAC TTT TAC AAG ACT TCA GAG T
	ncov _cap _165 _LEF T	ACG TGT GCT CTT CCG ATC TGG GTG TGG ACA TTG CTG CTA AT
	ncov _cap _167_ _LEF T	ACG TGT GCT CTT CCG ATC TAA CAA GCT AGT CTT AAT GGA GTC ACA
	ncov _cap _169 _LEF T	ACG TGT GCT CTT CCG ATC TTC TAC TGA



		TTG GAC TAG CTA AAC GTT T
	ncov _cap _171_ LEF T	ACG TGT GCT CTT CCG ATC TTG GTG TAA AGA TGG CCA TGT AGA A
	ncov _cap _173_ LEF T	ACG TGT GCT CTT CCG ATC TCA TTT TGG TGC TGG TTC TGA TAA AGG
	ncov _cap _175_ LEF T	ACG TGT GCT CTT CCG ATC TTT CAC TTA CAT TTG TGG GTT TAT ACA ACA
	ncov _cap _177_ LEF T	ACG TGT GCT CTT CCG ATC



		TAC AAA TCC AAT TCA GTT GTC TTC CTA TTC
	ncov _cap _179 _LEF T	ACG TGT GCT CTT CCG ATC TAG TCA GTG TGT TAA TCT TAC AAC CAG A
	ncov _cap _181_ LEF T	ACG TGT GCT CTT CCG ATC TTT TGC TTC CAC TGA GAA GTC TAA CA
	ncov _cap _183_ LEF T	ACG TGT GCT CTT CCG ATC TAC CTT GAA GGA AAA CAG GGT AAT TTC A
	ncov _cap	ACG TGT



	_185 _LEF T	GCT CTT CCG ATC TTG CTG CAG CTT ATT ATG TGG GT
	ncov _cap _187_ LEF T	ACG TGT GCT CTT CCG ATC TAG TTT TTA ACG CCA CCA GAT TTG C
	ncov _cap _189_ LEF T	ACG TGT GCT CTT CCG ATC TGG CAA ACT GGA AAG ATT GCT GAT
	ncov _cap _191_ LEF T	ACG TGT GCT CTT CCG ATC TAC CCA CTA ATG GTG TTG GTT ACC
	ncov _cap _193	ACG TGT GCT CTT



	LEF T	CCG ATC TTC TTG ACA TTA CAC CAT GTT CTT TTG G
	ncov _cap _195 _LEF T	ACG TGT GCT CTT CCG ATC TAG GTA TAT GCG CTA GTT ATC AGA CTC A
	ncov _cap _197 _LEF T	ACG TGT GCT CTT CCG ATC TCA ATC TTT TGT TGC AAT ATG GCA GTT T
	ncov _cap _199 _LEF T	ACG TGT GCT CTT CCG ATC TAA GTG ACA CTT GCA GAT GCT GG



	ncov _cap _201 _LEF T	ACG TGT GCT CTT CCG ATC TTG GCT TAT AGG TTT AAT GGT ATT GGA GT
	ncov _cap _203 _LEF T	ACG TGT GCT CTT CCG ATC TAC AAA GTT GAG GCT GAA GTG CA
	ncov _cap _205 _LEF T	ACG TGT GCT CTT CCG ATC TAC TTA TGT CCC TGC ACA AGA AAA GA



	ncov _cap _207 _LEF T	ACG TGT GCT CTT CCG ATC TTG ATC CTT TGC AAC CTG AAT TAG ACT
	ncov _cap _209 _LEF T	ACG TGT GCT CTT CCG ATC TTG GTG ACA ATT ATG CTT TGC TGT ATG
	ncov _cap _211 _LEF T	ACG TGT GCT CTT CCG ATC TGT TCG CGC TAC TGC AAC GAT AC
	ncov _cap _213 _LEF T	ACG TGT GCT CTT CCG ATC TAT GCT TTA GTC TAC TTC TTG CAG AGT



	ncov _cap _215 _LEF T	ACG TGT GCT CTT CCG ATC TAC ACA GTT ACT TCA CTT CAG ACT ATT ACC
	ncov _cap _217 _LEF T	ACG TGT GCT CTT CCG ATC TAC AAG CTG ATG AGT ACG AAC TTA TGT
	ncov _cap _219 _LEF T	ACG TGT GCT CTT CCG ATC TTC TAG AGT TCC TGA TCT TCT GGT CT
	ncov _cap _221 _LEF T	ACG TGT GCT CTT CCG ATC TGG CCA GTA ACT TTA GCT TGT



		TTT GT
	ncov _cap _223 _LEF T	ACG TGT GCT CTT CCG ATC TGT AAT CGG AGC TGT GAT CCT TCG
	ncov _cap _225 _LEF T	ACG TGT GCT CTT CCG ATC TTG CTT GTA CAG TAA GTG ACA ACA GA
	ncov _cap _227 _LEF T	ACG TGT GCT CTT CCG ATC TAC AAC AGT ACT TTT AAA AGA ACC TTG CT
	ncov _cap _229 _LEF T	ACG TGT GCT CTT CCG ATC TTG CTT CAC ACT CAA AAG



		AAA GAC AGA
	ncov _cap _231 _LEF T	ACG TGT GCT CTT CCG ATC TTG TAG TTG ATG ACC CGT GTC CT
	ncov _cap _233 _LEF T	ACG TGT GCT CTT CCG ATC TAG AGT ATC ATG ACG TTC GTG TTG T
	ncov _cap _235 _LEF T	ACG TGT GCT CTT CCG ATC TAG GAA GAC CTT AAA TTC CCT CGA GG
	ncov _cap _237 _LEF T	ACG TGT GCT CTT CCG ATC TTG AGG GAG CCT TGA ATA



		CAC CA
	ncov _cap _239 _LEF T	ACG TGT GCT CTT CCG ATC TCT CTT GCT TTG CTG CTG CTT GA
	ncov _cap _241 _LEF T	ACG TGT GCT CTT CCG ATC TCA TTG GCC GCA AAT TGC ACA AT
	ncov _cap _243 _LEF T	ACG TGT GCT CTT CCG ATC TAG GCT GAT GAA ACT CAA GCC TT
	ncov _cap _245 _LEF T	ACG TGT GCT CTT CCG ATC TTG AAT TCT CGT AAC TAC ATA GCA



		CAA GT
	RNA se_P _LEF T	ACG TGT GCT CTT CCG ATC TGA TTT GGA CCT GCG AGC G

SARS-CoV-2 primer pool 2 [version 1]

	NAM E	SEQ UEN CE
	ncov _cap _2_L EFT	ACG TGT GCT CTT CCG ATC TTC AGC ACA TCT AGG TTT CGT CC
	ncov _cap _4_L EFT	ACG TGT GCT CTT CCG ATC TCA TCA AAC GTT CGG ATG CTC GA
	ncov _cap _6_L EFT	ACG TGT GCT CTT CCG ATC TTG GCA



		CTG ATC CTT ATG AAG ATT TTC A
	ncov _cap _8_L EFT	ACG TGT GCT CTT CCG ATC TGC ATG AAA TTG CTT GGT ACA CGG
	ncov _cap _10_ LEF T	ACG TGT GCT CTT CCG ATC TGC GTC ACC AAA TGA ATG CAA CC
	ncov _cap _12_ LEF T	ACG TGT GCT CTT CCG ATC TAC CTG AGC ATA GTC TTG CCG AA
	ncov _cap _14_ LEF T	ACG TGT GCT CTT CCG ATC TGT GAC TTT



		AAA CTT AAT GAA GAG ATC GCC
	ncov _cap _16_ LEF T	ACG TGT GCT CTT CCG ATC TTG TAC GAT CAA TTT TCT CCC GCA
	ncov _cap _18_ LEF T	ACG TGT GCT CTT CCG ATC TAA ACC CGT CCT TGA TTG GCT TG
	ncov _cap _20_ LEF T	ACG TGT GCT CTT CCG ATC TTT TGT CAC GCA CTC AAA GGG AT
	ncov _cap _22_ LEF T	ACG TGT GCT CTT CCG ATC TGC CCT TGC ACC



		TAA TAT GAT GGT
	ncov _cap _24_ LEF T	ACG TGT GCT CTT CCG ATC TTG CTG TCA TAA AAA CTT TGC AAC CA
	ncov _cap _26_ LEF T	ACG TGT GCT CTT CCG ATC TGG TAA ACC TTT GGA ATT TGG TGC C
	ncov _cap _28_ LEF T	ACG TGT GCT CTT CCG ATC TAC ATT AAA AAT GCA GAC ATT GTG GAA GA
	ncov _cap _30_ LEF T	ACG TGT GCT CTT CCG ATC TTG TCT TCA



		TGT TGT CGG CCC AA
	ncov _cap _32_ LEF T	ACG TGT GCT CTT CCG ATC TAA GAT CGC TGA GAT TCC TAA AGA GG
	ncov _cap _34_ LEF T	ACG TGT GCT CTT CCG ATC TGT GGG TGA TGT TGT TCA AGA GGG
	ncov _cap _36_ LEF T	ACG TGT GCT CTT CCG ATC TTG GAA CTG TTT CTT GGA ATT TGC G
	ncov _cap _38_ LEF T	ACG TGT GCT CTT CCG ATC TTG TAA CAC ATG



		GCT TAA ATT TGG AAG AA
	ncov _cap _40_ LEF T	ACG TGT GCT CTT CCG ATC TTC CTA CCA CAT TCC ACC TAG ATG G
	ncov _cap _42_ LEF T	ACG TGT GCT CTT CCG ATC TAA CCT CAT AAT TCA CAT GAA GGT AAA ACA
	ncov _cap _44_ LEF T	ACG TGT GCT CTT CCG ATC TAG AGT TGA AGT TTA ATC CAC CTG CT
	ncov _cap _46_ LEF T	ACG TGT GCT CTT CCG ATC TCC



		TTA AGG GTG TAG AAG CTG TTA TGT
	ncov _cap _48_ LEF T	ACG TGT GCT CTT CCG ATC TCC TCA GAA TAC AAA GGT CCT ATT ACG G
	ncov _cap _50_ LEF T	ACG TGT GCT CTT CCG ATC TAG AAA CCT GCT TCA AGA GAG CTT
	ncov _cap _52_ LEF T	ACG TGT GCT CTT CCG ATC TGG GAA TGG ATA ATC TTG CCT GCG
	ncov _cap _54_ LEF T	ACG TGT GCT CTT CCG ATC TTG



		CTA AGC CTT TTC TTA ACA AAG TTG T
	ncov _cap _56_ LEF T	ACG TGT GCT CTT CCG ATC TTG GTT TTT ACT ATT AAG TGT TTG CCT AGG T
	ncov _cap _58_ LEF T	ACG TGT GCT CTT CCG ATC TAA TGG GAT TTA ACT GCT TTT GGC TT
	ncov _cap _60_ LEF T	ACG TGT GCT CTT CCG ATC TGC ATG TTG TAG ACG GTT GTA ATT CA
	ncov _cap _62_ T	ACG TGT GCT



	LEF T	CTT CCG ATC TGA CTT GTC ACT ACA GTT TAA AAG ACC AA
	ncov _cap _64_ LEF T	ACG TGT GCT CTT CCG ATC TAG AAT CAT CTG CAA AAT CAG CGT CT
	ncov _cap _66_ LEF T	ACG TGT GCT CTT CCG ATC TCT ACT TTT ATT TCA GCA GCT CGG C
	ncov _cap _68_ LEF T	ACG TGT GCT CTT CCG ATC TGC GCA GGT AGC AAA AAG TCA CA
	ncov _cap	ACG TGT



	70 LEF T	GCT CTT CCG ATC TGC AGT TAA TTA AAG TTA CAC TTG TGT TCC
	ncov _cap _72_ LEF T	ACG TGT GCT CTT CCG ATC TTG CTG CAG TCA TAA CAA GAG AAG T
	ncov _cap _74_ LEF T	ACG TGT GCT CTT CCG ATC TAC CAA TGT ACT AGA AGG TTC TGT TGC
	ncov _cap _76_ LEF T	ACG TGT GCT CTT CCG ATC TGA TCT TTA CCA GGA GTT TTC TGT GGT



	ncov _cap _78_ LEF T	ACG TGT GCT CTT CCG ATC TAC TGT ACT CTG TTT AAC ACC AGT TTA CT
	ncov _cap _80_ LEF T	ACG TGT GCT CTT CCG ATC TCC TTT AGT ACT TTT GAA GAA GCT GCG
	ncov _cap _82_ LEF T	ACG TGT GCT CTT CCG ATC TCC ACA AAC CTC TAT CAC CTC AGC
	ncov _cap _84_ LEF T	ACG TGT GCT CTT CCG ATC TGG CTG GTA ATG TTC AAC TCA GGG



	ncov _cap _86_ LEF T	ACG TGT GCT CTT CCG ATC TTG GTT CAT GTG GTA GTG TTG GTT
	ncov _cap _88_ LEF T	ACG TGT GCT CTT CCG ATC TAC CAC AAC TCT TAA TGA CTT TAA CCT TGT
	ncov _cap _90_ LEF T	ACG TGT GCT CTT CCG ATC TAC AAT GCT CAG GTG TTA CTT TCC A



	ncov _cap _92_ LEF T	ACG TGT GCT CTT CCG ATC TTT ACC TTC TCT TGC CAC TGT AGC T
	ncov _cap _94_ LEF T	ACG TGT GCT CTT CCG ATC TAC ACT TAT GAA TGT CTT GAC ACT CGT
	ncov _cap _96_ LEF T	ACG TGT GCT CTT CCG ATC TTG TAC TTG TTA CTT TGG CCT CTT TTG T
	ncov _cap _98_ LEF T	ACG TGT GCT CTT CCG ATC TAG TCT TAC TCT CAG TTT TGC AAC



		AAC T
	ncov _cap _100 _LEF T	ACG TGT GCT CTT CCG ATC TGT TCC CTT CCA TCA TAT GCA GCT
	ncov _cap _102 _LEF T	ACG TGT GCT CTT CCG ATC TGG GCA AAA GTT ACT AGT GCT ATG CA
	ncov _cap _104 _LEF T	ACG TGT GCT CTT CCG ATC TCA GGT TGT AGA TGC AGA TAG TAA AAT TGT
	ncov _cap _106 _LEF T	ACG TGT GCT CTT CCG ATC TTG CAC TGT TAT CCG ATT



		TAC AGG A
	ncov _cap _108_ LEF T	ACG TGT GCT CTT CCG ATC TAA GTG CCT GCC AAT TCA ACT GT
	ncov _cap _110_ LEF T	ACG TGT GCT CTT CCG ATC TAA AAA CAC AGT CTG TAC CGT CTG C
	ncov _cap _112_ LEF T	ACG TGT GCT CTT CCG ATC TAA TTC CTA AAA ACT AAT TGT TGT CGC TTC C
	ncov _cap _114_ LEF T	ACG TGT GCT CTT CCG ATC TTC AAC GTC TTA



		CTA AAT ACA CAA TGG CA
	ncov _cap _116_ LEF T	ACG TGT GCT CTT CCG ATC TCG AAA TGC TGG TAT TGT TGG TGT
	ncov _cap _118_ LEF T	ACG TGT GCT CTT CCG ATC TTG ACT TCA CGG AAG AGA GGT TAA AAC
	ncov _cap _120_ LEF T	ACG TGT GCT CTT CCG ATC TCC ACT TCA GAG AGC TAG GTG TTG
	ncov _cap _122_ LEF T	ACG TGT GCT CTT CCG ATC TCA CTT CTT CTT



		TGC TCA GGA TGG T
	ncov _cap _124 _LEF T	ACG TGT GCT CTT CCG ATC TCG CAT ATA CAA AAC GTA ATG TCA TCC C
	ncov _cap _126 _LEF T	ACG TGT GCT CTT CCG ATC TAT GGG TTG GGA TTA TCC TAA ATG TGA
	ncov _cap _128 _LEF T	ACG TGT GCT CTT CCG ATC TTC AGG AGA TGC CAC AAC TGC TT
	ncov _cap _130 _LEF T	ACG TGT GCT CTT CCG ATC TTA CTC TCT



		GAC GAT GCT GTT GTG
	ncov _cap _132 _LEF T	ACG TGT GCT CTT CCG ATC TAG ATC CAT CAA GAA TCC TAG GGG C
	ncov _cap _134 _LEF T	ACG TGT GCT CTT CCG ATC TGG AAC CTG AGT TTT ATG AGG CTA TGT
	ncov _cap _136 _LEF T	ACG TGT GCT CTT CCG ATC TAC TCA ACT TTA CTT AGG AGG TAT GAG CT
	ncov _cap _138 _LEF T	ACG TGT GCT CTT CCG ATC TAC GCT



		CAA AGC TAC TGA GGA GA
	ncov _cap _140 _LEF T	ACG TGT GCT CTT CCG ATC TAC CGA GGT ACA ACA ACT TAC AAA TTA AAT
	ncov _cap _142 _LEF T	ACG TGT GCT CTT CCG ATC TCT CTA CTA CCC TTC TGC TCG CA
	ncov _cap _144 _LEF T	ACG TGT GCT CTT CCG ATC TCT TTG ATG AAA TTT CAA TGG CCA CAA A
	ncov _cap _146 _LEF T	ACG TGT GCT CTT CCG ATC TGC



		ACA TAA AGA CAA ATC AGC TCA ATG C
	ncov _cap _148 _LEF T	ACG TGT GCT CTT CCG ATC TGT TGA TTC ATC ACA GGG CTC AGA
	ncov _cap _150 _LEF T	ACG TGT GCT CTT CCG ATC TTG GGT TAC ATC CTA CAC AGG CA
	ncov _cap _152 _LEF T	ACG TGT GCT CTT CCG ATC TGT CGA GGG GTG TCA TGC TAC TA
	ncov _cap _154 _LEF T	ACG TGT GCT CTT CCG ATC TAG TGA



		CAC ACT TAA AAA TCT CTC TGA CA
	ncov _cap _156 _LEF T	ACG TGT GCT CTT CCG ATC TCA GGT AAC CTA CAA AGC AAC CAT GA
	ncov _cap _158 _LEF T	ACG TGT GCT CTT CCG ATC TCC CAG TTC TTC ACG ACA TTG GT
	ncov _cap _160 _LEF T	ACG TGT GCT CTT CCG ATC TAA CTT GCC TGG TTG TGA TGG TG
	ncov _cap _162 _LEF T	ACG TGT GCT CTT CCG ATC TTG TAG



		ACA TCA TGC TAA TGA GTA CAG AT
	ncov _cap _164 _LEF T	ACG TGT GCT CTT CCG ATC TCA CAA AAG TTG ATG GTG TTG ATG TAG A
	ncov _cap _166 _LEF T	ACG TGT GCT CTT CCG ATC TAC TGA AAC GAT TTG TGC ACC ACT
	ncov _cap _168 _LEF T	ACG TGT GCT CTT CCG ATC TCC AAC AAT TAC CTG AAA CTT ACT TTA CTC A
	ncov _cap _170	ACG TGT GCT CTT



	LEF T	CCG ATC TGG TTC ATC TAA GTG TGT GTG TTC TGT
	ncov _cap _172 _LEF T	ACG TGT GCT CTT CCG ATC TAG AAT GCT ATT AGA AAA GTG TGA CCT TCA
	ncov _cap _174 _LEF T	ACG TGT GCT CTT CCG ATC TTG ATG CAG ATT CAA CTT TGA TTG GTG
	ncov _cap _176 _LEF T	ACG TGT GCT CTT CCG ATC TGC TGA TCT TTA TAA GCT CAT GGG ACA
	ncov _cap	ACG TGT



	_178 _LEF T	GCT CTT CCG ATC TAG GGG TAC TGC TGT TAT GTC TTT AAA
	ncov _cap _180 _LEF T	ACG TGT GCT CTT CCG ATC TAG GAC TTG TTC TTA CCT TTC TTT TCC A
	ncov _cap _182 _LEF T	ACG TGT GCT CTT CCG ATC TGT TGG ATG GAA AGT GAG TTC AGA GT
	ncov _cap _184 _LEF T	ACG TGT GCT CTT CCG ATC TTC GGC TTT AGA ACC ATT GGT AGA TT



	ncov _cap _186 _LEF T	ACG TGT GCT CTT CCG ATC TAG TGT ACG TTG AAA TCC TTC ACT GT
	ncov _cap _188 _LEF T	ACG TGT GCT CTT CCG ATC TAC TTT TAA GTG TTA TGG AGT GTC TCC T
	ncov _cap _190 _LEF T	ACG TGT GCT CTT CCG ATC TAG ATT GTT TAG GAA GTC TAA TCT CAA ACC T
	ncov _cap _192 _LEF T	ACG TGT GCT CTT CCG ATC TAT GGT TTA ACA GGC ACA



		GGT GT
	ncov _cap _194 _LEF T	ACG TGT GCT CTT CCG ATC TAT CAA CTT ACT CCT ACT TGG CGT G
	ncov _cap _196 _LEF T	ACG TGT GCT CTT CCG ATC TAA CTC TAT TGC CAT ACC CAC AAA TTT T
	ncov _cap _198 _LEF T	ACG TGT GCT CTT CCG ATC TAG TCA AAC AAA TTT ACA AAA CAC CAC CA
	ncov _cap _200 _LEF T	ACG TGT GCT CTT CCG ATC TTT GCC ACC TTT



		GCT CAC AGA TG
	ncov _cap _202 _LEF T	ACG TGT GCT CTT CCG ATC TGC ACT TGG AAA ACT TCA AGA TGT GG
	ncov _cap _204 _LEF T	ACG TGT GCT CTT CCG ATC TTG CTG CTA CTA AAA TGT CAG AGT GT
	ncov _cap _206 _LEF T	ACG TGT GCT CTT CCG ATC TTG TAA CAC AAA GGA ATT TTT ATG AAC CAC A
	ncov _cap _208 _LEF T	ACG TGT GCT CTT CCG ATC TGT AAA



		CAT TCA AAA AGA AAT TGA CCG CC
	ncov _cap _210 _LEF T	ACG TGT GCT CTT CCG ATC TGC TCA AAG GAG TCA AAT TAC ATT ACA CA
	ncov _cap _212 _LEF T	ACG TGT GCT CTT CCG ATC TCC AAG GGT GTT CAC TTT GTT TGC
	ncov _cap _214 _LEF T	ACG TGT GCT CTT CCG ATC TGC CAA CTA TTT TCT TTG CTG GCA
	ncov _cap _216 _LEF T	ACG TGT GCT CTT CCG ATC TTT



		CAC ACA ATC GAC GGT TCA TCC
	ncov _cap _218 _LEF T	ACG TGT GCT CTT CCG ATC TCC TTA CTG CGC TTC GAT TGT GT
	ncov _cap _220 _LEF T	ACG TGT GCT CTT CCG ATC TAA GCT CCT TGA ACA ATG GAA CCT
	ncov _cap _222 _LEF T	ACG TGT GCT CTT CCG ATC TTT CTT TCA GAC TGT TTG CGC GT
	ncov _cap _224 _LEF T	ACG TGT GCT CTT CCG ATC TAA CGC TTT CTT



		ATT ACA AAT TGG GAG C
	ncov _cap _226 _LEF T	ACG TGT GCT CTT CCG ATC TTC TCA ATT AGA TGA AGA GCA ACC AAT GG
	ncov _cap _228 _LEF T	ACG TGT GCT CTT CCG ATC TCG TCT ATC AGT TAC GTG CCA GAT
	ncov _cap _230 _LEF T	ACG TGT GCT CTT CCG ATC TAT TAT CTT TTG GTT CTC ACT TGA ACT GC
	ncov _cap _232 _LEF T	ACG TGT GCT CTT CCG ATC TAC



		ATC GAT ATC GGT AAT TAT ACA GTT TCC T
	ncov _cap _234 _LEF T	ACG TGT GCT CTT CCG ATC TAC CCT CAG ATT CAA CTG GCA GT
	ncov _cap _236 _LEF T	ACG TGT GCT CTT CCG ATC TAA TGA AAG ATC TCA GTC CAA GAT GGT
	ncov _cap _238 _LEF T	ACG TGT GCT CTT CCG ATC TGC AGT CAA GCC TCT TCT CGT TC
	ncov _cap _240 _LEF T	ACG TGT GCT CTT CCG ATC



		TGC AAA AAC GTA CTG CCA CTA AAG C
	ncov _cap _242 _LEF T	ACG TGT GCT CTT CCG ATC TAC ACA GGT GCC ATC AAA TTG GA
	ncov _cap _244 _LEF T	ACG TGT GCT CTT CCG ATC TTG CTG ACT CAA CTC AGG CCT AA
	RNA se_P _LEF T	ACG TGT GCT CTT CCG ATC TGA TTT GGA CCT GCG AGC G

Preparation of barcoded beads for RNA capture

1

Note

The first stage in our workflow is the preparation of bead-bound DNA probes that will facilitate subsequent viral RNA capture, cDNA synthesis and multiplex PCR to generate ~200 nt amplicons covering the entire SARS-CoV-2 genome. We use two distinct sets of probes for this, requiring two independent bead preparations for each sample to be processed. Ultimately, this enables library preparation to be performed in two parallel reactions, designed such that neighbouring amplicons do not overlap within the same library PCR. This strategy protects against preferential amplification of short overlap products, and has been successfully developed and employed elsewhere (Quick *et al.* [2017] *Nat Protoc.* 12: 1261; <https://www.protocols.io/view/ncov-2019-sequencing-protocol-bbmui6w/abstract>). Throughout the protocol we use the generic terms “pool 1” and “pool 2” to refer to our two capture probe sets, as well as the corresponding primer sets needed for library amplification.

We envisage that capture beads could be prepared in bulk in advance and supplied to testing labs in 96-well or 384-well format. Each uniquely barcoded bead preparation could thus be used for multiple sample batches, provided each batch is sequenced independently. The protocol below is for a small-scale preparation of up to 10 uniquely barcoded bead sets, to be used in development and evaluation of a prototype. We also provide an alternative, lower-cost bead preparation method at the end of this page.

Wash  200 μ L of streptavidin-coated Dynabeads™ MyOne™ C1 magnetic beads as per the manufacturer's instructions:

- 1.1 Resuspend the Dynabeads™ magnetic beads in the vial (vortex for >  00:00:30).
- 1.2 Transfer  200 μ L of Dynabeads™ magnetic beads to a tube.
- 1.3 Add an equal volume of 1X B&W buffer and resuspend.
- 1.4 Place the tube on a magnet for  00:01:00 and discard the supernatant.
- 1.5 Remove the tube from the magnet and resuspend the washed magnetic beads in  200 μ L of 1X B&W buffer.

- 1.6 Repeat steps 1.4–1.5 twice, for a total of 3 washes.
- 1.7 Pellet beads once more, and resuspend in  400 µL of 2X B&W buffer.
- 2 For each unique barcode, take an aliquot of the prepared beads and immobilise an oligonucleotide of the structure 5'-[biotin][spacer][partial Illumina p5 adapter][sample barcode][generic stub]-3' onto the beads as follows:
 - 2.1 Add  40 µL of 5 µM oligonucleotide ([Biotin]-barcoded bead tether) to  40 µL of washed beads.
 - 2.2 Incubate for  00:15:00 at room temperature using gentle rotation.
 - 2.3 Separate the biotinylated DNA coated beads with a magnet for 2–3 min.
 - 2.4 Wash the coated beads 2–3 times with  100 µL 1X B&W buffer.
 - 2.5 Resuspend in  120 µL of nuclease-free water.
- 3 For each barcoded sample, generate two distinct libraries of bead-bound SARS-CoV-2 probes by annealing separate oligonucleotide pools with the general structure 3'-[generic stub complement][SARS-CoV-2 target]-5' and performing strand extension:
 - 3.1 Mix the following components in 0.2 mL PCR tubes:

Component	Pool 1	Pool 2
Bead-bound template (previous step)	 60 µL	 60 µL
5X PrimeSTAR GXL buffer	 20 µL	 20 µL
dNTP mixture (2.5 mM each)	 8 µL	 8 µL
SARS-CoV-2 probe pool 1 or 2 (10 µM)	 10 µL	 10 µL

PrimeSTAR GXL DNA polymerase

 2 μL  2 μL

3.2 Incubate the reaction as follows:

Step	Temperature	Time	Cycles
Denaturation	 98 °C	 00:00:10	1
Annealing	 55 °C	 00:00:15	1
Extension	 68 °C	 00:02:00	1
Hold	 4 °C	Indefinite	1

4 For each bead sample, dissociate and remove the non-biotinylated strand as follows:

4.1 Pellet the beads on a magnet and wash once in  50 μL 1X SSC.4.2 Resuspend the beads in  20 μL of freshly prepared 0.15 M NaOH and incubate at room temperature for  00:10:00 .

4.3 Pellet the beads on a magnet and discard the supernatant containing the non-biotinylated strand.

4.4 Wash the beads once with  50 μL 0.1M NaOH, once with  50 μL of 1X B&W buffer and once with  50 μL TE buffer.4.5 Resuspend beads in  50 μL TE buffer.**Note**

The above steps are designed to produce two libraries of single-stranded, bead-bound SARS-CoV-2 probes for each sample barcode, with the general structure 5'-[biotin][spacer][partial Illumina p5 adapter][sample barcode][generic stub][SARS-CoV-2 target]-3'. Once the method has been suitably tested and is ready for deployment at scale, we envisage that these steps could be simplified substantially by synthesising the biotinylated probe sets in full for each sample barcode and immobilising these directly onto the streptavidin-coated beads.



Sample collection and preparation

- 5 Collect nasopharyngeal swabs in a  1.5 mL volume of transport medium.

Note

We are currently testing a range of media with the aim of identifying a formulation that optimally preserves viral RNA and supports hybridisation with our bead-bound probes in the following step. These include UTM, Lysis/binding buffer and a combination of the two (see Materials). Alternatively, kit-extracted RNA or synthetic control RNA can be used.

- 6 Quantify  10 µL of swab diluent using the Qubit RNA HS assay or similar.

RNA capture

- 7 For each sample to be tested, transfer 2 x  25 µL of unique barcoded beads (one for each probe set) into separate 0.2 mL PCR tubes.

Note

Alternatively, use multi-well PCR plates. For 384-well plates, sample volumes will need to be scaled down in subsequent steps.

- 8 Pellet the beads on a magnet and discard the supernatant.
- 9 Add  50 µL of swab diluent directly to each bead pellet (i.e. 2 x  50 µL per swab sample) and resuspend completely by pipette mixing. Alternatively, use purified or synthetic control RNA in a buffer suitable for probe hybridisation e.g. 1X SSC.
- 10 Place reactions in a preheated thermal cycler with a heated lid and incubate at  75 °C for  00:02:00, then ramp to  25 °C at a rate of -0.1°C/s and proceed immediately to the next step.



- 11 Pellet the beads on a magnet and discard the supernatant.
- 12 Resuspend each reaction in  20 μL of Washing buffer A.
- 13 Combine all “pool 1” samples together into a fresh tube and repeat this process for “pool 2”.

Note

For large sample batches it should only be necessary to pool a small aliquot of each reaction. The combined samples will be used as input for cDNA synthesis. We are currently testing whether there is a limit to the amount of beads that can be tolerated in the reverse transcription reaction, as well as the subsequent library amplification.

cDNA synthesis

- 14 Wash each pooled bead sample once with an equivalent volume of Washing buffer A and twice with the same volume of Washing buffer B, using a magnet to separate the beads from the solution between each washing step. After the final wash, leave the beads suspended in Washing buffer B at room temperature and proceed immediately to the next step.
- 15 Prepare a reverse transcriptase master mix by combining the following components from the SuperScript™ IV First-Strand Synthesis System:

Component	Volume
5X SSIV buffer	20 μL
10 mM dNTP mix (10 mM each)	5 μL
100 mM DTT	5 μL
Ribonuclease inhibitor (40 U/ μL)	5 μL
DEPC-treated water	60 μL



- 16 Pellet beads from step 14 on a magnet and discard the supernatant. Pulse spin tubes for 5 s, return to the magnet and discard any remaining liquid to ensure all the Washing buffer B has been removed (residual buffer may inhibit the reverse transcriptase reaction).
- 17 Perform one final wash in  25 μL of reverse transcriptase master mix. Again pellet beads on a magnet and discard the supernatant.
- 18 Resuspend each pooled bead sample in  19 μL of reverse transcriptase master mix and add  1 μL of SuperScript™ IV Reverse Transcriptase (200 U/ μL). Mix well by pipetting and briefly centrifuge.
- 19 Incubate the reactions at  55 °C for  00:10:00 .
- 20 Inactivate by incubating at  80 °C for  00:10:00 .
- 21 Pellet beads by placing each tube on a magnet. Discard the supernatant and wash once in TE buffer.

Library amplification

22

Note

Here we generate amplicon libraries from each pooled cDNA sample by multiplex PCR. We also incorporate Illumina sequencing adapters and additional barcodes at this stage, producing sequence-ready samples that could potentially be combined with additional library pools for further scalability. Sequencing adapters are added over two sequential rounds of PCR, although we envisage that a single step protocol may also be feasible.

Pellet the “pool 1” and “pool 2” cDNA-coated bead samples from the previous step on a magnet, discard the supernatant and resuspend each in  28 μL of nuclease-free water.

- 23 Combine the *first-round* PCR components in separate 0.2 mL PCR tubes:

Component	Pool 1	Pool 2
-----------	--------	--------

cDNA-coated beads (pool 1 or 2)

26.7 µL

26.7 µL

5X Q5 reaction buffer

10 µL

10 µL

10 mM dNTPs

1 µL

1 µL

Illumina P5 index primer (10 µM)

2.5 µL

2.5 µL

SARS-CoV-2 primer pool 1 or 2 (10 µM)

9.3 µL

9.3 µL

Q5 Hot Start High-Fidelity Polymerase

0.5 µL

0.5 µL

Note

The SARS-CoV-2 primer pools contain 124 (pool 1) and 123 (pool 2) individual primer sequences. The concentration of each primer in the final reaction is 0.015 µM, as recommended in previous studies (Quick *et al.* [2017] *Nat Protoc.* 12: 1261).

24 Mix well by pipetting and briefly centrifuge.

25 Incubate the *first-round* PCR reactions as follows:

Step	Temperature	Time	Cycles
Heat activation	98 °C	00:00:30	1
Denaturation	98 °C	00:00:15	
25-35			
Annealing/extension	65 °C	00:05:00	
25-35			
Hold	4 °C	Indefinite	1

26 For each reaction, pellet the cDNA-coated beads on a magnet and transfer the supernatant containing the amplified product to a new tube.



- 27 Clean up each reaction by performing a left side size selection with a 1.2X volume of SPRIselect beads as per the manufacturer's instructions (Beckman Coulter). Elute the sample in  40 μL of nuclease-free water.
- 28 Quantify  1 μL of each purified sample using the Qubit dsDNA HS assay or similar.
- 29 Combine the *second-round* PCR components in separate 0.2 mL PCR tubes:

Component	Pool 1	Pool 2
Purified template (previous step)	 10 ng	 10 ng
5X PrimeSTAR GXL buffer	 10 μL	 10 μL
dNTP mixture (2.5 mM each)	 4 μL	 4 μL
Illumina P5 flow cell primer (10 μM)	 1.25 μL	 1.25 μL
Illumina P7 index primer (10 μM)	 1.25 μL	 1.25 μL
PrimeSTAR GXL DNA polymerase	 1 μL	 1 μL
Nuclease-free water	to  50 μL	to
 50 μL		

- 30 Mix well by pipetting and briefly centrifuge.

- 31 Incubate the *second-round* PCR reactions as follows:

Step	Temperature	Time	Cycles
Denaturation	 98 $^{\circ}\text{C}$	 00:00:10	6
Annealing	 60 $^{\circ}\text{C}$	 00:00:15	6
Extension	 68 $^{\circ}\text{C}$	 00:01:00	6



Hold

4 °C

Indefinite

1

- 32 Clean up each reaction by performing a left side size selection with a 0.8X volume of SPRIselect beads as per the manufacturer's instructions (Beckman Coulter). Elute the sample in 40 µL of nuclease-free water.

Library quantitation and sequencing

- 33 Determine the concentration of the final "pool 1" and "pool 2" multiplexed libraries in ng/µL using the Qubit dsDNA HS assay or similar.
- 34 Assess the fragment size distribution of each library and determine the average fragment length using an Agilent 2100 Bioanalyzer or similar.
- 35 Calculate sample molarity using the following formula:

$$\text{Concentration in nM} = [(\text{concentration in ng/}\mu\text{L}) / (660 \text{ g/mol} \times \text{average size in bp})] \times 10^6$$

- 36 Combine an equimolar amount of each library and mix.

Note

If sequencing additional samples, quantitate them as above and combine at the desired relative molarity.

- 37 Sequence the final pooled library using an Illumina instrument, generating 2 × 150 nt paired end reads

Alternative bead preparation method

38

Note

This alternative method enables the production of multiple uniquely barcoded bead preparations from a single biotinylated template oligo. This reduces oligo synthesis costs, but does introduce an additional level of protocol complexity.

Wash  200 μL of streptavidin-coated Dynabeads™ MyOne™ C1 magnetic beads as per the manufacturer's instructions:

38.1 Resuspend the Dynabeads™ magnetic beads in the vial (vortex for >  00:00:30).

38.2 Transfer  200 μL of Dynabeads™ magnetic beads to a tube.

38.3 Add an equal volume of 1X B&W buffer and resuspend.

38.4 Place the tube on a magnet for  00:01:00 and discard the supernatant.

38.5 Remove the tube from the magnet and resuspend the washed magnetic beads in  200 μL of 1X B&W buffer.

38.6 Repeat steps 38.4–38.5 twice, for a total of 3 washes.

38.7 Pellet beads once more, and resuspend in  400 μL of 2X B&W buffer.

39 Immobilise a generic oligonucleotide of the structure 5'-[biotin][spacer][partial Illumina p5 adapter]-3' onto the beads as follows:

39.1 Add  400 μL of 5 μM oligonucleotide ([Biotin]-bead tether) to  400 μL of washed beads.

39.2 Incubate for  00:15:00 at room temperature using gentle rotation.

39.3 Separate the biotinylated DNA coated beads with a magnet for 2–3 min.

39.4 Wash the coated beads 2–3 times with 1X B&W buffer.

39.5 Resuspend in  1.200 mL of nuclease-free water.



40 For each unique barcode, generate two libraries of bead-bound SARS-CoV-2 probes via a two-stage annealing and extension process using a barcoded linker oligonucleotide of the structure 3'-[partial Illumina p5 adapter complement][sample barcode][generic stub complement]-5' and two separate oligonucleotide pools with the general structure 3'-[generic stub complement][SARS-CoV-2 target]-5':

40.1 Mix the following components in a 0.2 mL PCR tube:

Component	Pool 1	Pool 2
Bead-bound template (previous step)	 60 µL	 60 µL
5X PrimeSTAR GXL buffer	 20 µL	 20 µL
dNTP mixture (2.5 mM each)	 8 µL	 8 µL
Barcoded linker (10 µM)	 5 µL	 5 µL
SARS-CoV-2 probe pool 1 or 2 (10 µM)	 5 µL	 5 µL
PrimeSTAR GXL DNA polymerase	 2 µL	 2 µL

40.2 Incubate the reaction as follows:

Step	Temperature	Time	Cycles
Denaturation	 98 °C	 00:00:10	4
Annealing	 55 °C	 00:00:15	4
Extension	 68 °C	 00:02:00	4
Hold	 4 °C	Indefinite	4

41 For each bead sample, dissociate and remove the non-biotinylated strand as follows:

41.1 Pellet the beads on a magnet and wash once in  50 µL 1X SSC.



- 41.2 Resuspend the beads in  20 μL of freshly prepared 0.15 M NaOH and incubate at room temperature for  00:10:00 .
- 41.3 Pellet the beads on a magnet and discard the supernatant containing the non-biotinylated strand.
- 41.4 Wash the beads once with  50 μL 0.1M NaOH, once with  50 μL of 1X B&W buffer and once with  50 μL TE buffer.
- 41.5 Resuspend beads in  50 μL TE buffer.