

Jan 15, 2024

a-Synuclein protein expression and purification

DOI

dx.doi.org/10.17504/protocols.io.3byl4brwrvo5/v1



andrew.west west¹, Arpine Sokratian¹

¹Duke University



andrew.west west

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4brwrvo5/v1>

Protocol Citation: andrew.west west, Arpine Sokratian 2024. a-Synuclein protein expression and purification. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4brwrvo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: May 21, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 62989

Keywords: ASAPCRN, synuclein purification, synuclein protein expression, synuclein, purification protocol, purification, protein, rt-quic

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020527

Abstract

Protocol for recombinant a-synuclein purification, useful as monomer template for seeded-amplification assays (like RT_QUIC or PMCA) . Recommendations are to store the protein always on ice while not running and do not stop purification when it is started.

Materials

1. BL21-CodonPlus (DE3)-RIL Chemical Competent Cells (Agilent #230245-41)
2. Thermo Scientific™ Low Protein Binding Collection Tubes (1.5 mL) PI90411

☒ Nalgene™ Rapid-Flow™ Sterile Single Use Vacuum Filter Units **Thermo Fisher Scientific Catalog #565-0010**

☒ SnakeSkin Dialysis Tubing 3.5K MWCO 35 **Thermo Fisher Scientific Catalog #88244**

☒ Endotoxin detection kit LAL **Genscript Catalog #95045-024**

☒ ToxinEraser™ Endotoxin Removal Kit **Genscript Catalog #89233-330**

☒ ToxinEraser™ Endotoxin Removal Resin **Genscript Catalog #L00402**

☒ HiPrep Q HP anion exchange chromatography column **Cytiva Catalog #29018182**

☒ MilliporeSigma™ Amicon™ Ultra-15 Centrifugal Filter Units **Catalog #MilliporeSigma™ UFC901024**



Protocol materials

☒ Nalgene™ Rapid-Flow™ Sterile Single Use Vacuum Filter Units **Thermo Fisher Scientific Catalog #565-0010**

☒ SnakeSkin Dialysis Tubing 3.5K MWCO 35 **Thermo Fisher Scientific Catalog #88244**

☒ Endotoxin detection kit LAL **Genscript Catalog #95045-024**

☒ ToxinEraser™ Endotoxin Removal Kit **Genscript Catalog #89233-330**

☒ ToxinEraser™ Endotoxin Removal Resin **Genscript Catalog #L00402**

☒ HiPrep Q HP anion exchange chromatography column **Cytiva Catalog #29018182**

☒ MilliporeSigma™ Amicon™ Ultra-15 Centrifugal Filter Units **Catalog #MilliporeSigma™ UFC901024**

☒ Nalgene™ Rapid-Flow™ Sterile Single Use Vacuum Filter Units **Thermo Fisher Scientific Catalog #565-0010**

☒ SnakeSkin Dialysis Tubing 3.5K MWCO 35 **Thermo Fisher Scientific Catalog #88244**

☒ HiPrep Q HP anion exchange chromatography column **Cytiva Catalog #29018182**

☒ MilliporeSigma™ Amicon™ Ultra-15 Centrifugal Filter Units **Catalog #MilliporeSigma™ UFC901024**

☒ ToxinEraser™ Endotoxin Removal Kit **Genscript Catalog #89233-330**

☒ ToxinEraser™ Endotoxin Removal Resin **Genscript Catalog #L00402**

☒ Endotoxin detection kit LAL **Genscript Catalog #95045-024**

Troubleshooting

Transformation

1d

- 1 Thaw down an aliquot of plasmid construct (pRK172) encoding WT-human-a-synuclein

[M] 0.3 mg/mL

🌡️ On ice

15m

- 2 Thaw down 🌡️ On ice an aliquot of BL21 (DE3) RIL competent E Coli cells

15m

- 3 Add 🧴 1 µL of plasmid construct to the thawed competent cells and gently mix by flicking the bottom of the tube with a finger a few times



Safety information

do not resuspend

- 3.1 Incubate the reaction mix 🌡️ On ice

15m



- 4 Perform heat-shock transformation 🌡️ 42 °C in water bath incubator with manually shaking at 🌀 100 rpm, 00:00:45

1m



Equipment

Precision™ General Purpose Water Bath

NAME

Water Bath

TYPE

Thermo Scientific

BRAND

TSGP10

SKU



5 Immediately transfer the tube on ice and incubate for 1 min.

1m

6 Add  1000 μL of SOC media to a chilled reaction

10s

7 Incubate the bacteria  200 rpm, 37°C, 00:30:00

30m

Equipment

ThermoMixer® C

NAME

ThermoMixer


TYPE



Eppendorf

BRAND


5382000023

SKU

7.1 Prepare sterile 10cm LB agar plate containing  0.1 mg/mL of ampicillin

8 Centrifuge at  500 x g, 10°C, 00:03:00 and discard the supernatant leaving  50 μL of media

3m

8.1 Spread 50 μL of cell suspension onto a selection plate and incubate overnight at  37 °C in bacterial incubator

10m



Equipment

Isotemp™ Microbiological Incubator, 178 L, Stainless Steel^{NAME}



Microbiological Incubator^{TYPE}

Fisherbrand^{BRAND}

15-103-0513^{SKU}

Protein expression

12h

- 9 Pick one colony and transfer into  10 mL LB media with  0.1 mg/mL of ampicillin start in the morning (9:00 am)

- 9.1 Incubate the bacteria  250 rpm, 37°C, 05:00:00 until it reaches OD 0.2-0.3

5h

Equipment

Natural convection incubator^{NAME}

Bacterial shaker^{TYPE}

Innova^{BRAND}

M1335-0000^{SKU}



10 Transfer a starter culture to 2X2L flasks filled with 0.5L LB media with **[M] 0.1 mg/mL** of ampicillin
5 mL to each flask

11 Incubate the culture at the same conditions until it reaches OD 0.8 (use nanodrop or cuvette) (reaches optimal density at 6-7 pm)

5h

12 Induce protein expression by adding **[M] 0.05 millimolar (mM)** IPTG, incubate at **🔥 18 °C** for **🕒 12:00:00** overnight

12h



Note

To cool down the grown culture, transfer the flasks into ice-bath and incubate until it reaches desired temperature

Cell lysis

11m 45s

13 Collect the pellets by centrifugation (JA14 rotor) at **🌀 5000 x g** x g at **🔥 4 °C** for **🕒 00:10:00** . Used 250 ml Beckman tubes
Usually get 10-12 g from 2L

10m

14 Add to pellets 80 ml of lysis buffer (total): **[M] 10 millimolar (mM)** ThisHCl **📏 7.6** ,
[M] 750 millimolar (mM) NaCl, **[M] 1 millimolar (mM)** EDTA, **[M] 1 millimolar (mM)** PMSF (add just before using, have aliq frozen **[M] 0.1 Molarity (M)**), protease inhibitors (use MAXI version, need only one tablet);

15 Carefully resuspend the pellets to homogenize the solution

15.1 Heat up **🔥 1 L** of water in a high temperature resistant glass beaker (turn heat to the max on the magnetic stirrer)

16 While waiting on water to get to the boiling point sonicate the lysates (use thick prob-tip) for **🕒 00:01:00** , 30%, **🕒 00:00:15** ON **🕒 00:00:30** OFF of amplitude then go to next falcon, had 3 falcons (repeat 3 times, avoid overheating)

10m



- 17 After sonication samples need to get boiled thereby put the falcon tubes into glass beaker and boil for 00:25:00 . Use tweezers to pull out the tubes 25m
- 18 Transfer boiled homogenates into new 50 mL falcon tubes; chill down suspensions at room temperature for 20 min 20m
- 19 Prepare 4 L of buffer 10 millimolar (mM) TrisHCl 7.6 , 50 millimolar (mM) NaCl, 1 millimolar (mM) EDTA, 1 millimolar (mM) PMSF for dialysis
- 20 Centrifuge the homogenates at 20000 x g for 01:00:00 at 4 °C 1h
- 21 Filter the supernatant using 0.45 um filter unit
 Nalgene™ Rapid-Flow™ Sterile Single Use Vacuum Filter Units **Thermo Fisher Scientific Catalog #565-0010**
- 22 Transfer filtered supernatant into dialysis bag which is: SnakeSkin Dialysis tubing, 3.5K MWCO, 35 mm dry I.D., 35 feet.
Measure the dialysis tube taking into consideration that 5 cm length of tube holds 48 mL of the sample (plus 2.5cm at each end for closure). Clip the tube using green clips, make sure it does not leak.
Place the dialysis bag into 4 L plastic beaker filled with dialysis buffer, incubate overnight on magnetic plate on the slow mode (Chromatography fridge)
 SnakeSkin Dialysis Tubing 3.5K MWCO 35 **Thermo Fisher Scientific Catalog #88244**



Equipment

ÄKTA pure 25 L1

NAME

ÄKTA pure chromatography system

TYPE


Cytiva

BRAND

29018225

SKU

Protein purification (anion-exchange chromatography)

23 After a night of dialysis ( 4 °C slow mixing) collect the suspension into 100 mL glass bottle (filter the sample before running on the column, 0.22 um filter).


24 Column - HiPrep Q HP 16/10 column 1×20 ml (stored in 70% ethanol);

 HiPrep Q HP anion exchange chromatography column **Cytiva Catalog #29018182**

24.1 Wash the column 2V of miliQ degassed water

24.2 Wash the column with 2V of **STARTING BUFFER** [M] 10 millimolar (mM) TrisHCl


 7.6 , [M] 50 millimolar (mM) NaCl

24.3 Activate with 1V of [M] 10 millimolar (mM) TrisHCl  7.6 , [M] 1 Molarity (M) NaCl

24.4 Equilibrate with 3V of starting buffer


25 Load  80 mL of suspension and then washed with 100 ml [M] 50 millimolar (mM)





NaCl [M] 10 millimolar (mM) TrisHCL,  300 mL of gradient elution (0-100%), 2 ml/min flow rate. Collected samples using fraction collector 2, every fraction 4 ml (use 10 ml glass tubes)

26 Place supernatant into channel A1 (was previously use for starting buffer, do not generate bubbles)



27 Place starting buffer in channel A2 (clean the tubing using the program mode)

28 Place elution buffer in channel B1 ([M] 10 millimolar (mM) TrisHCl  7.6 ,
[M] 1 Molarity (M) NaCl)
Collected samples using fraction collector 2, every fraction 4 ml (use 10 ml glass tubes);

29 Analyze the fractions eluted at 250-350 mM salt (20 RFU conductivity) though SDS-PAGE (stain with Coomassie).
Combine  10 μ L of each fraction with  10 μ L of 2X laemmli buffer and analyze fractions by SDS-PAGE with 4–20% gradient gels, followed by coomassie staining/destaining


30 Measure A280/260 for the fractions containing single a-syn band, avoid collecting samples with A280/260 > 0.85

31 Combine the evaluated factions and measure total protein concentration using nanodrop.

32 Dialyze with  4 L of [M] 10 millimolar (mM) TrisHCl  7.6 ,
[M] 50 millimolar (mM) NaCl (follow instruction for dialysis)



Further purification

33 Repeat section 'Protein purification (anion-exchange chromatography)' for the further fractionation of the purified preparation  [go to step #23](#)



Protein concentration

10m

34 Concentrate dialyzed protein sample to approximately [M] 30 mg/mL aliquot



Prepare the ultra-concentration system

35 Use 50 mL ultra centrifugation units with 3K cutoff



MilliporeSigma™ Amicon™ Ultra-15 Centrifugal Filter
Units **Catalog #**MilliporeSigma™ UFC901024

36 Wash off the unit with miliQ water through centrifugation at 5000 x g at 4 °C
for 00:05:00 , JA10 rotor

5m

37 Load first 15 mL of the sample into ultracentrifugation unit (max load of the unit is approx. 15 mL)

38 Centrifuge at 5000 x g at 4 °C for 00:05:00 , JA10 rotor

5m

39 Resuspend concentrated sample, add more of protein sample and concentrate until the total volume is ~ 5 mL

40 Store at -80 °C . Yield should be approximately 80 mg per 2 L culture



Endotoxin removal

41 Follow instructions for



ToxinEraser™ Endotoxin Removal Kit **Genscript Catalog #**89233-330 with
modifications

For a more successful endotoxin removal, add 1 mL of



ToxinEraser™ Endotoxin Removal Resin **Genscript Catalog #**L00402 before the
regeneration process

42 Collect the eluate into 5 mL endotoxin-free tube and save 2 aliquots (10 µL and
 50 µL) for protein concentration and endotoxin measurements

Endotoxin quantification



43 Follow instructions for Endotoxin detection kit LAL Genscript

 Endotoxin detection kit LAL **Genscript Catalog #95045-024**