


Feb 12, 2025

🌐 Absolute Proteome Quantification of Extracellular Proteins

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygxwnyov8j/v1>



Borja Ferrero Bordera^{1,2}, Sandra Maaß³

¹LMU Munich; ²University Greifswald; ³University of Greifswald



Borja Ferrero Bordera

LMU Munich

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygxwnyov8j/v1>

External link: <https://journals.asm.org/doi/full/10.1128/spectrum.02616-23#sec-4>

Protocol Citation: Borja Ferrero Bordera, Sandra Maaß 2025. Absolute Proteome Quantification of Extracellular Proteins. protocols.io <https://dx.doi.org/10.17504/protocols.io.kxygxwnyov8j/v1>

**Manuscript citation:**

Ferrero-Bordera B, Bartel J, van Dijk JM, Becher D, Maaß S. 2024. From the outer space to the inner cell: deconvoluting the complexity of *Bacillus subtilis* disulfide stress responses by redox state and absolute abundance quantification of extracellular, membrane, and cytosolic proteins. *Microbiol Spectr* 12:e02616-23. <https://doi.org/10.1128/spectrum.02616-23>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 20, 2024

Last Modified: February 12, 2025

Protocol Integer ID: 116532

Keywords: Proteomics, Mass-Spectrometry, Absolute Quantification, Bioinformatics, Python, Exoproteome, Extracellular proteins, absolute proteome quantification of extracellular protein, extracellular proteins in bacteria, proteins in the culture medium, extracellular protein, absolute proteome quantification, protein, standard protein, secreted protein, bacteria

Funders Acknowledgements:

People Programme (Marie Skłodowska-Curie Actions) of the European Union's Horizon 2020 Programme

Grant ID: 813979

Abstract

To quantify secreted proteins in the culture medium, we described a simple and straightforward protocol for the absolute quantification of extracellular proteins in bacteria. We concentrated extracellular proteins, which are highly diluted in the medium, using StrataClean resin along with a set of standard proteins to determine the extent of the concentration step.



Materials

Equipment

- Low-binding tubes
- Centrifuge
- Overhead shaker
- SpeedVac centrifuge

Common Reagents & Buffers

- Tris-EDTA (TE) buffer (10 mM Tris-HCl, 1 mM EDTA, pH 7.5)
- StrataClean resin (Agilent)
- 37 % HCl
- Loading buffer (4% (w/v) SDS, 20% (v/v) glycerol, 40 µg/mL bromophenol blue, 125 mM Tris-HCl pH 6.8 with 10% (v/v) β-mercaptoethanol)
- Porcine Trypsin

Protein Standards

- Anchor proteins (e.g. Universal Protein Standards -UPS2- from Sigma)
- Enrichment standards (described in this protocol)

Extracellular protein extraction

1d 1h

- 1 Pellet by centrifuging at 4 °C an aliquot of the culture to be analysed.

10m

 8500 x g, 4°C, 00:10:00

- 2 Recover the supernatant into a new tube.

10m

Note

Low-binding tubes are recommended when preparing proteomics samples.

- 3 Add the Enrichment Standards to the raw supernatant. **Important!** This step is crucial to quantify the total protein in the raw supernatant allowing to estimate the protein binding efficiency.

10m



Standards were added based on the supernatant's protein concentration in a ratio of 1:35 (ng of standard mixture: ng of estimated sample protein). Below is the table of proteins used as standards in our paper.

	Protein	Organism	Accession	Supplier	Concentration (ng/μL)	Protein amount (fmol)
	alpha-Lactalbumin	Bos taurus (Bovine)	P00711	Sigma Merck	174.07	5000
	Glycerol-3-Phosphate dehydrogenase	Oryctolagus cuniculus (Rabbit)	P46406	Sigma Merck	460.05	5000
	Alcohol dehydrogenase	Saccharomyces cerevisiae (Yeast)	P00330	Sigma Merck	45.76	500



	Protein	Organism	Accession	Supplier	Concentration (ng/ μ L)	Protein amount (fmol)
	Soybean trypsin inhibitor	Glycine max (Soybean)	P01071	Gibco	24.87	500
	Lyszyme	Gallus gallus (Chicken)	P00698	Pierce	17.78	500
	Bovine serum albumin	Bos taurus (Bovine)	P02769	Pierce	8.60	50
	Carbonic anhydrase	Bos taurus (Bovine)	Q1LZA1	Sigma Merck	3.61	50

- 4 Add 20 μ L HCl-primed StrataClean resin (Agilent; REF 400714) into the raw supernatant and incubate at 4 °C overnight with overhead shaking.  800 rpm, 4°C, 16:00:00

1d

Note

Other protocols (TCA/Acetone precipitation) can be used as they showed comparable efficiency to the one described here. Nonetheless, StrataClean beads reduce the usage of organic solvent

- 5 Centrifuge at maximum speed at 4 °C and discard the supernatant to recover the resin with the bound proteins. Wash the resin with TE buffer and repeat the centrifugation step.
- 6 Vacuum-dry the resin containing the bound protein.

30m

30m



Protein Electroelution

1h 40m

- 7 Resuspend the resin-bound protein in 20 μ L loading buffer (consider the gel pocket size to adjust the final volumes, as it will be combined with the anchor proteins in the well).

10m

Note

It is recommended to resuspend the dried resin in 2X the original amount to ease the pipetting into the gel as the resin is very viscous. For example, 20 μ L initial resin can be divided into 2 eppis with 10 μ L each, dried and resuspended in 20 μ L loading buffer.

- 8 Load the resin with the bound protein into the gel pocket. After sample loading, add the anchor proteins used for quantification (e.g. UPS2) into the gel pocket on top.

10m

- 9 Run the SDS-PAGE gel to the desired time. Lower times are recommended for samples with low complexity and/or highly diluted.

50m

- 10 Proceed with Coomassie staining to visualize the eluted protein in the gel lanes.

30m

In-Gel protein digestion & desalting

1d 1h 5m

- 11 Fractionate the gel lanes in squares of 2×2 mm

20m

- 12 Destain the gel with 200 mM ammonium bicarbonate buffer in 30% (vol/vol) acetonitrile with  00:15:00 incubations at  37 °C . Repeat until completely destained.

45m

- 13 Vacuum-dry the destained gel pieces and resuspend them in 2 μ g/ μ L trypsin. Incubate overnight at 37 °C.

1d



- 14 To elute the digested peptides, resuspend the digested gel pieces in MilliQ-water and sonicate for 15 minutes. Centrifuge at maximum speed and transfer the supernatant to a low-binding eppi.

- 15 Vacuum-dry and resuspend in 0.1 % MS-grade acetic acid.

- 16 Proceed with desalting. We recommend ZipTip (Millipore) desalting because of the simplicity of the protocol. Nonetheless, other protocols like StageTipping can be



performed.

- 17 Vacuum-dry and resuspend with Buffer A for MS measurement of the samples.