


Jun 01, 2023

🌐 Abundance of fungal hyphae in seawater by epifluorescence microscopy using Calcofluor White stain method

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgq933klk5/v1>



M.H. Gutiérrez¹

¹COPAS Sur-Austral, Universidad de Concepción, Chile



Amy Gladfelter

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgq933klk5/v1>

Protocol Citation: M.H. Gutiérrez 2023. Abundance of fungal hyphae in seawater by epifluorescence microscopy using Calcofluor White stain method. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgq933klk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 02, 2018

Last Modified: June 01, 2023

Protocol Integer ID: 17442

Keywords: fungal hyphae from seawater, fungal hyphae, fungal, calcofluor white, hyphae, seawater, calcofluor, epifluorescence, stain

Abstract

Protocol to stain fungal hyphae from seawater with Calcofluor White.

Attachments



Abundance of fungi i...

30KB

Guidelines

References

Damare S, Raghukumar C (2008) Fungi and macroaggregation in deep-sea sediments. *Microb Ecol* 56:168–177

Jebarag CS, Raghukumar C (2009) Anaerobic denitrification in fungi from the coastal marine sediments off Goa, India. *Mycol Res* 113:100–109

Newell, S.Y., Statzell-Tallman, A. (1982) Factors for conversion of fungal biovolume values to biomass, carbon and nitrogen: variation with mycelial ages, growth conditions, and strains of fungi from a salt marsh. *Oikos* **39**: 261–268

Rasconi S, Jobard M, Jouve L, Sime-Ngando T (2009) Use of calcofluor white for detection, identification, and quantification of phytoplanktonic fungal parasites. *Appl Environ Microbiol* 75:2545–2553

Safety warnings

! See SDS (Safety Data Sheet) for hazards and safety guidelines.

Collection

- 1 Collect samples in sterile 50 mL polypropylene tubes and fix with formaldehyde or glutaraldehyde (2% final concentration). Store samples at  4 °C in the dark.

Filtration

- 2 Filter  5 mL to  30 mL (depending on the environment) of seawater by 0.22-µm mesh black 25 mm diameter polycarbonate filters (Millipore Corp.).

Stain

- 3 Stain filters with the retained material from filtration by applying  600 µL of aqueous 0.1% Calcofluor White directly to sample and filter, making sure to cover the entire area of the filter.

Allow  00:05:00 to  00:10:00 then use vacuum to remove the excess stain from the filter. **Avoid any light source.**

Abundance

- 4 Place the filter (sample side up) onto a slide and add 1 drop of non-fluorescent immersion oil to the top of the filter, then cover with a cover slip.
- 5 Count slides immediately on epifluorescent microscope or store at  -20 °C .
- 6 Use an epifluorescence microscope equipped with UV filter (used for DAPI, eg. filter set 49 Carl Zeiss Ltd., 365 nm excitation and 445-450 nm emission band pass) to examine at 1000X the entire effective area of the filters.
- 7 Count all hyphae identified and record their length and width. Use cylinder volume as a morphological approximation to estimate the biovolume of fungal filaments.



Note

Biomass can be estimated based from biomass: biovolume ratios described from fungi (Newell and Statzell-Tallman, 1982).