



Dec 22, 2025

Accurate detection of somatic mutations in single cells by scNanoSeq

DOI

dx.doi.org/10.17504/protocols.io.x54v9be5ml3e/v1

Muchun Niu¹, Yang Zhang¹, Jiayi Luo¹, Chenghang Zong¹

¹Baylor College of Medicine



Muchun Niu

Baylor College of Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9be5ml3e/v1>

External link: <https://www.biorxiv.org/content/10.1101/2023.10.09.561356v1>

Protocol Citation: Muchun Niu, Yang Zhang, Jiayi Luo, Chenghang Zong 2025. Accurate detection of somatic mutations in single cells by scNanoSeq. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9be5ml3e/v1>

Manuscript citation:

Niu, M., Zhang, Y., Luo, J., Sinson, J. C., Thompson, A. M., & Zong, C. (2023). Characterization of cancer evolution landscape based on accurate detection of somatic mutations in single tumor cells. *bioRxiv*, 2023-10.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 19, 2025

Last Modified: December 22, 2025

Protocol Integer ID: 235394

Keywords: cell whole genome amplification method, cell genome amplification method, whole genome amplification method, accurate detection of somatic mutation, somatic mutations in individual human cell, detecting somatic mutation, errors per human genome, somatic mutation, single cell, human genome, cell expansion experiment, sequencing scheme, mutation, individual human cell, nanoseq chemistry, based nanoseq chemistry

Funders Acknowledgements:

NIH Director's New Innovator Award

Grant ID: 1DP2EB020399

NIH SMaHT Program

Grant ID: 1UG3NS132132

McNair Scholarship

NIH SMaHT Program

Grant ID: 1UM1DA058229-01

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Accurately detecting somatic mutations in individual human cells poses a significant challenge for single-cell whole genome amplification methods. In this study, we introduce a novel single-cell genome amplification method based on the restriction enzyme-based NanoSeq chemistry, achieving an error rate of less than 10^{-10} in mutation calling (equivalent to 0.2 errors per human genome). We term this methodology as scNanoSeq. This precision is derived from theoretical error rate analyses within the duplexing sequencing scheme, which aligns with upper bound estimations derived from single-cell expansion experiments.



Protocol materials

- ✕ 1M KCl **bioworld Catalog #40120947-1**
- ✕ NP-40 Surfact-Amps™ Detergent Solution **Thermo Fisher Scientific Catalog #85124**
- ✕ 0.5M EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889-100ML**
- ✕ QIAGEN Protease (7.5 AU) **Qiagen Catalog #19155**
- ✕ Molecular grade water **Catalog #BP561-1 1L**
- ✕ 1M Tris-HCl, pH 8.0 **Thermo Fisher Scientific Catalog #15568025**
- ✕ 10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #93443-100ML**
- ✕ rCutSmart Buffer **New England Biolabs Catalog #B6004S**
- ✕ HpyCH4V - 500 units **New England Biolabs Catalog #R0620L**
- ✕ Nuclease free water
- ✕ CutSmart® Buffer **New England Biolabs Catalog #B7204S**
- ✕ NEBuffer 4 - 5.0 ml **New England Biolabs Catalog #B7004S**
- ✕ Klenow Fragment (3'-5' exo-) - 1,000 units **New England Biolabs Catalog #M0212L**
- ✕ Adenosine 5-Triphosphate (ATP) **New England Biolabs Catalog # P0756L**
- ✕ 15 µM xGen CS adapter **Integrated DNA Technologies, Inc. (IDT) Catalog #1080799**
- ✕ T4 DNA Ligase **New England Biolabs Catalog #M0202L**
- ✕ Ampure XP beads **Beckman Catalog #A63881**
- ✕ iTaq Universal SYBR Green Supermix **Bio-Rad Laboratories Catalog #172-5112**
- ✕ NEBNext Ultrall Q5 Master Mix **New England Biolabs Catalog #M0544X**

Troubleshooting



1. Single cell lysis

3h 55m

1 Prepare  1000 μL **Protease Lysis Mix**.

5m


 30 μL of  1M Tris-HCl, pH 8.0 Thermo Fisher Scientific Catalog #15568025


 4 μL of


 0.5M EDTA Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889-100ML



 20 μL of  1M KCl bioworld Catalog #40120947-1

 10 μL of

 10% Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #93443-100ML

 36 μL of

 NP-40 Surfact-Amps™ Detergent Solution Thermo Fisher Scientific Catalog #85124

 400 μL of [M] 20 mg/mL  QIAGEN Protease (7.5 AU) Qiagen Catalog #19155

 500 μL of  Molecular grade water Catalog #BP561-1 1L

2 Aliquot  2.5 μL of **Protease Lysis Mix** in each low-binding PCR tubes.

20m

3 FACS sort individual cells/nuclei into PCR tubes containing **Protease Lysis Mix**.

1h

4 Brief centrifugation, and run the following cell lysis program on a thermo-cycler.


2h 30m

 50 °C  02:00:00

 75 °C  00:20:00

 85 °C  00:05:00

 4 °C Hold

5 Store the lysed cells at  -80 °C .

Genome fragmentation

50m

6 Prepare the following **Fragmentation Mix** (2.5 uL per cell).



5m

- 0.5 μ L of CutSmart® Buffer **New England Biolabs Catalog #B7204S** or
 rCutSmart Buffer **New England Biolabs Catalog #B6004S**
0.15 μ L of HpyCH4V - 500 units **New England Biolabs Catalog #R0620L**
1.85 μ L of Nuclease free water

- 7 Add 2.5 μ L of **Fragmentation Mix** to each tube along the wall, tap to mix, and then centrifuge.

45m

Run the following program on a thermo-cycler.

- 37 °C 00:15:00
 65 °C 00:20:00

End repair

45m

- 8 Prepare the following **End Repair Mix** (10 μ L per cell).

5m

- 1 μ L of NEBuffer 4 - 5.0 ml **New England Biolabs Catalog #B7004S**
7.35 μ L of Nuclease free water
0.15 μ L of
 Klenow Fragment (3'-5' exo-) - 1,000 units **New England Biolabs Catalog #M0212L**

- 9 Add 10 μ L of **End Repair Mix** to each tube along the wall, tap to mix, and then centrifuge.

40m

Run the following program on a thermo-cycler.

- 37 °C 00:30:00

Y-shape adapter ligation

1h 7m

- 10 Prepare the following **Ligation Mix** (22.4 μ L per cell).

5m

- 2.24 μ L of NEBuffer 4 - 5.0 ml **New England Biolabs Catalog #B7004S**
15.53 μ L of Nuclease free water

3.74 μL of

Adenosine 5-Triphosphate (ATP) **New England Biolabs Catalog # P0756L**

0.33 μL of

15 μM xGen CS adapter **Integrated DNA Technologies, Inc. (IDT) Catalog #1080799**

0.56 μL of T4 DNA Ligase **New England Biolabs Catalog #M0202L**

- 11 Add 22.4 μL of **Ligation Mix** to each tube along the wall, tap to mix, and then centrifuge.

32m

Run the following program on a thermo-cycler.

20 $^{\circ}\text{C}$ 00:22:00

- 12 Immediately after ligation, purify the ligation product with 0.95X Ampure XP beads **Beckman Catalog #A63881**, and eluted in 15.5 μL of Nuclease free water.

30m

To avoid sample loss, flick the tubes for all the mixing steps instead of pipetting.

Library QC and amplification

4h

- 13 Use 0.5 μL of Purified ligation product for the following qPCR yield test.

2h

5 μL of

iTaq Universal SYBR Green Supermix **Bio-Rad Laboratories Catalog #172-5112**

0.25 μL of [M] 10 micromolar (μM) **Truseq5 primer**

(ACACTCTTTCCCTACACGAC)

0.25 μL of [M] 10 micromolar (μM) **Truseq7 primer**

(GTGACTGGAGTTCAGACGTGT)

4 μL of Nuclease free water

0.5 μL Purified ligation product

Run qPCR on a Roche LightCycler 96 machine following the program:

94 $^{\circ}\text{C}$ for 00:02:00

30 cycles of

- 94 °C for 00:00:20
- 58 °C for 00:00:20
- 72 °C for 00:01:00

Melting Curve

For a typical diploid cells, qPCR Ct value is expected to be 17.5~18.5.

- 14 For the cells with anticipated yield, use the remaining 15 µL of
Purified ligation product for PCR amplification.

2h

25 µL of

NEBNext Ultrall Q5 Master Mix **New England Biolabs Catalog #M0544X**

5 µL of [M] 10 micromolar (µM) **P5_index5_Truseq5 primer**

(AATGATACGGCGACCAACGAGATCTACAC[8-base-index]ACACTCTTTCCCTACACGACGCTCTTCCGATCT)

5 µL of [M] 10 micromolar (µM) **P7_index7_Truseq7 primer**

(CAAGCAGAAGACGGCATACGAGAT[8-base-index]GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT)

15 µL Purified ligation product

Run the following PCR program on a thermo-cycler:

98 °C for 00:00:30

17 cycles of

▪ 98 °C for 00:00:10

▪ 73 °C for 00:01:15

73 °C for 00:05:00

Purify the PCR product twice with 0.7X

Ampure XP beads **Beckman Catalog #A63881** , and eluted in 20 µL of

Nuclease free water .

Pooled libraries are ready for sequencing on an Illumina NovaSeq 6000 / NovaSeq X / NovaSeq XPlus platform with 2×150 bp paired-end mode. The target depth is 5~6 reads per fragment based on pre-amplification qPCR quantification.



Protocol references

1. Niu, M., Zhang, Y., Luo, J., Sinson, J. C., Thompson, A. M., & Zong, C. (2023). Characterization of cancer evolution landscape based on accurate detection of somatic mutations in single tumor cells. *bioRxiv*, 2023-10.
2. Abascal, F., Harvey, L. M., Mitchell, E., Lawson, A. R., Lensing, S. V., Ellis, P., ... & Martincorena, I. (2021). Somatic mutation landscapes at single-molecule resolution. *Nature*, 593(7859), 405-410.