


Aug 19, 2022

🌐 ACE-inhibitory activity assay: IC50

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g74q5kgwz/v1>



Ezequiel Coscueta¹, & Maria Manuela Pintado²

¹Universidade Catolica Portuguesa;

²Centro de Biotecnologia e Química Fina (CBQF)-Laboratório Associado, Escola Superior de Biotecnologia, Universidade Católica Portuguesa, Rua Diogo Botelho 1327, 4169-005 Porto, Portuga



Ezequiel R. Coscueta

Universidade Catolica Portuguesa

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g74q5kgwz/v1>

Protocol Citation: Ezequiel Coscueta, & Maria Manuela Pintado 2022. ACE-inhibitory activity assay: IC50. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g74q5kgwz/v1>

Manuscript citation:

Coscueta, E. R., Campos, D. A., Osório, H., Nerli, B. B., & Pintado, M. (2019). Enzymatic soy protein hydrolysis: A tool for biofunctional food ingredient production. *Food Chemistry: X*, 1(100006). <https://doi.org/10.1016/j.fochx.2019.100006>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 17, 2022

Last Modified: August 19, 2022

Protocol Integer ID: 68779

Keywords: ACE, angiotensin-converting enzyme, antihypertensive, IC₅₀, bioactivity, ngiotensin-I converting enzyme, angiotensin, inhibitory activity assay, ic₅₀ in inhibition, nitrophenylalanylproline, enzyme, terminal dipeptide, peptide substrate, action of ace, converting enzyme, assay, ace

Abstract

This protocol describes the procedure for the determination of the IC₅₀ in inhibition on angiotensin-I converting enzyme (ACE) activity. ACE is also known as peptidyl dipeptidase A because it removes C-terminal dipeptides from a wide variety of peptide substrates. In the assay described here, the chosen substrate is the intramolecularly quenched fluorescent tripeptide o-aminobenzoylglycyl-p-nitrophenylalanylproline (Abz-Gly-Phe(NO₂)-Pro). Hydrolysis of this substrate by the action of ACE generates the fluorescent product o-aminobenzoylglycine (Abz-Gly).

Guidelines

ACE-inhibitory activity is measured by fluorescence using the method of Sentandreu & Toldrá (2006) with some modifications (Coscueta, Campos, Osório, Nerli, & Pintado, 2019). A total of 40 µL of ultrapure water or ACE working solution are added to each microtiter-plate well, then adjusted to 80 µL by adding ultrapure water to blank, control or samples. For direct or 1/2 diluted samples a sample blank is also made. The enzyme reaction is started by adding 160 µL of substrate solution and the mixture is incubated at 37 °C. Serial dilutions are made of each sample, usually 1/1 to 1/32. The fluorescence generated is measured at 30 min using a Multidetector plate reader (Synergy H1, Vermont, USA). The assay is performed in a black 96-well microplate (Nunc, Denmark). Excitation and emission wavelengths are 350 and 420 nm, respectively.

Note: In some cases it is necessary to extend the range to 1/5 or 1/10 serial dilutions.

Note: Generally the most concentrated dilutions are the ones that can fluoresce as the sample itself, so in the example sample blanks (SPLB) are used only for the two most concentrated dilutions. In the case of highly fluorescent samples, SPLB should be used for all dilutions.

Citation

Sentandreu, M. Á., & Toldrá, F. (2006)
. A rapid, simple and sensitive fluorescence method for the assay of angiotensin-I converting enzyme.
Food Chemistry.

<https://doi.org/10.1016/j.foodchem.2005.06.006>

LINK

Citation

Coscueta, E. R., Campos, D. A., Osório, H., Nerli, B. B., & Pintado, M. (2019)
. Enzymatic soy protein hydrolysis: A tool for biofunctional food ingredient production. Food Chemistry: X.

<https://doi.org/10.1016/j.fochx.2019.100006>

LINK

Materials

- ⊗ Angiotensin Converting Enzyme from rabbit lung **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A6778**
- ⊗ Tris Base **Fisher Scientific Catalog #BP152-1** or ⊗ Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**
- ⊗ Ultrapure Water
- ⊗ ZnCl₂ anhydrous **Thermo Fisher Scientific Catalog #11497737**
- ⊗ Sodium Chloride **Fisher Scientific Catalog #S271**
- ⊗ Abz-Gly-p-nitro-Phe-Pro-OH Trifluoroacetate **bachem Catalog #4003531.0050**

Equipment

Synergy H1M

NAME

Multidetetection plate reader

TYPE

BioTek

BRAND

679SH1M-SN

SKU

Detection modes: UV-Vis absorbance; fluorescence intensity; luminescence ^{SPECIFICATIONS}
 Read methods: endpoint, kinetic, spectral scanning, well area scanning
 Microplate types: 6- to 384-well plates
 Temperature control: to 45 °C
 Shaking: linear, orbital, double orbital

Equipment

96 Well Black Plate	NAME
microplate	TYPE
Nunc	BRAND
237108	SKU
96-Well black polystyrene microplate, flat bottom, 50 to 250 µL, non-sterile, non-treated	SPECIFICATIONS
	

Before start

Prepare the necessary reagents carefully.

[M] 1 U/mL ACE stock solution

- Dissolve the

 Angiotensin Converting Enzyme from rabbit lung **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A6778**

(peptidyl-dipeptidase A, EC 3.4.15.1) in a solution of 50% glycerol in ultrapure water, to obtain a final concentration of 1 U/mL.

- Make aliquots of  200 μ L of the solution and store at  -20 $^{\circ}$ C .

[M] 0.150 Mass Percent Tris buffer 8.3

- Dissolve  1.817 g of Tris base (MW = 121.14) or  2.364 g Tris-HCl (MW = 157.60) in approx.  90 mL ultrapure water .
- Titrate to  8.77 at the lab temperature of  Room temperature with monovalent strong base or acid as needed.
- Make up volume to  100 mL with ultrapure water.

Buffer will be  8.3 at  37 $^{\circ}$ C .

[M] 42 mU/mL ACE work solution

- Prepare a solution of [M] 1 millimolar (mM) ZnCl_2 , dissolving  1.4 mg ZnCl_2 (MW = 136.28) in  10 mL ultrapure water . Store at  -20 $^{\circ}$ C .
- Prepare [M] 0.1 millimolar (mM) ZnCl_2 [M] 0.150 Mass Percent Tris buffer  8.3 (Enzyme buffer). Dilute 1/10 the previous solution ([M] 0.1 millimolar (mM)) and add  25 μ L of this solution to  25 mL of [M] 0.150 Mass Percent Tris buffer  8.3 . Store at  4 $^{\circ}$ C for a maximum of one week, or at  -20 $^{\circ}$ C for a maximum of six months.
- Dilute 1/24 the [M] 1 U/mL ACE stock solution with the Enzyme buffer. Prepare daily.

[M] 0.45 millimolar (mM) Substrate solution

- Prepare [M] 1.125 Mass Percent NaCl [M] 0.150 Mass Percent Tris buffer  8.3 (Substrate buffer). Dissolve  3.2872 g NaCl (MW = 58.44) in  50 mL of [M] 0.150 Mass Percent Tris buffer  8.3 . Store at  4 $^{\circ}$ C for a maximum of one week, or at  -20 $^{\circ}$ C for a maximum of six months.
- Dissolve 3.6 mg of substrate  Abz-Gly-p-nitro-Phe-Pro-OH Trifluoroacetate **bachem Catalog #4003531.0050** in



16 mL Substrate buffer (for 96 wells). Prepare at the moment.

Analysis

30m

1

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	BLK	BLK	SPLB1:1 1	SPLB1:1 1	SPLB2:1 1	SPLB2:1 1	SPLB3:1 1	SPLB3:1 1	SPLB4:1 1	SPLB4:1 1	
B	CTL	CTL	CTL	SPLB1:2 2	SPLB1:2 2	SPLB2:2 2	SPLB2:2 2	SPLB3:2 2	SPLB3:2 2	SPLB4:2 2	SPLB4:2 2	
C	SPL1:1 1	SPL1:1 1	SPL1:1 1	SPL2:1 1	SPL2:1 1	SPL2:1 1	SPL3:1 1	SPL3:1 1	SPL3:1 1	SPL4:1 1	SPL4:1 1	SPL4:1 1
D	SPL1:2 2	SPL1:2 2	SPL1:2 2	SPL2:2 2	SPL2:2 2	SPL2:2 2	SPL3:2 2	SPL3:2 2	SPL3:2 2	SPL4:2 2	SPL4:2 2	SPL4:2 2
E	SPL1:3 4	SPL1:3 4	SPL1:3 4	SPL2:3 4	SPL2:3 4	SPL2:3 4	SPL3:3 4	SPL3:3 4	SPL3:3 4	SPL4:3 4	SPL4:3 4	SPL4:3 4
F	SPL1:4 8	SPL1:4 8	SPL1:4 8	SPL2:4 8	SPL2:4 8	SPL2:4 8	SPL3:4 8	SPL3:4 8	SPL3:4 8	SPL4:4 8	SPL4:4 8	SPL4:4 8
G	SPL1:5 16	SPL1:5 16	SPL1:5 16	SPL2:5 16	SPL2:5 16	SPL2:5 16	SPL3:5 16	SPL3:5 16	SPL3:5 16	SPL4:5 16	SPL4:5 16	SPL4:5 16
H	SPL1:6 32	SPL1:6 32	SPL1:6 32	SPL2:6 32	SPL2:6 32	SPL2:6 32	SPL3:6 32	SPL3:6 32	SPL3:6 32	SPL4:6 32	SPL4:6 32	SPL4:6 32

BLK Blank CTL Control SPLB Sample blank SPL Sample

Microplate outline example for IC₅₀ determination.

- 1.1 Add  80 µL ultrapure water to blank (BLK)
- 1.2 Add  40 µL ultrapure water to control (CTL) and sample blank (SPLB)
- 1.3 Add  40 µL sample dilution to sample (SPL) and sample blank (SPLB)
- 1.4 Add  40 µL ACE work solution to control (CTL) and sample (SPL)
- 1.5 Add  160 µL Substrate solution to control all wells
- 1.6 Incubate at  37 °C  00:30:00 and read fluorescence with **350 nm excitation** wavelength and **420 nm emission** wavelength


30m

Result treatment

- 2 Inhibitory activity is expressed as the peptide concentration required to inhibit the original ACE activity by 50% (IC_{50}). The formula applied to calculate the percentage of ACE-inhibitory is:

$$iACE (\%) = ((F_{CTL} - F_{BLK}) - (F_{SPL} - F_{SPLB})) * \frac{100}{F_{CTL} - F_{BLK}}$$

Non-linear fitting to the data is performed to calculate the IC_{50} values, using the 5 Parameter curve fit method and then Interpolating to 50.



Example of a typical inhibition curve as a function of dilution factor.



Citations

Sentandreu, M. Á., & Toldrá, F.. A rapid, simple and sensitive fluorescence method for the assay of angiotensin-I converting enzyme

<https://doi.org/10.1016/j.foodchem.2005.06.006>

Coscueta, E. R., Campos, D. A., Osório, H., Nerli, B. B., & Pintado, M.. Enzymatic soy protein hydrolysis: A tool for biofunctional food ingredient production

<https://doi.org/10.1016/j.fochx.2019.100006>