



Apr 15, 2019

aCGH Hybridization & wash

 [PLOS One](#)

DOI

dx.doi.org/10.17504/protocols.io.zj3f4qn

Hendrik F HFvE Van Essen¹

¹Amsterdam UMC



Hendrik F HFvE Van Essen

Amsterdam UMC

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.zj3f4qn>

External link: <https://doi.org/10.1371/journal.pone.0223827>

Protocol Citation: Hendrik F HFvE Van Essen 2019. aCGH Hybridization & wash. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.zj3f4qn>

Manuscript citation:

Vincenten JPL, Essen HFv, Lissenberg-Witte BI, Bulkman NWJ, Krijgsman O, Sie D, Eijk PP, Smit EF, Ylstra B, Thunnissen E (2019) Clonality analysis of pulmonary tumors by genome-wide copy number profiling. PLoS ONE 14(10): e0223827. doi: [10.1371/journal.pone.0223827](https://doi.org/10.1371/journal.pone.0223827)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 27, 2019

Last Modified: April 15, 2019

Protocol Integer ID: 21851

Keywords: acgh hybridization, wash

Troubleshooting

aCGH Hybridization

- 1 Turn on hybridization oven  65 °C
Turn on heat block  95 °C
Turn on heat block  37 °C
- 2 Mix Cy3 and Cy5 samples
- 3 Mix the following components in 1.5 centrifuge tubes
 - * Cy3/Cy5 gDNA mixture 39 µl
 - * Cot-1 DNA (1.3 µg/µl) 5 µl
 - * Agilent 10x blocking agent 11 µl
 - * Agilent 2x Hybridization buffer 55 µl
 - Total 110 µl
- 4 Mix sample by pipetting, by careful not to create any bubbles.
- 5 Spin down with maximum speed for a few seconds.
- 6 Heat for  00:03:00 at  95 °C and incubate immediately thereafter for  00:30:00 at  37 °C .
- 7 Leave samples in heat block, spin down for  00:01:00 at maximum speed.
- 8 Load gasket slide in hybridization chamber and dispense 100 µl hyb mix in a drag and drop manner (use a 200 µl pipette, be careful not to create any bubbles)
- 9 Place microarray slide – active side down – and assemble chamber.



- 10 Check for immobile air bubbles by rotating the assembled hybridization chamber. If necessary, slam side of the chamber on table to loosen the air bubbles.
- 11 Place assembled chamber in the hybridization oven ( 65 °C at 20 rpm) and make sure to balance with empty chambers.
- 12 Hybridize for  24:00:00 ,  65 °C at 20 rpm.

aCGH wash

- 13 Turn on Agilent Slide Scanner
- 14 Wash buffer dishes with distilled water and blow dry with pressurized air.
- 15 Add buffers to dishes and have tin foil covering ready.-
Wash a maximum of 5 Agilent slides at the same time.

Dish	Method	Buffer	Temp (°C)	Time (min)
#1	Disassembly	Buffer 1	RT	-
#2	1st wash step	Buffer 1	RT	5 (stirrer)
#3	2nd wash step	Buffer 2	37	1 (stirrer & hot plate)
#4	Dry step	Acetonitrile	RT	1
- 16 Separate gasket slide from microarray slide in dish #1 (buffer 1)
- 17 Place slides in holding rack and was in dish #2 for  00:05:00 , rotation of stirrer is switched ON!
- 18 Cover dish with tin foil.
- 19 Start the "*Scan Control Software*"
For Agilent 180k, select the "AgilentG3_CGH" Option



- 20 Move holding rack to dish #3 (buffer 2) and wash for  00:01:00 , rotation of stirrer and heating is switched ON!
- 21 Move holding rack to dish #4 (acetonitrile) leave for  00:01:00 .
- 22 Remove slide gently with tweezers, thus preventing the formation of droplets.
- 23 Place slide in holder and load into the carousel of the scanner
Agilent Slide Scanner scans through the back of the slide (this side must contain the numeric barcode)
- 24 Click on "Set Values" to confirm all settings
- 25 Click Scan Slot m-n on the Scan Control main window
The letter m represents the Start slot where your first slide is located, and the letter n represents the End slot where your last slide is located.
- 26 Start the scan
- 27 Open "feature extraction" and check TIFF file