

Jul 03, 2023

Actin flow

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8pmb8g2w/v1>

Juan.Gonzalez¹

¹University of Glasgow

Centre for the Cellular M...



Juan.Gonzalez

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8pmb8g2w/v1>

Protocol Citation: Juan.Gonzalez 2023. Actin flow . protocols.io <https://dx.doi.org/10.17504/protocols.io.5jyl8pmb8g2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: July 01, 2023

Last Modified: July 03, 2023

Protocol Integer ID: 84352

Keywords: actin flow, actin

Abstract

Actin flow



Experimental procedure

- 1 Transfect the cells with the LifeAct-GFP plasmid following the protocol "Cells electroporation for cell transfection with NEON system"
- 2 The next day, trypsinize the cells as usual
- 3 Seed the cells on glass (i.e. gels on glass bottom petri dishes or the glass itself)

Note

The cell density depends on the cell type used. i.e. for C2C12 myoblasts, we seed 10000 cells/cm²

- 4 After  03:00:00 , image the cells with a confocal microscope. Use the 40x objective. If the CO₂ cannot be controlled, use CO₂-independent medium. The channel would be eGFP. 3h
- 5 Focus the cell in a focal plane where the actin cytoskeleton can be seen
- 6 We perform a time series of 120 images every  00:00:02 ( 00:04:00 in total per cell) 4m 2s

Data analysis

- 7
Open the image sequence with Fiji software

8

Expected result



9

Analyze, Set measurements, select "Bounding rectangle" and "Min & Max gray value"

10

Check where the cell cytoskeleton moves



11

Draw a line with the wand tool crossing the movement

12

Analyze, Multi kimograph, multi kymograph, linewidth: 1

13


A black & white gradient pops up. White lines indicate cell movement. The higher angle indicates more movement

14

Draw a line with the wand tool in the white line seen in the kymograph

15

Expected result



Data analysis from kymograph with Fiji software

- 16 Press "M" and copy & paste the measurements in an Excel file
- 17 In the Excel file, split the width by the length. This new value means flow in "pixels per frame"
- 18 Divide the pixels/frame value obtained by 2. This new value means flow in "pixels per second"
- 19 Split the "pixels/s" value by the "pixels/micron" ratio of the objective used (i.e. 3.2055 in the LSM900 confocal microscope) and multiply by 1000. This final value is the actin flow in nm/s