

Oct 22, 2021

Version 9

Adult mouse kidney dissociation (on ice) V.9

DOI

<https://dx.doi.org/10.17504/protocols.io.bzd9p296>



Andrew Potter¹

¹Cincinnati Children's Hospital Medical Center



Andrew Potter

CCHMC

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bzd9p296>

Protocol Citation: Andrew Potter 2021. Adult mouse kidney dissociation (on ice). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bzd9p296> Version created by **Andrew Potter**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: October 22, 2021

Last Modified: January 07, 2026

Protocol Integer ID: 54433

Keywords: kidney, dissociation, single cell, adult mouse kidney dissociation, mouse kidney dissociation, development for p1 mouse kidney, p1 mouse kidney, major kidney cell type, ice, additional enzymatic digestion, artifact gene expression, cell viability

Abstract

Protocol for adult (8-10 week) mouse kidney dissociation performed on ice to reduce artifact gene expression. The protocol is based on our protocol published in *Development* for P1 mouse kidney, however this protocol includes two layers to provide additional enzymatic digestion. Cell viability is ~80% with major kidney cell types represented.

Attachments



Adult mouse kidney c...

17KB

Guidelines

***Bacillus Licheniformis* Enzyme Mix (2 × 1 mL per 25 mg tissue):**

100 µL b. lich 100 mg/mL (10 mg/mL final) (Sigma, P5380)

5 µL 1 M CaCl₂ (5 mM)

5 µL DNase (125 U) (StemCell, #07469)

890 µL DPBS (no Ca, Mg) ThermoFisher (cat. #14190)

Preparing enzymes:

The enzymes are made up in DPBS (#14190). They are aliquoted and stored at -80 °C. FBS (for making the FBS/PBS) is heat-inactivated and sterile-filtered.

Required reagents:

Red Blood Cell Lysis Buffer - Sigma (R7757)

Required Equipment & Consumables:

Thermomixer

Centrifuges for 1.5 mL and 15 mL conicals (MLS)

Pipettes and pipet tips (MLS)

15 ml Conicals (MLS)

1.5 mL tubes (MLS)

30 µM filter (MLS)

Petri dishes (MLS)

Razor blades (MLS)

Ice bucket w/ice (MLS)

Hemocytometers - InCyto Neubauer Improved (DHC-NO1-5)

The protocol workflow is as follows:

A. Isolate Kidney

B. First layer

C. Second layer

D. Preparing cells for single cell analysis



Materials

MATERIALS

⊗ DPBS (no Ca, no Mg) **Thermofisher** Catalog #14190144

⊗ RBC Lysis Buffer **Merck MilliporeSigma (Sigma-Aldrich)** Catalog #R7757

Before start

- Prepare enzyme mixes and leave on ice.
- Cool centrifuges to 4 °C.
- Isolate and transport tissue in ice-cold DPBS.



Isolate kidney

1 Perfuse kidneys to remove RBC. Extract and isolate kidneys in ice-cold PBS. Leave kidneys in ice-cold PBS until ready to dissociate (Up to 1 hr).

2 Coarsely mince tissue in PBS.

00:02:00 mince on ice

Layer 1

3 Weigh out 25 mg coarsely minced tissue for each set of kidneys (remove PBS before weighing).

25 μ L minced kidney tissue

4 Continue mincing kidneys on top of petri dish, on ice, using razor blade in small vol. (~50 μ L) PBS. (1-2 min) until fine paste.

50 μ L PBS

5 Prepare a separate 1 mL aliquot of B. Lich enzyme mix for each set of adult kidneys (prepare on ice). Use p200 w/cut tip to transfer minced kidney tissue from petri dish to tube of enzyme mix.

1 mL B. Lich enzyme mix

6 Incubate digest mix for 10 min on ice with trituration and shaking. Triturate 15 strokes using 1 mL pipet set to 600 μ L every 2 min; shake vigorously every min.

13m

00:01:00 shake vigorously

00:02:00 triturate 15X

00:10:00 incubate on ice

7 After 10 min, let tissue chunks settle on ice for 1 min. Save supernatant (700 μ L, 70%) and apply to 30 μ M filter on 15 mL conical. Rinse filter w/5 mL 10% FBS/PBS. Leave filter on 15 mL conical for the next layer.

1m



⌚ 00:01:00 let tissue chunks settle

🧴 700 μ L save supernatant

🧴 5 mL rinse filter with 5 mL 10% FBS/PBS

Layer 2

- 8 Add additional 700 μ L B. Lich enzyme mix to residual tissue chunks. Continue incubating on ice 15-20 min with shaking and trituration, until tubules and glomeruli are fully broken up.

23m

🧴 700 μ L B. lich enzyme mix

⌚ 00:01:00 shake vigorously

⌚ 00:02:00 triturate 15X

⌚ 00:20:00 incubate on ice

- 9 Once digestion is adequate (tubules / glomeruli are broken up), triturate and add entire digest mix to same 30 μ M filter as used in previous step on 15 mL conical.

- 10 Spin 15 mL conical with combined flow-through from layer 1 and layer 2 (isolated cells) 300 g for 5 min at 4 °C.

5m

⌚ 00:05:00 300 g at 4 °C

Prepare cells for 10X Chromium / scRNA-Seq

- 11 Discard supernatant. If necessary, perform RBC removal according to manufacturers instructions.

- 12 Re-suspend cells in 10 mL 10% FBS/PBS. Spin 300 g for 5 min at 4 °C.

🧴 10 mL 10% FBS/PBS

⌚ 00:05:00 spin 300 g

- 13 Discard supernatant. Re-suspend cells in 500 μ L 10% FBS/PBS (or other compatible buffer). Analyze cell viability and concentration using trypan blue dye exclusion with a hemocytometer. Adjust cell concentration to 700-1,200 cells/ μ L for 10X single cell 3'v3.1.



🧪 500 μ L 10% FBS/PBS