

Feb 01, 2023

Version 3

Adult mouse skin dissociation protocol (on ice) V.3

DOI

dx.doi.org/10.17504/protocols.io.n92ld391og5b/v3

Andrew Potter¹

¹CCHMC

Human Cell Atlas Metho...



Andrew Potter

CCHMC

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ld391og5b/v3>

Protocol Citation: Andrew Potter 2023. Adult mouse skin dissociation protocol (on ice). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ld391og5b/v3> Version created by **Andrew Potter**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 01, 2023

Last Modified: February 01, 2023

Protocol Integer ID: 76216

Keywords: CAP, skin, dissociation, single cell, adult mouse skin dissociation protocol, mouse skin, bacillus licheniformis protease cocktail, ice, layers of digestion, protease, mouse, digestion, cell, protocol

Abstract

This protocol was developed to dissociate adult (8-10 wk) mouse skin "on ice". It utilizes two layers of digestion with a *Bacillus licheniformis* protease cocktail, combined with mechanical disruption from a dounce homogenizer. The cell yield is 3000 cells/mg.

Attachments



[Skin adult mouse dis...](#)

99KB



Guidelines

***Bacillus Licheniformis* Enzyme Mix (1 mL per 23 mg tissue):**

100 μ L b. lich 100 mg/mL (10 mg/mL final) (Sigma, P5380)

1 μ L 0.5 M EDTA (Sigma, A8806)

899 μ L DPBS (no Ca, Mg) ThermoFisher (cat. #14190)

Preparing enzymes:

The enzyme is made up in DPBS (#14190). It is aliquoted and stored at -80 °C at 100 mg/mL in 100 μ L aliquots..

Reagents

Enzymes and BSA are made up in DPBS (no Ca, no Mg) from Thermo Fisher (14190).

Bovine Serum Albumin - Sigma (A8806).

Hypothermosol FRS

Required supplies:

2 mL dounce homogenizer – Bellco (1984-10002)

Centrifuge for 1.5 mL, 15 mL conicals

Pipettes and pipet tips

15 ml Conicals (MLS)

1.5 mL tubes (MLS)

30 μ M filters - Miltenyi (130-098-458)

Petri dishes (MLS)

Razor blades (MLS)

Ice bucket w/ice


Hemocytometers - InCyto Neubauer Improved (DHC-NO1-5)

Troubleshooting




Isolating tissue


- 1 After euthanizing mouse, remove hair using Nair: dab with Nair, wait 30 secs, wipe with wet paper towel.
- 2 Isolate tissue and place in ice-cold hypothermosol.
- 3 Scrape off underlying layer of fatty / connective tissue using scalpel before proceeding.
- 4 Mince skin tissue thoroughly on petri dish on ice for 3-4 min on ice into 1-mm³ pieces using razor blade while manipulating tissue with forceps – you will need to use grinding motion and vigorously break up tissue.

 00:04:00 mincing

1st digest layer

- 5 Place 23 mg minced tissue into 1 mL B. Lich enzyme cocktail. Incubate on ice.

 23 mg minced tissue
- 6 Shake every min; triturate 10x every 2 min with p1000 w/tip cut (start triturating at 2 min) for 20 min.

 00:20:00 digest

 00:01:00 shake  00:02:00 triturate 10X

- 7 After 20 mins of triturating on ice, use pipet to transfer digest mix to 2 mL dounce homogenizer. Use 10 strokes of Pestle A every 2 min (4 series total, 8 min). Digest mix should become turbid.

 00:02:00 dounce homogenize  00:08:00 digest using dounce

- 8 Transfer back to 1.5 mL tube using 1 mL serological pipet. Mix thoroughly and allow to settle on ice 2 min.



00:02:00 settle on ice

- 9 Save 70% (700 μ L) of supernatant, leaving chunks at the bottom of the tube; apply to 30 μ M filter on 15 mL conical. Rinse filter w/5 mL ice-cold PBS/BSA 0.04%. Save flow through on ice and keep filter on tube for 2nd layer.

 700 μ L save supernatant

5 mL rinse filter w/ice-cold PBS/BSA 0.04%

2nd digest layer

- 10 Add additional 1 mL b. Lich enzyme mix to residual tissue chunks.

1 mL b. lich enzyme mix

- 11 Triturate 10x every 2 min, shake every min while incubating on ice for 20 additional mins. (50 min. total digest time).

00:20:00 additional digest time

00:02:00 triturate 10x

00:01:00 shake every min.

- 12 Transfer entire volume to same 30 μ M filter on 15 mL conical. Rinse with additional 5 mL ice-cold PBS/BA 0.04%.

5 mL ice-cold PBS/BSA 0.04%

Preparing cells for single cell analysis

- 13 Centrifuge at 300 g for 5 min at 4 °C. Remove supernatant & re-suspend in 100 μ L PBS/BSA 0.04%. Examine using hemocytometer with trypan blue.

00:05:00 centrifuge at 300 g

 100 μ L re-suspend in ice-cold PBS/BSA 0.04%

- 14 Adjust concentration to 1,000 cells/ μ L for Chromium or 100 cells/ μ L for DropSeq.

