



May 10, 2018 Version 1

# Adult Mouse Spleen Dissociation (On ice) V.1

DOI

[dx.doi.org/10.17504/protocols.io.p2ddqa6](https://dx.doi.org/10.17504/protocols.io.p2ddqa6)

Andrew Potter<sup>1</sup>

<sup>1</sup>CCHMC

Human Cell Atlas Metho...



Andrew Potter

CCHMC

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.p2ddqa6](https://dx.doi.org/10.17504/protocols.io.p2ddqa6)

**Protocol Citation:** Andrew Potter 2018. Adult Mouse Spleen Dissociation (On ice). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.p2ddqa6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 09, 2018

**Last Modified:** May 10, 2018

**Protocol Integer ID:** 12069

**Keywords:** spleen

## Abstract

Protocol used to dissociate adult (8-10 wk) mouse spleen into single cells. Attained >95% viability, a variety of cell sizes, and ~10 million cells from 12 mg tissue.



## Guidelines

### **Collagenase Enzyme Mix (two tubes, 1 mL each)**

7.5 mg/mL Collagenase A (Sigma, 10103578001)  
7.5 mg/mL Collagenase Type 4 (Worthington, CLS-4)  
100 µg/mL soybean trypsin inhibitor (Sigma, 10109886001)  
125 U DNase (Applichem, A3778)  
5 mM CaCl<sub>2</sub>  
740 µL DPBS (no Ca, Mg)

---

+12 mg chopped spleen / tube

## Materials

### MATERIALS


 Red Blood Cell Lysis Buffer Hybri-Max Merck MilliporeSigma (Sigma-Aldrich) Catalog #R7757


## Before start


- Set centrifuges to 4° C.
- Make two tubes of 1 mL enzyme mix.
- Make ~25 mL of DPBS/0.04% BSA



- 1 Chop tissue coarsely (~30 secs) using razor blade on petri dish, on ice.
- 2 Add 12 mg chopped tissue to 1 mL enzyme mix.
- 3 Incubate tube on ice for 10 minutes. Triturate 10X every 2 mins and shake every min.  

 00:10:00 Incubate on ice
- 4 After 10 mins of digestion, let tissue chunks settle for 1 min on ice & remove 80% of supernatant with released cells & filter using 70  $\mu$ M filter on 50 mL conical, on ice. Rinse filter with 5 mL ice-cold PBS/0.04% BSA. Leave filter and 50 mL conical on ice, it will be used for the steps as well.
- 5 Add additional 1 mL enzyme mix to tissue chunks.
- 6 Continue to triturate 10x every 2 minutes and shake every minute while incubating on ice, for 10 additional minutes.  

 00:10:00 Incubate on ice
- 7 Triturate and add entire volume of cell digestion to 70  $\mu$ M filter on 50 mL conical. Rinse w/5 mL ice-cold PBS/0.04% BSA.
- 8 Transfer flow-through to 15 mL conical. Spin 650 g for five minutes at 4 °C. After spin, remove supernatant (down to ~100  $\mu$ L).
- 9 Perform RBC lysis: add 1 mL RBC lysis buffer to cells and triturate 10X. Let sit 3 minutes on ice. Add 10 mL ice-cold PBS/BSA 0.04%  

 00:03:00 incubate on ice
- 10 Spin 650 g for 5 mins at 4 °C. Remove supernatant and re-suspend in 1 mL PBS/BSA 0.04%.